






























Gift Set Student Microscope							
No.	Cata. No.	Product	Specification	Package			Picture
1	A11.1513	Toy Biological Microscope Gift Set	Monocular, 10x, Triple Nosepiece, 10x,40x,90x, Clip Stage, LED+Mirror, With Plastic Suitcase & Accessories Kit For Experiment Kit Include: Packing Box Size 35x42x12cm, Shrimp Hatchery, Brine ShrimpEggs, Specimen Slicer, Petri Dish, Plastic Tweezers, Plastic Scalpel, Plastic Spatula, Needle, Plastic Stirring Rod, Graduated Cylinder, Eosin, Gum Media (bottle), Sea Salt, Prepared Slides*5, Blank Slides*7, Slide Covers*7, Statistical Slide Covers*7, Slide Labels*7, Collecting Vials*2, Spare LED Bulb	5 pcs / 1 ctn 55*40*36cm, 10.5/8Kgs			
2	A11.1124	Student Microscope, Cordless Dual Light, Plastic Body	Monocular, WF10x, Triple, Achromatic 4x10x40x, Plain Stage 90x90mm, Disc Diaphragm, Up/Down LED, 3xAAA Batteries	10 pcs / 1 ctn 81*45*37cm, 16/15Kgs			
3	A11.1526	Student Biological & Stereo Microscope	Biological & Stereo Microscope, Monocular, WF10x/18mm, Zoom Lens 2-10x, Round Stage Dia.82mm, Disc Diaphragm, Up/Bottom LED, Cordless 3xAAA Batteries	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
4	A11.1528	Student Microscope	Monocular, WF10x, WF25x, Triple, Achromatic 4x10x40x, Plain Stage 90x90mm, Disc Diaphragm, LED, Slide Box 5 Blank + 5 Prepared	10 pcs / 1 ctn 74*47*34cm, 12/11Kgs			
5	A11.1529	Student Microscope, Cordless Dual Light, Plastic Body	Monocular, WF10x, Triple, Achromatic 4x10x40x, Plain Stage 90x90mm, Disc Diaphragm, Up/Down LED, 3xAAA Batteries And DC In Power Supply	10 pcs / 1 ctn 74*47*34cm, 12/11Kgs			
6	A11.1509-M	Digital Student Biological Microscope Gift Set	Microscope + Foam Carton Monocular, Eyepiece WF10x, 16x, Triple, Achromatic 4x10x40x, Clip Stage, 6 Holes Disc Diaphragm With 5 Filters Built-in, Up + Down LED, Foam+Carton Package	10 pcs / 1 ctn 68*39*34cm, 14/11Kgs			
7	A11.1509	Digital Student Biological Microscope Gift Set	Microscope + Plastic Suitcase + Acce. Kit +100K Camera Monocular, Eyepiece WF10x, 16x, 2x Barlow Lens, Triple, Achromatic 4x10x40x, Clip Stage, 6 Holes Disc Diaphragm With 5 Filters Built-in, Up + Down LED, Foam+Carton Package	4 pcs / 1 ctn 57*42*29cm, 12/11Kgs			
8	A11.1510-M	Digital Student Biological Microscope Gift Set	Microscope + Foam Carton Monocular, Eyepiece WF10x, 16x, Triple, Achromatic 4x10x40x, Clip Stage, 6 Holes Disc Diaphragm With 5 Filters Built-in, Up + Down LED, Foam+Carton Package	10 pcs / 1 ctn 68*39*34cm, 14/0Kgs			









No.	Cata. No.	Product	Specification	Package			Picture
9	A11.1510	Digital Student Biological Microscope Gift Set	Microscope + Plastic Suitcase + Acce. Kit +100K Camera Monocular, Eyepiece WF10x, 16x, 2x Barlow Lens, Triple, Achromatic 4x10x40x, Clip Stage, 6 Holes Disc Diaphragm With 5 Filters Built-in, Up + Down LED	4 pcs / 1 ctn 57*42*29cm, 12/0Kgs			
10	A11.1511-M	Digital Student Biological Microscope Gift Set	Microscope + Foam Carton Monocular, WF10x,16xx, Triple, 4x10x40x, Clip Stage 90x90mm, Bottom LED	8 pcs / 1 ctn 71*49*34cm, 16/0Kgs			
11	A11.1511	Digital Student Biological Microscope Gift Set	Monocular, WF10x,16xx, Triple, 4x10x40x, Clip Stage 90x90mm, Up/Bottom LED, Including Gift Box, Accessories, Including Digital Eyepiece Camera	4 pcs / 1 ctn 64*28*42cm, 16/0Kgs			
12	A11.1512-M	Digital Student Biological Microscope Gift Set	Microscope + Foam Carton Monocular, Eyepiece WF10x, 16x, Triple, Achromatic 4x10x40x, Clip Stage, 6 Holes Disc Diaphragm With 5 Filters Built-in, Up + Down LED, Foam+Carton Package	8 pcs / 1 ctn 71*49*34cm, 16/15Kgs			
13	A11.1512	Digital Student Biological Microscope Gift Set	Microscope + Plastic Suitcase + Acce. Kit Monocular, Eyepiece WF10x, 16x, 2x Barlow Lens, Triple, Achromatic 4x10x40x, Clip Stage, 6 Holes Disc Diaphragm With 5 Filters Built-in, Up + Down LED, 5.5V Transformer, Plastic	2 pcs / 1 ctn 42*36*28cm, 9/8Kgs			
Primary Level Student Microscope							
14	A11.1505	Student Biological Microscope, Round Stage	Monocular Straight Tube, H5x,10x,15x, Triple, 5x10x45x, Round Stage Dia.120mm+ Clip, Disc Diaphragm, Coarse & Fine Focusing, Plan-Concave Mirror Dia.50mm	4 pcs / 1 ctn 60*36*39cm, 14.5/11.4Kgs			
15	A11.1508-01	Student Biological Microscope, Vertical Monocular Head	Monocular Vertical,10X,12.5X, Triple 4x10x40x, Clip Stage110x120mm, 5 Holes Disc Diaphragm, Coarse & Fine Focusing, Mirror	4 pcs / 1 ctn 53*43*39cm, 14.5/13Kgs			
16	A11.1508-02	Student Biological Microscope, Vertical Monocular Head	Monocular Vertical,10X,16X, Triple 4x10x40x, Clip Stage110x120mm, 5 Holes Disc Diaphragm, Coarse & Fine Focusing, Mirror	4 pcs / 1 ctn 53*43*39cm, 15/14Kgs			
17	A11.1508-03	Student Biological Microscope, Vertical Monocular Head	Monocular Straight Tube, H5x,10x, Triple, 10x40x100, Vernial Scale Stage 120x120mm, Coarse & Fine Focusing, N.A.1.25 Condenser + Iris Diaphragm, Plan-Concave Mirror Dia.50mm	4 pcs / 1 ctn 53*43*39cm, 15/14Kgs			
18	A11.1508-04	Student Biological Microscope, Vertical Monocular Head	Monocular Straight Tube, H5x,10x12.5x, Triple, 10x40x100, Vernial Scale Stage 120x120mm, Coarse & Fine Focusing, N.A.1.25 Condenser + Iris Diaphragm, Plan-Concave Mirror Dia.50mm	4 pcs / 1 ctn 53*43*39cm, 15/14Kgs			






No.	Cata. No.	Product	Specification	Package			Picture
19	A11.1508-05	Student Biological Microscope, Vertical Monocular Head	Monocular Straight Tube, H5x,10x15x, Triple, 10x40x100, Vernial Scale Stage 120x120mm, Coarse & Fine Focusing, N.A.1.25 Condenser + Iris Diaphragm, Plan-Concave Mirror Dia.50mm	4 pcs / 1 ctn 53*43*39cm, 15/14Kgs			
20	A11.1508-06	Student Biological Microscope, Vertical Monocular Head	Monocular Straight Tube, H5x,10x16x, Triple, 10x40x100, Vernial Scale Stage 120x120mm, Coarse & Fine Focusing, N.A.1.25 Condenser + Iris Diaphragm, Plan-Concave Mirror Dia.50mm	4 pcs / 1 ctn 53*43*39cm, 15/14Kgs			
21	A11.1321	Student Biological Microscope	Monocular, WF10x, Triple, Achromatic 4x10x40, Clip Stage 95x95mm, Single Lens Condenser + Disc Diaphragm, Coarse & Fine Focusing, LED Light Brightness Adjustable, External AD Adapter (or 3xAA Battery)	1 pcs / 1 ctn 25*20*50cm, 2/1.5Kgs			
22	A11.1322	Student Biological Microscope, Binocular, Dual LED Light	Binocular, WF10x, Triple, Achromatic 4x10x40, Vernial Scale Stage 95x95mm/50x30mm, Disc Diaphragm, Coarse & Fine Focusing, Dual LED Light Brightness Adjustable, External AD Adapter (or 3xAA Battery)	1 pcs / 1 ctn 25*20*50cm, 2/1.5Kgs			
23	A11.1326	Student Biological Microscope, Round Stage	Monocular, WF10x, Triple, Achromatic 4x10x40x, Round Quick Move Stage, Disc Diaphragm, Coarse Focusing, Up/Down LED, Adapter/Battery	10 pcs / 1 ctn 88*49*36cm, 17/15Kgs			
24	A11.1518-M	Student Biological Microscope	Monocular, WF10x, Triple, Achromatic 4x10x40x, Plain Stage 90x90mm, Disc Diaphragm, LED, 3xAA Batteries And DC In Power Supply	8 pcs / 1 ctn 71*49*34cm, 16/14Kgs			
25	A11.1518-B	Student Biological Microscope	Binocular, WF10x, Triple, Achromatic 4x10x40x, Plain Stage 90x90mm, Disc Diaphragm, LED, 3xAA Batteries And DC In Power Supply	8 pcs / 1 ctn 71*49*34cm, 16/14Kgs			
26	A11.1143	Student Microscope	Monocular, WF10x Pointer, Triple, Achromatic 4x,10x,40x, Clip Stage 95x95mm, Coaxial Coarse & Fine Focusing, Uppler LED, Bottom LED Brightness Adjustable, Rechargeable	1 pcs / 1 ctn 28.5*21*41cm, 2.5/1.8Kgs			
27	A11.1123-M	Biological Microscope	Monocular, WF10x Pointer, Quarduple, Achromatic 4x10x40x100x, Double Layer 100x90mm, N.A.1.25, LED, Incl. Rechargeable Batteries	1 pcs / 1 ctn 31*21*42cm, 4/3Kgs			
28	A11.1123	Biological Microscope	Binocular, WF10x Pointer, Quarduple, Achromatic 4x10x40x100x, Double Layer 100x90mm, N.A.1.25, LED, Incl. Rechargeable Batteries	1 pcs / 1 ctn 31*21*42cm, 4/3Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
29	A11.1133-A	Biological Microscope	Monocular, WF10x, Triple, Achromatic 4x10x40x, Plain Stage 100*90mm, LED	2 pcs / 1 ctn 48*32*45cm, 16/14Kgs			
30	A11.1133-B	Biological Microscope	Monocular, WF10x, Quaduple, Achromatic 4x10x40x100x, Plain Stage With Vernial Scale, 100*90mm, LED	2 pcs / 1 ctn 48*32*45cm, 16/14Kgs			
31	A11.1133-C	Biological Microscope	Sliding Binocular, WF10x, Quaduple, Achromatic 4x10x40x100x, Plain Stage With Vernial Scale, 100*90mm, LED	2 pcs / 1 ctn 48*32*45cm, 16/14Kgs			
32	A11.1133-D	Biological Microscope	Sedentopf Binocular, WF10x, Quaduple, Achromatic 4x10x40x100x, Mechanical Stage 115*125mm, LED	2 pcs / 1 ctn 48*32*45cm, 16/14Kgs			
33	A11.1009-D	Student Biological Microscope	Monocular, WF10x, Triple, Achromatic 4x10x40x, Mechanical Stage 110x120mm/70x20mm, Abbe 1.2 Rack & Pinion, Coaxial Coarse&Fine, Incandescent 15W	1 pcs / 1 ctn 28*21*39cm, 3/0Kgs			
34	A11.1009-E	Student Biological Microscope	Binocular, WF10x, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 110x120mm/70x20mm, Abbe 1.2 Rack & Pinion, Coaxial Coarse&Fine, LED 1W	1 pcs / 1 ctn 28*21*39cm, 3/0Kgs			
35	A11.1020-M	Biological Microscope	Monocular, WF10x/16, WF25x/9, Triple, Achromatic 4x10x40x, N.A. 0.65+ Disc Diaphragm, Plain Stage 95*95mm, Coarse Coaxial & Fine Focusing, Up LED, Bottom LED Brightness Adjustable, 3xAAA Batteries/AC Power	1 pcs / 1 ctn 28*19*39cm, 3/2.5Kgs			
36	A11.1020-B	Biological Microscope	Binocular, WF10x/16, WF25x/9, Triple, Achromatic 4x10x40x, N.A. 0.65+ Disc Diaphragm, Plain Stage 95*95mm, Coarse Coaxial & Fine Focusing, Up LED, Bottom LED Brightness Adjustable, 3xAAA Batteries/AC Power	1 pcs / 1 ctn 28*19*39cm, 3/2.5Kgs			
37	A11.1303-B	Biological Microscope	Binocular, WF10x, Quadruple, Achromatic 4x10x40x100x, Double Layer 125x130/60x30mm, N.A.1.25, Halogen 6V20W	4 pcs / 1 ctn 69*45*47cm, 20/0Kgs			
38	A11.1303-T	Biological Microscope	Trinocular, WF10x, Quadruple, Achromatic 4x10x40x100x, Double Layer 125x130/60x30mm, N.A.1.25, Halogen 6V20W, Incl. 1x C-mount, Eyepiece Adapter	4 pcs / 1 ctn 69*45*47cm, 20/0Kgs			






No.	Cata. No.	Product	Specification	Package			Picture
39	A11.1313-B	Biological Microscope	Seidentopf Binocular, WF10x/18mm, Quarduple, Achromatic 4x10x40x100x, Double Layer Mechanical 125x130mm/30x60mm, Coaxial Coarse & Fine Focusing, 6V20W Halogen	2 pcs / 1 ctn 50*34*49cm, 12/0Kgs			
40	A11.1313-T	Biological Microscope	Seidentopf Trinocular, WF10x/18mm, Quarduple, Achromatic 4x10x40x100x, Double Layer Mechanical 125x130mm/30x60mm, Coaxial Coarse & Fine Focusing, 6V20W Halogen	2 pcs / 1 ctn 50*34*49cm, 13/0Kgs			
41	A11.1521-M1	Student Biological Microscope, Monocular	Monocular, WF10x/16mm, Triple, Ach. 4x10x40x, Coaxial, Plain Stage 120x120mm, N.A.0.65 Condenser + Disc Diaphragm, LED	2 pcs / 1 ctn 45*27*41cm, 7.5/6Kgs			
42	A11.1521-M	Microscope	Monocular, WF10x/16mm, Quarduple, Ach. 4x10x40x100x, Mechanical Stage 115x125mm, Halogen 6V20W/LED	2 pcs / 1 ctn 45*27*41cm, 7.5/5Kgs			
43	A11.1521-B	Microscope	Binocular, WF10x/16mm, Quarduple, Ach. 4x10x40x100x, Mechanical Stage 115x125mm, Halogen 6V20W/LED	2 pcs / 1 ctn 45*27*41cm, 7.5/6.4Kgs			
44	A11.1521-T	Microscope	Trinocular, WF10x/16mm, Quarduple, Ach. 4x10x40x100x, Mechanical Stage 115x125mm, Halogen 6V20W/LED, 1x C-Mount	2 pcs / 1 ctn 45*27*41cm, 7.5/6.5Kgs			
45	A11.5121-M	Microscope	Monocular Teaching, WF10x, WF16x, Quarduple, Ach. 4x10x40x100x, Mechanical Stage 115x125mm, Coaxial Focusing, LED, 0.3x C-Mount, With 10 Prepared Slides, With Experiment Kit, Pack In Hand Carry Case+Carton	1 pcs / 1 ctn 30*20*40cm, 4.5/4Kgs			
46	A11.5121-B	Microscope	Binocular, WF10x, WF16x, Quarduple, Ach. 4x10x40x100x, Mechanical Stage 115x125mm, Coaxial Focusing, LED, With 10 Prepared Slides, With Experiment Kit, Pack In Hand Carry Case+Carton	1 pcs / 1 ctn 30*20*40cm, 4.5/4Kgs			
47	A11.5121-T	Microscope	Trinocular, WF10x, WF16x, Quarduple, Ach. 4x10x40x100x, Mechanical Stage 115x125mm, Coaxial Focusing, LED, 0.3x C-Mount, With 10 Prepared Slides, With Experiment Kit, Pack In Hand Carry Case+Carton	1 pcs / 1 ctn 30*20*40cm, 4.5/4Kgs			
48	A11.1007-17W	Laboratory Biological Microscope XSZ107BN,	Sliding Binocular White, WF10x16x, Quaduple, Achromatic 4x10x40x100x, Double Layer Mechanical 140x140mm /50x75mm, Coaxial Coarse&Fine, Abbe 1.25, LED, Foam +Carton Packing	1 pcs / 1 ctn 33*28*43cm, 6/5Kgs			 HOT

No.	Cata. No.	Product	Specification	Package			Picture
49	A11.1007-17WT	Laboratory Biological Microscope XSZ107BN,	Sliding Trinocular White, WF10x16x, Quaduple, Achromatic 4x10x40x100x, Double Layer Mechanical 140x140mm /50x75mm, Coaxial Coarse&Fine, Abbe 1.25, 3W LED, Eyepiece Adapter + 1x C-Mount, Foam +Carton Packing	1 pcs / 1 ctn 33*28*43cm, 6/5Kgs			
50	A11.1535-M	Microscope	Monocular, WF10x/18mm, Quaduple, Achromatic 4x10x40x100x, Abbe 1.25, Double Layer Mechanical Stag 142*132mm, Coaxial Focusing, LED	1 pcs / 1 ctn 40*27*48cm, 6/5Kgs			
51	A11.1535-B	Microscope	Binocular Head, WF10x/18mm, Quaduple, Achromatic 4x10x40x100x, Abbe 1.25, Double Layer Mechanical Stag 142*132mm, Coaxial Focusing, LED	1 pcs / 1 ctn 40*27*48cm, 6/5Kgs			
52	A11.1535-T	Microscope	Trinocular, WF10x/18mm, Quaduple, Achromatic 4x10x40x100x, Abbe 1.25, Double Layer Mechanical Stag 142*132mm, Coaxial Focusing, LED, With 1x C-Mount	1 pcs / 1 ctn 40*27*48cm, 6/5Kgs			
53	A11.1535-V	Microscope	Teaching Head, WF10x/18mm, Quaduple, Achromatic 4x10x40x100x, Abbe 1.25, Double Layer Mechanical Stag 142*132mm, Coaxial Focusing, LED, With 1x C-Mount	1 pcs / 1 ctn 40*27*48cm, 6/5Kgs			
54	A11.1531-M	Microscope	Monocular, WF10x, Quaduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 110x125mm, 3W LED Adjustable, Power AC100V-240V+AA1.5V Battery 3pcs	1 pcs / 1 ctn 33*25.5*42cm, 5/4.2Kgs			
55	A11.1531-V	Microscope	Dual-View Head, WF10x, Quaduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 110x125mm, 3W LED Adjustable, Power AC100V-240V+AA1.5V Battery 3pcs	1 pcs / 1 ctn 33*25.5*42cm, 5/4.2Kgs			
56	A11.1531-B	Microscope	Binocular, WF10x, Quaduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 110x125mm, 3W LED Adjustable, Power AC100V-240V+AA1.5V Battery 3pcs	1 pcs / 1 ctn 33*25.5*42cm, 5/4.2Kgs			
57	A11.1531-T	Microscope	Trinocular, WF10x, Quaduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 110x125mm, 3W LED Adjustable, Power AC100V-240V+AA1.5V Battery 3pcs	1 pcs / 1 ctn 33*25.5*42cm, 5/4.2Kgs			
58	A11.1551-M	Student Biological Microscope	Monocular, WF10x, Quaduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 142x132mm, 6V20W Halogen	1 pcs / 1 ctn 34*25.5*49cm, 6/5.3Kgs			





No.	Cata. No.	Product	Specification	Package			Picture
59	A11.1551-B	Student Biological Microscope	Binocular, WF10x, Quarduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 142x132mm, 6V20W Halogen	1 pcs / 1 ctn 34*25.5*49cm, 6/5.3Kgs			
60	A11.1551-T	Student Biological Microscope	Trinocular, WF10x, Quarduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 142x132mm, 6V20W Halogen	1 pcs / 1 ctn 34*25.5*49cm, 6/5.3Kgs			
61	A11.1536-B	Student Biological Microscope	Binocular, WF10x/18mm, Quarduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 115x125mm, 3W LED, Rechargeable Battery Box	1 pcs / 1 ctn 34*25.5*49cm, 6/5.3Kgs			
62	A11.1536-T	Student Biological Microscope	Trinocular, WF10x/18mm, Quarduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 115x125mm, 3W LED, Rechargeable Battery Box, Digital Eyepiece Adapter	1 pcs / 1 ctn 34*25.5*49cm, 6/5.3Kgs			
63	A11.1537-B	Student Biological Microscope	Binocular, WF10x/18mm, Quarduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 115x125mm, 3W LED	1 pcs / 1 ctn 34*25.5*49cm, 6/5.3Kgs			
64	A11.1537-T	Student Biological Microscope	Trinocular/18mm, WF10x, Quarduple, Achromatic 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 115x125mm, 3W LED, Digital Eyepiece Adapter	1 pcs / 1 ctn 34*25.5*49cm, 6/5.3Kgs			
65	A11.1170-B	Biological Microscope	Binocular, WF10x, WF25x, Quarduple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage 130x130mm, Abbe N.A.1.25, LED	1 pcs / 1 ctn 41*27*17cm, 5/4Kgs			
66	A11.1170-T	Biological Microscope	Trinocular, WF10x, WF25x, Quarduple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage 130x130mm, Abbe N.A.1.25, LED	1 pcs / 1 ctn 41*27*17cm, 5/4Kgs			
67	A11.1171-B	Biological Microscope	Binocular, WF10x, WF25x, Quarduple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage 130x130mm, Abbe N.A.1.25, LED	1 pcs / 1 ctn 41*27*17cm, 5/4Kgs			
68	A11.1171-T	Biological Microscope	Trinocular, WF10x, WF25x, Quarduple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage 130x130mm, Abbe N.A.1.25, LED	1 pcs / 1 ctn 41*27*17cm, 5/4Kgs			











No.	Cata. No.	Product	Specification	Package			Picture
69	A12.1303-M	Laboratory Microscope	Monocular, WF10x/18, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 125*125mm, LED Brightness Adjustable	2 pcs / 1 ctn 44*37*47cm, 8/0Kgs			
70	A12.1303-S	Laboratory Microscope	Dual Viewing 45°+45°, WF10x/18, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 125*125mm, LED Brightness Adjustable	2 pcs / 1 ctn 44*37*47cm, 9/0Kgs			
71	A12.1303-V	Laboratory Microscope	Dual Teaching 30°+90°, WF10x/18, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 125*125mm, LED Brightness Adjustable	2 pcs / 1 ctn 44*37*47cm, 9/0Kgs			
72	A12.1303-B	Laboratory Microscope	Binocular, WF10x/18, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 125*125mm, LED Brightness Adjustable	2 pcs / 1 ctn 44*37*47cm, 10/0Kgs			
73	A12.1303-T	Laboratory Microscope	Trinocular, WF10x/18, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 125*125mm, LED Brightness Adjustable	2 pcs / 1 ctn 44*37*47cm, 10.5/0Kgs			











High Level Student Microscope











74	A11.0913-M	Biological Microscope	Monocular, WF Plan10X/18mm, Quadruple, Achromatic Objective 4x10x40x, Fixed Single Layer Stage 110x120mm, Coaxial Coarse & Fine Focus, N.A.0.65, 0.2W LED	1 pcs / 1 ctn 44*47*24cm, 9/6.8Kgs			
75	A11.0913-B	Biological Microscope	Binocular, WF Plan10X/18mm, Quadruple, Achromatic Objective 4x10x40x100x, Double Layer Mechanical Stage 115x125mm, N.A.1.25, Coaxial Coarse & Fine Focus, 0.2W LED	1 pcs / 1 ctn 44*47*24cm, 9/6.8Kgs			
76	A11.0913-T	Biological Microscope	Trinocular, WF Plan10X/18mm, Quadruple, Achromatic Objective 4x10x40x100x, Double Layer Mechanical Stage 115x125mm, N.A.1.25, Coaxial Coarse & Fine Focus, 0.2W LED	1 pcs / 1 ctn 44*47*24cm, 9/6.8Kgs			
77	A11.2601-M	Student Biological Microscope	Monocular, WF10x/18, Quadruple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage 125x115/75x35mm, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 34*24*42cm, 7/5.5Kgs			
78	A11.2601-B	Student Biological Microscope	Binocular, WF10x/18, Quadruple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage 125x115/75x35mm, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 34*24*42cm, 7/5.5Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
79	A11.2601-T	Student Biological Microscope	Trinocular, WF10x/18, Quaduple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage 125x115/75x35mm, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 34*24*42cm, 7/5.5Kgs			
80	A11.2601-BP	Student Biological Microscope	Binocular, WF10x/18, Quaduple, Plan 4x10x40x100x, Double Layer Mechanical Stage 125x115/75x35mm, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 34*24*42cm, 7/5.5Kgs			
79	A11.2601-T	Student Biological Microscope	Trinocular, WF10x/18, Quaduple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage 125x115/75x35mm, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 34*24*42cm, 7/5.5Kgs			
80	A11.2601-BP	Student Biological Microscope	Binocular, WF10x/18, Quaduple, Plan 4x10x40x100x, Double Layer Mechanical Stage 125x115/75x35mm, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 34*24*42cm, 7/5.5Kgs			
81	A11.2601-TP	Student Biological Microscope	Trinocular, WF10x/18, Quaduple, Plan 4x10x40x100x, Double Layer Mechanical Stage 125x115/75x35mm, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 34*24*42cm, 7/5.5Kgs			
82	A11.0214-A	Biological Microscope	Binocular, WF10x, Quadruple, Achromatic 4x10x40x100x, Mechanical, Abbe NA1.25, Coaxial Coarse & Fine Focusing, Halogen 6V20W	1 pcs / 1 ctn 40*22*31cm, 11/9Kgs			
83	A11.0214-AT	Biological Microscope	Trinocular, WF10x, Quadruple, Achromatic 4x10x40x100x, Mechanical, Abbe NA1.25, Coaxial Coarse & Fine Focusing, Halogen 6V20W	1 pcs / 1 ctn 40*22*31cm, 11/9Kgs			
84	A11.0216-B	Student Biological Microscope	Binocular, WF10x/18, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 115x125mm, Abbe NA1.25, Coaxial Coarse & Fine Focusing, 1W LED Low Voltage	1 pcs / 1 ctn 40*22*31cm, 11/9Kgs			
85	A11.0216-T	Student Biological Microscope	Trinocular, WF10x/18, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 115x125mm, Abbe NA1.25, Coaxial Coarse & Fine Focusing, 1W LED Low Voltage	1 pcs / 1 ctn 40*22*31cm, 11/9Kgs			
86	A12.0909-A	Biological Microscope	Monocular, WF10x/18, Quaduple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage, Size 132*140mm, Abbe N.A.1.25 condenser, Coaxial Coarse & Fine Focusing, 3W LED	1 pcs / 1 ctn 50.5*29.5*57cm, 10/8Kgs			
87	A12.0909-B	Biological Microscope	Binocular, WF10x/18, Quaduple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage, Size 132*140mm, Abbe N.A.1.25 condenser, Coaxial Coarse & Fine Focusing, 3W LED	1 pcs / 1 ctn 50.5*29.5*57cm, 10/8Kgs			
88	A12.0909-C	Biological Microscope	Trinocular, WF10x/18, Quaduple, Achromatic 4x10x40x100x, Double Layer Mechanical Stage, Size 132*140mm, Abbe N.A.1.25 condenser, Coaxial Coarse & Fine Focusing, 3W LED	1 pcs / 1 ctn 50.5*29.5*57cm, 10/8Kgs			


No.	Cata. No.	Product	Specification	Package			Picture
89	A12.0909-S-C	Biological & Stereo Microscope	Trinocular, PL10x/18, Quaduple, Achromatic 4x10x40x, LWD Stereo 2x/0.04, Round Stage, Up/Down 3W LED, Built-in 4 Rechargeable Batteries	1 pcs / 1 ctn 50.5*29.5*57cm , 8/8Kgs			
90	A12.3601-B	Microscope	Binocular, WF10X/18, Quaduple, Achromatic 4x10x40x100x, Coaxial, Mechacnical Stage 140x132mm, LED	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
91	A12.3601-T	Microscope	Trinocular, WF10X/18, Quaduple, Achromatic 4x10x40x100x, Coaxial, Mechacnical Stage 140x132mm, LED, 0.5x C-Mount	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
92	A11.5101	Heating Stage Biological Microscope	Monocular, H10x,H16x, Triple Nosepiece, Achroamtic 4x10x40x, Heated Stage 38C 110x120mm, Coarse Focusing, LED, Brightness Adjustable, Aluminum Case	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
93	A11.5101-B	Heating Stage Biological Microscope	Binocular, H10x,H16x, Triple Nosepiece, Achroamtic 4x10x40x, Heated Stage 10-50°C Adjustable, Size 110x120mm, Coarse Focusing, LED, Brightness Adjustable, Aluminum Case	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
94	A12.0806-PB	Heating Stage Biological Microscope	Binocular, WF10x/20mm, WF16x/15, Quaduple, Plan 4x10x40x100x, Mechanical Stage 160x142/76x52mm, With Heater Temprature 20-50°, Heating Electric Power 12V36W, Abbe N.A.1.25, 3W LED	1 pcs / 1 ctn 42.5*28.5*51cm , 8.5/7Kgs			
95	A12.0806-PT	Heating Stage Biological Microscope	Trinocular, WF10x/20mm, WF16x/15, Quaduple, Plan 4x10x40x100x, Mechanical Stage 160x142/76x52mm, With Heater Temprature 20-50°, Heating Electric Power 12V36W, Abbe N.A.1.25, 3W LED, With 1.0x C Mount	1 pcs / 1 ctn 42.5*28.5*51cm , 8.5/7Kgs			
96	A33.5100	7" LCD Digital Heating Stage Biological Microscope	7" LCD Resolution 800x480, 0.6M 1/3" CMOS, 976x592 Pixels, AV Output, Power by 5V DC USB Adapter, Dual Viewing Head, H10x,H16x, Triple Nosepiece, Achroamtic 4x10x40x, Heated Stage 37-38C, LED, Coaxial	1 pcs / 1 ctn 22*34*45cm, 6/5Kgs			
97	A12.0805-B	Biological Microscope	Binocular, WF10x/20, WF16x/15, Achromatic 4x10x40x100x, Mechanical Stage 160x142/76x52mm, Abbe N.A.1.25, Kohler 3W LED	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
98	A12.0805-T	Biological Microscope	Binocular, WF10x/20, WF16x/15, Achromatic 4x10x40x100x, Mechanical Stage 160x142/76x52mm, Abbe N.A.1.25, Kohler 3W LED	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
Entry Level Laboratory Microscope							
99	A11.1551-IB	Student Biological Microscope	Binocular, WF10x, Quaduple, Infinity Plan 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 142x132mm, 6V20W Halogen	1 pcs / 1 ctn 34*25.5*49cm, 6/5.3Kgs			
100	A11.1551-IT	Student Biological Microscope	Trinocular, WF10x, Quaduple, Infinity Plan 4x10x40x100x, Coarse & Fine Focusing, Double Layer Mechanical Stage 142x132mm, 6V20W Halogen	1 pcs / 1 ctn 34*25.5*49cm, 6/5.3Kgs			
101	A12.0921-B	Biological Microscope	Binocular, PL10x/20mm Diopter Adjustable, High Eye Point, Quaduple, Infinity Plan Achromatic Objective 4x10x40x100x, Double Stage 150x140mm, N.A.1.25, Coaxial Coarse & Fine Focus, 3W LED	1 pcs / 1 ctn 52*49*27cm, 9.5/7.2Kgs			
102	A12.0921-T	Biological Microscope	Trinocular, WF10x/20mm Diopter Adjustable, High Eye Point, Quaduple, Infinity Plan Achromatic Objective 4x10x40x100x, Stage 150x140mm, N.A.1.25, Coaxial Coarse & Fine Focus, 3W LED	1 pcs / 1 ctn 52*49*27cm, 9.5/7.2Kgs			
103	A12.0911-B	Biological Microscope	Binocular, PL10x/20mm, Coded Quintuple, Infinity Plan 4x10x20x40x100x, Double Layer Mechanical Stage 150x162mm /76x 50mm, N.A.1.25 Condenser With DF/PH Slot, 3W LED, LCD Screen On	1 pcs / 1 ctn 57*56*28cm, 10.5/9Kgs			
104	A12.0911-T	Biological Microscope	Trinocular, PL10x/20mm, Coded Quintuple, Infinity Plan 4x10x20x40x100x, Double Layer Mechanical Stage 150x162mm /76x 50mm, N.A.1.25 Condenser With DF/PH Slot, 3W LED, LCD Screen On	1 pcs / 1 ctn 57*56*28cm, 10.5/9Kgs			
105	A12.6603-B	Biological Microscope	Binocular, WF10x/20, Quaduple, Infinity E-Plan 4x10x40x100x, Mechanical Stage 145X140/76x52mm, Coaxial Coarse & Fine Focusing, Abbe Condenser N.A.1.25, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 26*35*46cm, 7/5.5Kgs			
106	A12.6603-T	Biological Microscope	Trinocular, WF10x/20, Quaduple, Infinity E-Plan 4x10x40x100x, Mechanical Stage 145X140/76x52mm, Coaxial Coarse & Fine Focusing, Abbe Condenser N.A.1.25, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 26*35*46cm, 7/5.5Kgs			
107	A12.2601-A	Laboratory Biological Microscope	Binocular, WF10x/18, Quaduple, Achromatic 4x10x40x100x, Mechanical Stage 145x140mm, Working Stage Moving 76x52mm, Abbe N.A.1.25, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 55*38*27cm, 8.5/7Kgs			
108	A12.2601-AT	Laboratory Biological Microscope	Trinocular, WF10x/18, Quaduple, Achromatic 4x10x40x100x, Mechanical Stage 145x140mm, Working Stage Moving 76x52mm, Abbe N.A.1.25, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 55*38*27cm, 8.5/7Kgs			




No.	Cata. No.	Product	Specification	Package			Picture
109	A12.2601-B	Laboratory Biological Microscope	Binocular,WF10x/20, Quarduple, Plan 4x10x40x100x, Mechanical Stage 145x140mm, Working Stage Moving 76x52mm, Abbe N.A.1.25, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 55*38*27cm, 8.5/7Kgs			
110	A12.2601-BT	Laboratory Biological Microscope	Trinocular,WF10x/20, Quarduple, Plan 4x10x40x100x, Mechanical Stage 145x140mm, Working Stage Moving 76x52mm, Abbe N.A.1.25, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 55*38*27cm, 8.5/7Kgs			
111	A12.2601-C	Laboratory Biological Microscope	Binocular,WF10x/20, Quarduple, Infinity E-Plan 4x10x40x100x, Mechanical Stage 145x140mm, Working Stage Moving 76x52mm, Abbe N.A.1.25, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 55*38*27cm, 8.5/7Kgs			
112	A12.2601-CT	Laboratory Biological Microscope	Trinocular,WF10x/20, Quarduple, Infinity E-Plan 4x10x40x100x, Double Layer Mechanical Stage 145x140mm,Working Stage Moving 76x52mm, Abbe N.A.1.25, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 55*38*27cm, 8.5/7Kgs			
113	A12.2601-D	Laboratory Biological Microscope	Binocular,WF10x/20, Quarduple, Infinity Plan 4x10x40x100x, Mechanical Stage 145x140mm, Working Stage Moving 76x52mm, Abbe N.A.1.25, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 55*38*27cm, 8.5/7Kgs			
114	A12.2601-DT	Laboratory Biological Microscope	Trinocular,WF10x/20, Quarduple, Infinity Plan 4x10x40x100x, Mechanical Stage 145x140mm, Working Stage Moving 76x52mm, Abbe N.A.1.25, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 55*38*27cm, 8.5/7Kgs			
115	A12.2602-B	Laboratory Biological Microscope	Binocular, WF10x/20, Quinduple, Infinity Plan 4x10x40x100x, Mechanical Stage 175x145mm, Abbe N.A. 1.25, Halogen 12V20W	1 pcs / 1 ctn 54.2*26.5*35.5cm, 11/0Kgs			
116	A12.2602-T	Laboratory Biological Microscope	Trinocular, WF10x/20, Quinduple, Infinity Plan 4x10x40x100x, Mechanical Stage 175x145mm, Abbe N.A. 1.25, Halogen 12V20W	1 pcs / 1 ctn 54.2*26.5*35.5cm, 11/0Kgs			
Middle Level Laboratory Microscope							
117	A12.0205-B	Laboratory Biological Microscope	Binocular,WF10x/22mm,Quarduple,Infinity Plan 4x10x40x100x,Double layer Mechanical 180*150mm/75*50mm ,Abbe NA1.25,Coaxial Coarse & Fine Focusing,3W LED	1 pcs / 1 ctn 36*27*50cm, 12/0Kgs			
118	A12.0205-T	Laboratory Biological Microscope	Trinocular,WF10x/22mm,Quarduple,Infinity Plan 4x10x40x100x,Double layer Mechanical 180*150mm/75*50mm, Abbe NA1.25,Coaxial Coarse & Fine Focusing,3W LED	1 pcs / 1 ctn 36*27*50cm, 12/0Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
119	A12.0810-B	Laboratory Microscope	Binocular, WF10x/20, WF16x/16, Quarduple, Infinity Plan HC ICCF 4x10x40x100x, Mechanical Stage 188x145mm/76x54mm, Abbe N.A.1.25, Kohler 3W LED	1 pcs / 1 ctn 74*46*31cm, 17/13Kgs			
120	A12.0810-T	Laboratory Microscope	Trinocular, WF10x/20, WF16x/16, Quarduple, Infinity Plan HC ICCF 4x10x40x100x, Mechanical Stage 188x145mm/76x54mm, Abbe N.A.1.25, Kohler 3W LED	1 pcs / 1 ctn 74*46*31cm, 17/13Kgs			
121	A12.0701-E	Biological Microscope	Seidentopf Binocular,30° Inclined,360° Rotatable,Interpupillary 48-75mm;Eyepiece WF10x/18;Achromatic Objective 4x10x40x100x;Quadruple Nosepiece;Coaxial Coarse And Fine Adjustment,Fine Division	1 pcs / 1 ctn 33*27*46cm, 7/6Kgs			
122	A12.0701-F	Biological Microscope	Seidentopf Trinocular,30° Inclined,360° Rotatable,Interpupillary 48-75mm;Eyepiece WF10x/18;Achromatic Objective 4x10x40x100x;Quadruple Nosepiece;Coaxial Coarse And Fine Adjustment,Fine Division	1 pcs / 1 ctn 34*29*66cm, 7/6Kgs			
123	A12.1010-A	Laboratory Microscope	Binocular, WF10x/20, Quadruple, Infinity Plan 4x10x40x100x, Mechanical Stage 230x150mm, Abbe 1.25, Coaxial Coarse&Fine, 12V20W	1 pcs / 1 ctn 27*38*42cm, 8.3/8Kgs			
124	A12.1010-B	Laboratory Microscope	Trinocular, WF10x/20, Quadruple, Infinity Plan 4x10x40x100x, Mechanical Stage 230x150mm, Abbe 1.25, Coaxial Coarse&Fine, 12V20W	1 pcs / 1 ctn 27*38*42cm, 8.3/8Kgs			
125	A12.1150-B	Laboratory Biological Microscope	Binocular, WF10x/22, Quintuple, Infinity Semi-Plan 4x10x40x100x, Double Layer Stage 180x160mm /50x75mm, Kohler Illumination 12V30W Halogen	1 pcs / 1 ctn 53*33*39cm, 10/9Kgs			
126	A12.1150-T	Laboratory Biological Microscope	Trinocular, WF10x/22, Quintuple, Infinity Semi-Plan 4x10x40x100x, Double Layer Stage 180x160mm /50x75mm, Kohler Illumination 12V30W Halogen	1 pcs / 1 ctn 53*33*39cm, 10/9Kgs			
127	A12.1502-AB	Laboratory Microscope	Binocular,WF10x/20, Quarduple, Achromatic HC 4x,10x,40x100x, Abbe N.A.1.25, Coaxial, Double Layer Mechanical Graphite 143*132mm, Moving 75x50mm, 3W LED	1 pcs / 1 ctn 53*33*39cm, 10/8Kgs			
128	A12.1502-AT	Laboratory Microscope	Trinocular,WF10x/20, Quarduple, Achromatic HC 4x,10x,40x100x, Abbe N.A.1.25, Coaxial, Double Layer Mechanical Graphite 143*132mm, Moving 75x50mm, 3W LED, 0.5x C Mount	1 pcs / 1 ctn 53*33*39cm, 10/8Kgs			






No.	Cata. No.	Product	Specification	Package			Picture
129	A12.1502-B	Laboratory Microscope	Binocular, WF10x/22, Quaduple, Infinity Plan HC 4x,10x,40x100x, Abbe N.A.1.25, Coaxial, Double Layer Mechanical Graphite 216x150mm, Moving 75x55mm, 3W LED Kohler	1 pcs / 1 ctn 53*33*39cm, 10/8Kgs			
130	A12.1502-T	Laboratory Microscope	Trinocular, WF10x/22, Quaduple, Infinity Plan HC 4x,10x,40x100x, Abbe N.A.1.25, Coaxial, Double Layer Mechanical Graphite 216x150mm, Moving 75x55mm, 3W LED Kohler, 0.5x C Mount	1 pcs / 1 ctn 53*33*39cm, 10/8Kgs			
131	A12.1503-A	Biological Microscope	Binocular, WF10/18, Quaduple, 4x10x40x100x, Abbe 1.25, Coaxial Coarse & Fine Focusing, Double Layer Mechanical Stage 142x132mm/ 75x50mm, 3W LED	1 pcs / 1 ctn 36*25*46cm, 7/0Kgs			
132	A12.1503-B	Biological Microscope	Trinocular, WF10/18, Quaduple, 4x10x40x100x, Abbe 1.25, Coaxial Coarse & Fine Focusing, Double Layer Mechanical Stage 142x132mm/ 75x50mm, 3W LED, 0.5x C-Mount	1 pcs / 1 ctn 36*25*46cm, 7/0Kgs			
133	A12.2701-M	Laboratory Microscope	Monocular, WF10x/18, Quaduple, Infinity Ach.4x10x40x100x, Coaxial Focusing, Mechanical Stage, Abbe NA1.25, Kohler 6V20W Halogen or LED	1 pcs / 1 ctn 48*47*28cm, 9/0Kgs			
134	A12.2701-A	Laboratory Microscope	Binocular, WF10x/18, Quaduple, Infinity Ach.4x10x40x100x, Coaxial Focusing, Mechanical Stage, Abbe NA1.25, Kohler 6V20W Halogen or LED	1 pcs / 1 ctn 29*47*48cm, 9/0Kgs			
135	A12.2701-B	Laboratory Microscope	Trinocular, WF10x/18, Quaduple, Infinity Ach.4x10x40x100x, Coaxial Focusing, Mechanical Stage, Abbe NA1.25, Kohler 6V20W Halogen or LED	1 pcs / 1 ctn 48*47*28cm, 9/0Kgs			
136	A12.2702-A	Laboratory Microscope	Binocular, WF10x/20, Quaduple, Infinity Plan 4x10x40x100x, Coaxial Focusing, Mechanical Stage 156x138mm, Abbe NA1.25, Kohler 6V20W Halogen or LED	1 pcs / 1 ctn 48*48*33cm, 12/10Kgs			
137	A12.2702-B	Laboratory Microscope	Trinocular, WF10x/20, Quaduple, Infinity Plan 4x10x40x100x, Coaxial Focusing, Mechanical Stage 156x138mm, Abbe NA1.25, Kohler 6V20W Halogen or LED, 1x C Mount	1 pcs / 1 ctn 48*48*33cm, 12/10Kgs			
High Level Laboratory Microscope							
138	A12.2603-B	Laboratory Biological Microscope	Binocular, WF10x/22 Dia.30mm, Quinduple, Infinity Plan 4x10x20x40x100x, Rackless Mechanical Stage 182x140mm, Travel Range 77x52mm, Swing Out Condenser, Halogen 6V30W	1 pcs / 1 ctn 52*47*28.5cm, 11/10Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
139	A12.2603-T	Laboratory Biological Microscope	Trinocular, WF10x/22 Dia.30mm, Quinduple, Infinity Plan 4x10x20x40x100x, Rackless Mechanical Stage 182x140mm, Travel Range 77x52mm, Swing Out Condenser, Halogen 6V30W	1 pcs / 1 ctn 52*47*28.5cm, 11/10Kgs			
140	A12.0204-A	Laboratory Biological Microscope	Trinocular, WF10x/22, Quadruple, Infinity Plan 4x10x40x100x, Double Layer Mechanical Stage 210*140mm/75*50mm, Abbe NA1.25, Coaxial Coarse & Fine Focusing, Halogen 6V30W Kohler	1 pcs / 1 ctn 37*33*55cm, 11/10Kgs			
141	A12.0204-B	Laboratory Biological Microscope	Trinocular, WF10x/22, Quintuple, Infinity Plan 4x10x20x40x100x, Double Layer Mechanical Stage 210*140mm/75*50mm, Abbe NA1.25, Coaxial Coarse & Fine Focusing, Halogen 6V30W Kohler	1 pcs / 1 ctn 37*33*55cm, 11/10Kgs			
142	A12.0204-DIC	DIC Laboratory Biological Microscope	Trinocular, WF10x/22, Quintuple, Infinity Plan LWD 40x, 100x, DIC 5x, 10x, 20x, Mechanical Stage 140*210/75x50mm, Condenser Combine Abbe, Polarizer, DIC All Together, Halogen 6V30W Kohler	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
143	A12.0908-A	Laboratory Microscope	Binocular, PL10x/22mm Dia.30mm, Quintuple, Infinity Plan 4x10x20x40x100x, Double Layer 175x145/76x50mm, Swing-out Condenser N.A.1.2/0.22, Halogen 6V30W	1 pcs / 1 ctn 43.9*41.6*55cm , 15.5/11Kgs			
144	A12.0908-B	Laboratory Microscope	Trinocular, PL10x/22mm Dia.30mm, Quintuple, Infinity Plan 4x10x20x40x100x, Double Layer 175x145/76x50mm, Swing-out Condenser N.A.1.2/0.22, Halogen 6V30W	1 pcs / 1 ctn 43.9*41.6*55cm , 15.5/11Kgs			
145	A12.0910	Research Scientific Laboratory Microscope, Semi-APO, Upgrade BF+DF+PL+PH+FL+DIC	Trinocular, PL10x/25, Sextuple BF Nosepiece With Slot For DIC/PL, Infinity Plan Semi-APO Fluorescent 4x, 10x, 20x, 40x, 100x, Double Mechanical Stage 187*166mm, Moving Range 80*55mm, Swing-Out	1 pcs / 1 ctn 89*51.3*51cm, 28/23Kgs			
146	A12.0910-LED	Research Scientific Laboratory Microscope, Semi-APO, Upgrade BF+DF+PL+PH+FL+DIC	Trinocular, PL10x/23, Sextuple BF Nosepiece With Slot For DIC/PL and intensity memory, Infinity Plan objectives 4x, 10x, 20x, 40x, 100x, Double Mechanical Stage 187*166mm, Moving Range 80*55mm, Non Swing-Out	1 pcs / 1 ctn 89*51.3*51cm, 28/23Kgs			
147	A12.0933	Biological Microscope	Trinocular, PL10x/22mm, Coded Quintuple, Infinity Plan 4x10x40x100x, Double Layer Mechanical Stage 210x171mm /78x51mm, N.A.1.25 Condenser With DF/PH Slot, 3W LED, LCD Screen, Color Temperature Adjustable 2800K-7000K, Moving Range 80*55mm	1 pcs / 1 ctn 57*56*28cm, 10.5/9Kgs			
148	A12.1018	Laboratory Microscope	Trinocular, EW10x/22mm, Quintuple, Infinity Plan Objective 4x10x40x100x, Mechanical Stage 185x142/ 55x75mm, Swingout NA0.9/0.25, Coaxial Coarse & Fine, External Halogen 6V30W, Kohler Illumination Aspherical Collector, Eyepiece	1 pcs / 1 ctn 43*39*61cm, 14/0Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
149	A12.1025-A	Laboratory Microscope	Titular Binocular, EW10x, Quintuple, Infinity Plan Objective 4x10x40x100x, Mechanical Stage, Swingout NA0.9-0.25, Coaxial Coarse & Fine, LED Kohler Illumination	1 pcs / 1 ctn 43*39*61cm, 14/0Kgs			
150	A12.1025-B	Laboratory Microscope	Titular Trinocular, EW10x, Quintuple, Infinity Plan Objective 4x10x40x100x, Mechanical Stage, Swingout NA0.9-0.25, Coaxial Coarse & Fine, LED Kohler Illumination	1 pcs / 1 ctn 43*39*61cm, 14/0Kgs			
151	A12.1030-B	Laboratory Microscope, Coded Nosepiece & LCD	Binocular, WF10x/20, Coded Quadruple, Infinity Plan NIS45 4x10x40x100x, Rackless Stage 180x130/ 74x30mm, N.A.1.25 Condenser With Objective Position Scales, Coaxial Focusing, 3W LED	1 pcs / 1 ctn 47*32*59cm, 12.1/11Kgs			 HOT
152	A12.1030-T	Laboratory Microscope, Coded Nosepiece & LCD	Trinocular, WF10x/20, Coded Quadruple, Infinity Plan NIS45 4x10x40x100x, Rackless Stage 180x130/ 74x30mm, N.A.1.25 Condenser With Objective Position Scales, Coaxial Focusing, 3W LED	1 pcs / 1 ctn 47*32*59cm, 12.1/11Kgs			 HOT
153	A12.1038-B	Laboratory Microscope, LCD Screen	Binocular, WF10x/20, Coded Quintuple, Infinity Semi-Plan 4x10x40x100x, Mechanical Stage 150x140/75x52mm, Abbe Condenser N.A.1.25, Coaxial Coarse & Fine Focusing, 3W LED, Green Filter, LCD Screen Base	1 pcs / 1 ctn 47*32*59cm, 12/11Kgs			
154	A12.1038-T	Laboratory Microscope, LCD Screen	Trinocular, WF10x/20, Coded Quintuple, Infinity Semi-Plan 4x10x40x100x, Mechanical Stage 150x140/75x52mm, Abbe Condenser N.A.1.25, Coaxial Coarse & Fine Focusing, 3W LED, Green Filter, LCD Screen Base	1 pcs / 1 ctn 47*32*59cm, 12/11Kgs			
155	A12.1062-B	Laboratory Microscope	Binocular, WF10x/22, Coding Quintuple, Infinity Plan NIS60 4x10x40x100x, Rackless Stage 230x150/ 78x54mm, N.A.1.25 Condenser With Slot, Coaxial, 3W LED, LCD Screen Base, Brightness Auto	1 pcs / 1 ctn 47*32*59cm, 12.1/11Kgs			 HOT
156	A12.1062-T	Laboratory Microscope	Trinocular, WF10x/22, Coding Quintuple, Infinity Plan NIS60 4x10x40x100x, Rackless Stage 230x150/ 78x54mm, N.A.1.25 Condenser With Slot, Coaxial, 3W LED, LCD Screen Base, Brightness Auto	1 pcs / 1 ctn 47*32*59cm, 12.1/11Kgs			 HOT
157	A12.1091-H	Research Scientific Laboratory Microscope, Manual	Trinocular, EW10x/25mm, Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler Illumination, 12V100W Halogen, ECO Function	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			
158	A12.1091-L	Research Scientific Laboratory Microscope, Manual	Trinocular, EW10x/25mm, Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler Illumination, 3W SLED, ECO Function	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			











No.	Cata. No.	Product	Specification	Package			Picture
159	A12.1092-H	Research Scientific Laboratory Microscope, Manual	Trinocular, EW10x/25mm, Coded Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler Illumination, 12V100W Halogen , ECO	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			
160	A12.1092-L	Research Scientific Laboratory Microscope, Manual	Trinocular, EW10x/25mm, Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler Illumination, 3W SLED, ECO Function,	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			
161	A12.1093-H	Research Scientific Laboratory Microscope, Semi Motorized	Tilting Trinocular 0-35°, EW10x/25mm, Coded Motorized Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Motorized Swing-Out Condenser NA0.9/0.25, Coaxial, Transmit	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			
162	A12.1093-L	Research Scientific Laboratory Microscope, Semi Motorized	Tilting Trinocular 0-35°, EW10x/25mm, Motorized Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Motorized Swing-Out Condenser NA0.9/0.25, Coaxial, Transmit 3W	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			
163	A12.1095	Research Scientific Laboratory Microscope, Motorized Nosepiece + XYZ Stage	Trinocular, EW10x/25mm, Motorized Coded Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Motorized XYZ Control Stage 190x152/125x75mm, With Joy Stick Controller, Swing-Out Condenser NA0.9/0.25, Coaxial Auto Focus, Transmit 3W SLED Kohler Illumination, Touch Screen, Motor Objective Switch, Brightness Memory, ECO Setup, With Software	1 pcs / 1 ctn 79*58*52cm, 28/26Kgs			

Nikon, Olympus

164	A12.0705-A	Nikon E100 Biological Microscope	CFI Optical System;E2-TB Binocular,30° Inclined;Eyepiece E1-CF1 10x/18;Quadruple Nosepiece;CFI BE Plan Achromatic Objective 4x10x40x100x; Abbe Condenser NA1.25;Coaxial Coarse And Fine	1 pcs / 1 ctn 51*30*33cm, 8/0Kgs			
165	A12.0705-B	Nikon E100 Biological Microscope	CFI Optical System;E2-TB Trinocular,30° Inclined;Eyepiece E1-CF1 10x/18;Quadruple Nosepiece;CFI BE Plan Achromatic Objective 4x10x40x100x; Abbe Condenser NA1.25;Coaxial Coarse And Fine	1 pcs / 1 ctn 51*30*33cm, 8/0Kgs			
166	A12.0706-A	Biological Microscope	E200 Binocular, Eyepiece CFI 10x/20;Quadruple Nosepiece;CFI E Plan Achromatic 4x10x40x100x;E2 Abbe Condenser NA1.25;Coaxial Coarse And Fine Adjustment,Rectangle Stage 216x150mm/78mmx54mm;	1 pcs / 1 ctn 33*37*55cm, 15/0Kgs			
167	A12.0706-B	Biological Microscope	E200 Trinocular, Eyepiece CFI 10x/20;Quadruple Nosepiece;CFI E Plan Achromatic 4x10x40x100x;E2 Abbe Condenser NA1.25;Coaxial Coarse And Fine Adjustment,Rectangle Stage 216x150mm/78mmx54mm;	1 pcs / 1 ctn 33*37*55cm, 15/0Kgs			
168	A12.0736-B	Olympus Biological Microscope, CX23	Binocular, WF10x/20, Quadruple, Infinity Plan 4x10x40x100x, Rackless Mechanical Stage 174x89mm, Moving 76x30mm, Abbe Condenser, LED 0.5W Light Source	1 pcs / 1 ctn 30*35*50cm, 6.5/6Kgs			


No.	Cata. No.	Product	Specification	Package			Picture
169	A12.0736-T	Olympus Biological Microscope, CX23	Trinocular, WF10x/20, Quaduple, Infinity Plan 4x10x40x100x, Rackless Mechanical Stage 174x89mm, Moving 76x30mm, Abbe Condenser, LED 0.5W Light Source, 0.5x C-Mount	1 pcs / 1 ctn 30*35*50cm, 6.5/6Kgs			
170	A12.0737-31B	OLYMPUS Microscope	Binocular, WF10x, Quaduple, Plan 4x10x40x100x, Abbe NA1.25, Mechanical Stage, 6V30W Halogen	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
171	A12.0737-31T	OLYMPUS Microscope	Trinocular, WF10x, Quaduple, Plan 4x10x40x100x, Abbe NA1.25, Mechanical Stage, 6V30W Halogen	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
172	A12.0738	Olympus Biological Microscope, CX33	Trinocular, WF10x/20mm, Quaduple, Plan Achromat (UIS2) 4x10x40x100x, Abbe Condenser, Wire Movement Working Stage 211x154mm, Moving 76x52mm, 2.4W LED BF,DF	1 pcs / 1 ctn 30*35*50cm, 6.5/6Kgs			
173	A12.0739-B	Olympus Biological Microscope, CX43	Binocular, WF10x/20mm, Quintuple, Plan Achromat (UIS2) 4x10x40x100x, Universal Condenser With 7 Turret Positions, Wire Movement Working Stage 211x154mm, Moving 76x52mm, 2.4W LED	1 pcs / 1 ctn 30*35*50cm, 6.5/6Kgs			
174	A12.0739-T	Olympus Biological Microscope, CX43	Trinocular, WF10x/20mm, Quintuple, Plan Achromat (UIS2) 4x10x40x100x, Universal Condenser With 7 Turret Positions, Wire Movement Working Stage 211x154mm, Moving 76x52mm, 2.4W LED	1 pcs / 1 ctn 30*35*50cm, 6.5/6Kgs			
175	A12.0740-53	Olympus Fluorescent Microscope	BX53F21BX53 Frame U-LHLEDC1001High Power LED Lamp Housing U-TR30-21Trinocular tube WHN10X-H1Widefield eyepiece 10X, focusable	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
176	A12.0740-43	Olympus Fluorescent Microscope	BX43F1Microscope frame for transmitted microscopy, including hexagonal driver ABxxxx 显微镜主机 U-TR30-21Trinocular tube 三目观察筒	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
177	A12.0740-CKX53SF	Olympus Fluorescent Microscope	CKX53SF1"CKX53 frame with binocular tube and eyepieces (2pcs), including L-wrench GA4720, clamp for camera cable DO123600 (2pcs) and AB4x8SA screw (2pcs) CKX53机架, 带三目观察筒"	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
Inverted Biological Microscope							
178	A14.0201	Inverted Biological Phase Contrast Microscope	Trinocular, WF10x/20, Quaduple, LWD Plan 10x25x40x, LWD Plan Phase Contrast 10x, Impellent Phase Contrast Set, Double Layer Mechanical Stage 200*152/75*30mm, LWD Condenser 30mm, Halogen 6V20W	1 pcs / 1 ctn 39*31*55cm, 12/10Kgs			











No.	Cata. No.	Product	Specification	Package			Picture
179	A14.0205	Inverted Biological Phase Contrast Microscope, DIC, PMC	Trinocular, WF10x/22mm, Centering Telescope 10x, Quaduple, Bright Field LWD Infinity Plan 10x20x40x, LWD Infinity Plan Phase Contrast 10x, Working Stage 224x205mm, Moving Stage 112x77mm, LWD Condenser	1 pcs / 1 ctn 64*43*60cm, 16/14Kgs			
180	A14.0900-A	Inverted Microscope, BF+PH, Kohler Illumination	Binocular, PL10x/22mm, Quintuple, LWD Infinity Plan Achromatic 4x10x40x, LWD Infinity Plan Contrast 20x, Fixed Stage 250x160mm, Glass Plate, Metal Plate, Petri Dish Holder 35mm, LWD Condenser N.A. 0.3, WD	1 pcs / 1 ctn 82.5*62*39.2cm , 14/11Kgs			
181	A14.0901-A	Inverted Microscope, BF+PH, Critical Illumination	Binocular, PL10x/22mm, Quintuple, LWD Infinity Plan Achromatic 4x10x40x, LWD Infinity Plan Contrast 20x, Fixed Stage 250x160mm, Glass Plate, Metal Plate, Petri Dish Holder 35mm, LWD Fixed Condenser N.A. 0.3,	1 pcs / 1 ctn 82.5*62*39.2cm , 13.8/11Kgs			
182	A14.0901-B	Inverted Microscope, BF+PH, Critical Illumination	Trinocular, PL10x/22mm, Quintuple, LWD Infinity Plan Achromatic 4x10x40x, LWD Infinity Plan Contrast 20x, Fixed Stage 250x160mm, Glass Plate, Metal Plate, Petri Dish Holder 35mm, LWD Fixed Condenser N.A. 0.3,	1 pcs / 1 ctn 82.5*62*39.2cm , 13.8/11Kgs			
183	A14.0912	Inverted Microscope, BF+PH, ECO	Trinocular, PL10x/22mm, Quintuple, Infinity LWD Plan 4x10x40x, LWD Postive Phase Contrast 20x, PH Slide 4x10x+20x40x, Centering Telescope, Fixed Stage 215*250mm, Glass Plate, Metal Plate, Petri Dish Holder 35mm,	1 pcs / 1 ctn 68*67*47cm, 21/19Kgs			
184	A14.1021	Inverted Microscope	Trinocular, EW10x/22, Quintuple, Infinity Plan LWD 4x40x, Infinity Plan Phase Contrast 10x20x, ELWD Condenser N.A.0.3, W.D.72mm, Max W.D.150mm, Fix Stage 160x250mm, Auxiliary Stage 70x180mm Phase	1 pcs / 1 ctn 46*41*67cm, 15/11Kgs			
185	A14.1063	Inverted Biological Microscope, Semi-APO, PH, DIC	Binocular (Photo Port On Side), WF10x/22, Quintuple, NIS Infinity Plan 4x, Infinity Plan Phase Contrast 10x20x40x, Plain Stage 170x250mm, Attachable Moving Ruler 128x80mm, Glass Plate, Auxilliary Plate, Universal	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
186	A14.1064	Inverted Biological Microscope, Semi-APO, PH, DIC	Binocular (Photo Port On Side), WF10x/22, Coded Quintuple Nosepiece, NIS Infinity Plan 4x, Infinity Plan Phase Contrast 10x20x40x, Plain Stage 170x250mm, Attachable Moving Ruler 128x80mm, Glass Plate, Auxilliary	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
187	A14.1065	Inverted Biological Microscope, Semi-APO, PH, DIC, LCD	Tilting Binocular (Photo Port On Side), WF10x/22, Coded Quintuple Nosepiece, NIS Infinity Plan 4x, Infinity Plan Phase Contrast 10x20x40x, Plain Stage 170x250mm, Attachable Moving Ruler 128x80mm, Glass Plate, Auxilliary	1 pcs / 2 ctn 69*39*65cm, 20/18Kgs			
188	A14.1096	Inverted Laboratory Microscope, Semi-Auto, Semi-APO, BF+PH	Trinocular, Built-in Bertrand Lens, SW10x/22mm, Sextuple Nosepiece With DIC Slot, Infintiy Plan Sem-APO PH10x20x40x, Triple Layer Working Stage Moving Range 130x85mm, With 3 Holder, LWD Condensor N.A.0.55	1 pcs / 1 ctn 120*80*70cm, 70/60Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
189	A14.1097	Inverted Biological Microscope, Semi-Auto, Semi-APO, BF+PH	Trinocular, Built-in Bertrand Lens, SW10x/22mm, LCD Touch Screen + Coded Sextuple Nosepiece With DIC Slot, Infinity Plan Semi-APO PH10x20x40x, Triple Layer Working Stage Moving Range 130x85mm, With	1 pcs / 1 ctn 120*80*70cm, 70/60Kgs			
190	A14.1501	Student Mini Cordless Inverted Microscope	Monocular, WF10x/14mm, WF16x, Triple, 4x10x20x, Coarse Focusing, Singler Layer Working Stage 105x118mm, LWD Condenser N.A.0.19, 0.5W LED, 3 AA 1.5V Batteries, Mobile Phone Holder	1 pcs / 1 ctn 26*21*37cm, 1.5/1.1Kgs			
191	A14.2603	Inverted Microscope, Transmit Light, Semi-APO	Trinocular, WF10x/22, Quintuple, Infinity Plan LWD Semi-APO FL10x20x40x, Infinity Phase Contrast LWD Semi-APO FL20x Objective, 20x/40x Phase Contrast Annular Spot, Halogen 6V30W	1 pcs / 1 ctn 66*60*33cm, 17/13.6Kgs			
192	A14.2603-TR	Inverted Microscope, Transmit & Relfect Light	Trinocular, WF10x/22, Quintuple, Infinity Plan LWD Semi-APO FL10x20x40x, Infinity Phase Contrast LWD Semi-APO FL20x, 20x/40x Phase Contrast Annular Spot, Incl. Attachable Stage For 96 Holes Plate, Transmit 5W	1 pcs / 1 ctn 66*60*33cm, 17/13.6Kgs			
193	A14.2605	Inverted Biological Phase Contrast Microscope	Trinocular + Dual Photo Output, WF10x/22, Quintuple, Infinity Plan Semi-APO LWD Fluorescent FL4x, Infinity Semi-APO Phase Contrast LWD Fluorescent FL10x20x40x, Turret Phase Contrast Condenser, Incl.	1 pcs / 1 ctn 66*60*33cm, 18/13.5Kgs			
194	A14.2701	Inverted Microscope	Binocular Head, Side Camera Port, Eyepiece WF10x/20, Quintuple, ing Telescope, Infinity Plan 10x20x40 & Infinity Phase Contract Plan 20x, ABBE N.A. 0.3 (W.D.72mm), Mechanical Stage 240x260/135x85mm, Coaxial Corse & Fine Focusing, 6V30W Halogen Illumination, PH Slider 10x20x40x, 0.5x C-Mount	1 pcs / 1 ctn 69*68*42cm, 32/0Kgs			
195	A14.2702	Inverted Microscope	UCIS Infinity Independent Achromatic Optical System, Trinocular, WF10X/20, Centering Telescope, Infinity Long Work Distance Plan Achromatic 10x, 20x, 40x, Infinity Plan Phase-Contrast Positive Objective	1 pcs / 1 ctn 30*53*66cm, 16.5/14Kgs			
Upright Metallurgical Microscope							
196	A13.2501-B	Portable Metallurgical Microscope	WF10x/18mm, LWD Plan Coverless Objective 10x, 40x, Rechargeable Batteries, Aluminum Case	1 pcs / 1 ctn 32*23*11cm, 3/2.5Kgs			
197	A13.0201-A2	Inverted Metallurgical Microscope, Reflect, Polarizing Optional	Trinocular, WF10x/18mm, Quadruple, LWD Plan Coverless 10x20x40x100x, Mechanical Stage 180x150mm/15x15mm, Coaxial Focus, 6V20W Halogen	1 pcs / 1 ctn 30*40*47cm, 11/10Kgs			
198	A13.0202-A	Metallurgical Microscope, Reflect, Polarizing	Trinocular, WF10x/18, Quadruple, LWD Plan Coverless 5x10x40x60x, Double Layer Mechanical Stage 185x140/75x40mm, Reflect Illumination 6V20W Halogen	1 pcs / 1 ctn 45*30*35cm, 9.2/0Kgs			





No.	Cata. No.	Product	Specification	Package			Picture
199	A13.0204-A	Metallurgical Microscope, Transmit & Reflect, Polarizing	Trinocular, WF10x/18, Quadruple, LWD Plan Coverless 5x10x40x60x, Double Layer Mechanical Stage 185x140/75x40mm, Transmit & Reflect Kohler Illumination 6V20W Halogen	1 pcs / 1 ctn 45*30*35cm, 10/9Kgs			
200	A13.0207	Upright Metallurgical Microscope, Reflect, BF/DF, PL, DIC	Trinocular, WF10x/22, Quintuple, LWD Infinity Plan PLL5x10x20x50x80x, Mechanical Stage 280x270/ 204x204mm, Epi Illumination 6V30W Halogen, Push/Pull Polaring Slide	1 pcs / 1 ctn 59*50*52cm, 22/20Kgs			
201	A13.0207-BD	Upright Metallurgical Microscope, Reflect, BF/DF, PL, DIC	Trinocular, WF10x/22, Quintuple, LWD Infinity Plan PLL BD 5x10x20x50x80x, Mechanical Stage 280x270/ 204x204mm, Epi Illumination 12V/50W Halogen, Push/Pull Polaring Slide	1 pcs / 1 ctn 59*50*52cm, 22/20Kgs			
202	A13.0207-DIC	Upright Metallurgical Microscope, Reflect, BF/DF, PL, DIC	Trinocular, WF10x/22, Quintuple, LWD Infinity Plan PLL BD 50x80x, LWD Infinity DIC Plan LMPlan 5x10x20x, DIC Slide For 5x10x20x, Mechanical Stage 280x270/ 204x204mm, Epi Illumination 12V/50W Halogen, Polaring Unit	1 pcs / 1 ctn 59*50*52cm, 22/20Kgs			
203	A13.0207-BD DIC	Upright Metallurgical Microscope, Reflect, BF/DF, PL, DIC	Trinocular, WF10x/22, Quintuple, LWD Infinity DIC BD Plan LMPlan 5x10x20x, LWD Infinity Plan PLL50x, DIC Slide For 5x10x20x, Mechanical Stage 280x270/ 204x204mm, Epi Illumination 12V/50W Halogen, Push/Pull Polaring Slide	1 pcs / 1 ctn 59*50*52cm, 22/20Kgs			
204	A13.0210	Metallurgical Microscope, BF/DF, DIC	Binocular + Side Photo Tube, WF10x/22, Quintuple, LWD Infinity Plan 10x20x50x100x, Double Layer Mechanical Stage242x200/30*30mm, Rotatable Stage Max Measurement Dia.130mm, Transmit 6V30W Halogen	1 pcs / 1 ctn 50*35*62cm, 17/14Kgs			
205	A13.0210-BD	Metallurgical Microscope, BF/DF, DIC	Binocular + Side Photo Tube, WF10x/22, Quintuple, LWD Infinity Plan BF/DF 5x10x20x50x, Double Layer Mechanical Stage242x200/30*30mm, Rotatable Stage Max Measurement Dia.130mm, Transmit 12V50W Halogen	1 pcs / 1 ctn 50*35*62cm, 17/14Kgs			
206	A13.0210-DIC	Metallurgical Microscope, BF/DF, DIC	Binocular + Side Photo Tube, WF10x/22, Quintuple, LWD Infinity Plan DIC 5x10x20x, Plan 50x, DIC Slide For 5x10x20x, Double Layer Mechanical Stage242x200/30*30mm, Rotatable Stage Max Measurement Dia.130mm, Transmit 12V50W Halogen	1 pcs / 1 ctn 50*35*62cm, 17/14Kgs			
207	A13.0210-BDDIC	Metallurgical Microscope, BF/DF, DIC	Binocular + Side Photo Tube, WF10x/22, Quintuple, LWD Infinity Plan BF/DF DIC 5x10x20x, BD Plan 50x, DIC Slide For 5x10x20x, Double Layer Mechanical Stage242x200/30*30mm, Rotatable Stage Max Measurement Dia.130mm, Transmit 12V50W Halogen	1 pcs / 1 ctn 50*35*62cm, 17/14Kgs			
208	A13.0211	Upright Metallurgical Microscope, Transmit & Reflect, BF/DF, DIC, Polarizing	Trinocular, WF10x/22, Quintuple, LWD Infinity Plan PLL5x10x40x60x, Double Layer Stage 210x140mm, Moving 63x50mm, Reflect Illumination 6V30W Halogen + Polarizing, Transmit Illumination 6V30W Halogen	1 pcs / 1 ctn 37*33*69cm, 17/15Kgs			











No.	Cata. No.	Product	Specification	Package			Picture
209	A13.0211-BD	Upright Metallurgical Microscope, Transmit & Reflect, BF/DF, DIC, Polarizing	Trinocular, WF10x/22, Quintuple, LWD Infinity Plan BF/DF PLL5x10x20x40x, Double Layer Stage 210x140mm, Moving 63x50mm, Transmit Illumination 12V30W Halogen + Polarizing, Reflect Illumination 12V50W Halogen	1 pcs / 1 ctn 37*33*69cm, 17/15Kgs			
210	A13.0211-DIC	Upright Metallurgical Microscope, Transmit & Reflect, BF/DF, DIC, Polarizing	Trinocular, WF10x/22, Quintuple, LWD Infinity Plan Coverless LMPlan DIC 5x10x20x, Plan 40x, 5X/10X/20X DIC Push-Pull Group, Double Layer Stage 210x140mm, Moving 63x50mm, Reflect Illumination 12V50W Halogen +	1 pcs / 1 ctn 37*33*69cm, 17/15Kgs			
211	A13.0211-BDDIC	Upright Metallurgical Microscope, Transmit & Reflect, BF/DF, DIC, Polarizing	Trinocular, WF10x/22, Quintuple, LWD Infinity Plan Coverless BF/DF LMPlan DIC 5x10x20x, BF/DF Plan 40x, 5X/10X/20X DIC Push-Pull Group, Double Layer Stage 210x140mm, Moving 63x50mm, Reflect Illumination 12V50W Halogen +	1 pcs / 1 ctn 37*33*69cm, 17/15Kgs			
212	A13.0212	Metallurgical Microscope, BF/DF, DIC	Trinocular, WF10x/22, Quintuple, Infinity Plan Coverless 5x10x20x50x, Double Layer Mechanical Stage185x140/35x30mm, Reflect 6V30W Halogen + Polarizing	1 pcs / 1 ctn 39*37*65cm, 20/18Kgs			
213	A13.0212-BD	Metallurgical Microscope, BF/DF, DIC	Trinocular, WF10x/22, Quintuple, LWD Infinity Plan Coverless BF/DF 5x10x20x50x, Double Layer Mechanical Stage185x140/35x30mm, Reflect 12V50W Halogen, Polarizing Set	1 pcs / 1 ctn 39*37*65cm, 20/18Kgs			
214	A13.0212-DIC	Metallurgical Microscope, BF/DF, DIC	Trinocular, WF10x/22, Quintuple, Infinity Plan Coverless DIC 5x10x20x, Plan 50x, DIC Slide For 5x10x20x, Double Layer Mechanical Stage185x140/35x30mm, Reflect 12V50W + Polarizing	1 pcs / 1 ctn 39*37*65cm, 20/18Kgs			
215	A13.0212-BDDIC	Metallurgical Microscope, BF/DF, DIC	Trinocular, WF10x/22, Quintuple, Infinity Plan Coverless BF/DF DIC5x10x20x, Plan 50x, DIC Slide For 5x10x20x, Double Layer Mechanical Stage185x140/35x30mm, Reflect 12V50W + Polarizing	1 pcs / 1 ctn 39*37*65cm, 20/18Kgs			
216	A13.0901N-R	Metallurgical Microscope, BF+DF+PL+DIC, Semi-APO, Reflect & Transmit Light	5°~35° Tilting Trinocular Erect Image, PL10x/22, Quintuple Nosepiece With DIC Slot, Infinity LWD Plan BF/DF Metallurgical DIC 5x,10x,20x, Semi-APO 50x, Three Layers Mechanical Stage 455x240mm/ Reflect	1 pcs / 1 ctn 71.2*61.9*48.2cm, 35/29Kgs			
217	A13.0901N-RT	Metallurgical Microscope, BF+DF+PL+DIC, Semi-APO, Reflect & Transmit Light	5°~35° Tilting Trinocular Erect Image, PL10x/22, Quintuple Nosepiece With DIC Slot, Infinity LWD Plan BF/DF Metallurgical DIC 5x,10x,20x, Semi-APO 50x, Three Layers Mechanical Stage 455x240mm/ Reflect	1 pcs / 1 ctn 71.2*61.9*48.2cm, 35/29Kgs			
218	A13.0909-R	Metallurgical Microscope, BF+PL, Transmit & Reflect Light	Trinocular, PL10x/22mm, Quintuple, LWD Infinity Plan Metallurgical 5x10x20x50x, Double Layer Stage 175x145mm, Moving 76x42mm, Reflect 5W LED With Oblique Light, Polarizing	1 pcs / 1 ctn 67*45.5*61.5cm, 17/14.5Kgs			








No.	Cata. No.	Product	Specification	Package			Picture
219	A13.0909-RT	Metallurgical Microscope, BF+PL, Transmit & Reflect Light	Trinocular, PL10x/22mm, Quintuple, LWD Infinity Plan Metallurgical 5x10x20x50x, Double Layer Stage 175x145mm, Moving 76x42mm, Reflect 5W LED With Oblique Light, Polarizing, Transmit 5W LED Swing-Out	1 pcs / 1 ctn 67*45.5*61.5cm 17/14.5Kgs			
220	A13.0910-R	Metallurgical Microscope, BF+DF+PL+DIC, Semi-APO, Reflect & Transmit Light	Trinocular Erect Image, PL10x/25, Sextuple BF/DF Nosepiece With DIC Slot, LWD Infinity Plan BF/DF DIC Semi-APO 5x10x20x50x100x, 4" Triple Layer Mechanical Stage 230x215mm, Moving Range 105*102mm, Reflect	1 pcs / 1 ctn 89.4*56.3*52.5cm 40/0Kgs			
221	A13.0910-RT	Metallurgical Microscope, BF+DF+PL+DIC, Semi-APO, Reflect & Transmit Light	Trinocular Erect Image, PL10x/25, Sextuple BF/DF Nosepiece With DIC Slot, LWD Infinity Plan BF/DF DIC Semi-APO 5x10x20x50x100x, 4" Triple Layer Mechanical Stage 230x215mm, Moving Range 105*102mm, Reflect	1 pcs / 1 ctn 89.4*56.3*52.5cm 40/0Kgs			
222	A13.1005-B	Metallurgical Microscope, Bright Field/Dark Field	Trinocular, Inclined 30° Interpupillary 48-75mm, EW10X/22 Φ30, Infinite Plan Achromatic Objective (Bright Field) 5x10x20x50x100x, Backward Quintuple Nosepiece, Size Of The Soleplate 396X276X22mm, C Mount 1X,	1 pcs / 2 ctn 38*39*56cm, 22.7/20Kgs			
223	A13.1013-B	Metallurgical Microscope, BF/DF, PL, DIC	Trinocular Head, EW10x/22, Quintuple, Infinity Plan BF/DF Objectives 5x,10x,20x,50x,100x, Infinity Plan 40x,100x, Double Layer Mechanical Stage 186x138mm/74x50mm, Max. Sample Height 30mm, Specimen	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
224	A13.1017-B	Metallurgical Microscope, Reflect Light	Trinocular, WF10x, Quadruple, LWD Infinity Plan 4x10x40x, Non-Porous Mechanical Stage, Coaxial, Reflected 6V20W Halogen, Kohler Illumination	1 pcs / 1 ctn 36*26*62cm, 12/0Kgs			
225	A13.1018-B	Metallurgical Microscope, Transmit & Reflect Light	Trinocular, WF10x, Quadruple, LWD Infinity Plan 4x10x40x, Non-Porous Mechanical Stage 150x140mm/ 75x 50mm, Coaxial, Reflected 6V20W Halogen, Kohler Illumination, Transmit 6V20W Halogen	1 pcs / 1 ctn 36*26*62cm, 12/0Kgs			
226	A13.2603-B	Metallurgical Microscope, Reflect Light	Binocular, WF10X/18mm/1pc, WF10X/18mm Retical/1pc, LWD Infinity Plan 10x20x50x, Working Stage 140x200mm, Reflect 5W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 57*31*41cm, 9/0Kgs			
227	A13.2603-T	Metallurgical Microscope, Reflect Light	Trinocular, WF10X/18mm/1pc, WF10X/18mm Retical/1pc, LWD Infinity Plan 10x20x50x, Working Stage 140x200mm, Reflect 5W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 57*31*41cm, 9/0Kgs			
228	A13.2604-B	Metallurgical Microscope, Reflect & Transmit	Trinocular, 10X/22mm/1pc, 10X/22mm Retical/1pc, Quintuple, LWD Infinity Plan Semi-APO Objective 5x10x20x50x, Working Stage 182x140mm, Transmit 5W LED, Reflect 5W LED With Polaring Set, Safety Low	1 pcs / 1 ctn 51.5*47*43cm, 14/13Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
229	A13.2602-M	Inverted Metallurgical Microscope	Monocular, WF10x/18mm, Quaduple, Infinity LWD Plan 5x10x50x, Working Stage 210x180mm/30x30mm, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 44*27*53cm, 8/7Kgs			
230	A13.2602-B	Inverted Metallurgical Microscope	Binocular, WF10x/18mm, Quaduple, Infinity LWD Plan 5x10x50x, Working Stage 210x180mm/30x30mm, 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 44*27*53cm, 8/7Kgs			
231	A13.2604-B-D	Metallurgical Microscope, Reflect & Transmit	Trinocular, 10X/22mm/1pc, 10X/22mmRetical/1pc, Quaduple, LWD Infinity Plan BD Objective 5x10x20x50x, Working Stage 182x140mm, Transmit 5W LED, Reflect 12V50W Halogen, Both With Polaring	1 pcs / 1 ctn 51.5*47*43cm, 14/13Kgs			
232	A13.2605-A	Metallurgical Microscope, Reflect Light	Binocular, WF10x/20mm 1pc, WF10x/20mm Retical 1pc, LWD Infinity Plan Coverless 10x20x50x, Working Stage 140x200mm, Reflect 5W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 57*31*41cm, 9/0Kgs			
Inverted Metallurgical Microscope							
233	A13.2606	Invered Metallurgical Microscope, Reflect Light	Trinoclular, WF10x/22mm, WF10x/22 Reticule 0.1mm, Quinuple, Metallurgical LWD Infinity Plan Semi-APO Objective (BF) 5x10x20x50x, Three Layer Stage 210x180mm/50x 50mm, Coaxial Coarse & Fine,	1 pcs / 1 ctn 66*59*32.5cm, 17/12.5Kgs			
234	A13.2606-BD	Invered Metallurgical Microscope, Reflect Light	Trinoclular, WF10x/22mm, WF10x/22 Reticule 0.1mm, Quaduple, Metallurgical LWD Infinity Plan Objective (BD) 5x10x20x50x, Three Layer Stage 210x180mm/50x 50mm, Coaxial Coarse & Fine,	1 pcs / 1 ctn 66*59*32.5cm, 17/12.5Kgs			
235	A13.2607	Invered Metallurgical Microscope, BF, PL	Binocular, WF10x/20mm, WF10x/20 Reticule 0.1mm, Quaduple, Metallurgical LWD Infinity Plan Objective (BF) 5x10x20x50x, Three Layer Stage 180x155mm/ 75x40mm, Coaxial Coarse & Fine, 5W LED, Safety	1 pcs / 1 ctn 0*0*0cm, 17/12.5Kgs			
236	A13.2704	Inverted Metallurgical Microscope	Trinocular, WF10X/20, WF10X/18 reticule eyepiece. Infinity long work distance plan objectives 10x, 20x, 50x and 80x, Quintuple, Mechanical stage, 12V50W halogen lamp	0 pcs / 0 ctn 58*49*28cm, 20/0Kgs			
237	A13.0912-A	Inverted Metallurgical Microscope, BF+DF+PL+DIC, Semi-APO, ECO	Trinoclular, WF10x/22mm, Coded Quintuple Nosepiece With DIC Slot, Infinity LWD BD DIC Plan Achromatic Metallurgical 5x, 10x, 20x, Semi-APO FL50x, 3 Layers Mechanical Stage 240x250mm/ 50x50mm, Metal Plate	1 pcs / 1 ctn 49*51*73cm, 22/20Kgs			
238	A13.0900-R	Metallurgical Microscope, 4", BF+PL+DIC, Semi-APO, Reflect & Transmit Light	Trinocular, PL10x/22, Quintuple Nosepiece with DIC Slot, Infinity LWD Plan Metallurgical DIC 5x, 10x, 20x, Semi-APO 50x, 4" Metal Stage, Moving 105x105mm, 90V-240V, Reflect 5W LED, With Oblique Light, Brightness	1 pcs / 1 ctn 60*53.4*48.1cm , 30/25Kgs			








No.	Cata. No.	Product	Specification	Package			Picture
239	A13.0900-RT	Metallurgical Microscope, 4", BF+PL+DIC, Semi-APO, Reflect & Transmit Light	Trinocular, PL10x/22, Quintuple Nosepiece with DIC Slot, Infinity LWD Plan Metallurgical DIC 5x,10x,20x, Semi-APO 50x, 4" Glass Stage, Moving 105x105mm, Transmit & Reflected 5W LED, With Oblique Light, Brightness	1 pcs / 1 ctn 60*53.4*48.1cm , 30/25Kgs			
240	A13.0911-R	Metallurgical Microscope, 8", BF+DF+PL+DIC, Semi-APO, Reflect & Transmit Light	Trinocular, Erect Image, 30° Tilttable, PL10x/25mm, IPD 50-76mm, Light Split E100:P0/E0:P100, 6 Holes Nosepiece, With DIC Slot, Infinity LWD BD DIC Semi-APO Metallurgical 5x10x20x, Semi-APO 50x, 8" 3 Layers Stage	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
241	A13.0911-RT	Metallurgical Microscope, 8", BF+DF+PL+DIC, Semi-APO, Reflect & Transmit Light	Trinocular, Erect Image, 30° Tilttable, PL10x/25mm, IPD 50-76mm, Light Split E100:P0/E0:P100, 6 Holes Nosepiece, With DIC Slot, Infinity LWD BD DIC Semi-APO Metallurgical 5x10x20x, Semi-APO 50x, 8" 3 Layers Stage	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
242	A13.0913-R	Metallurgical Microscope, 12", BF+DF+PL+DIC, Semi-APO, Motorized Nosepiece	Trinocular, Erect Image, 0-35° Tilttable, PL10x/25mm, IPD 50-76mm, Light Split E100:P0/E0:P100, 6 Holes Motorized Nosepiece, With DIC Slot, Infinity LWD BD DIC Semi-APO Metallurgical 5x10x20x, Semi-APO 50x, 12x14" 3	1 pcs / 2 ctn 0*0*0cm, 0/0Kgs			
243	A13.0921-21	Precision Measuring Metallurgical Microscope, BF+DF+PL+DIC	Trinocular, Erect Image, 5-35° Tilttable, PL10x/22mm, Sextuple BF, LWD Infinity Plan Metallurgical 50x, DIC5x10x20x, Manual XY Stage 258x358mm, Moving 200x100mm, Marble Base 500x700x960mm/155Kgs,	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
244	A13.1091-TR	Metallurgical Microscope, Reflect, Transmit/Reflect, Semi-APO, BF/DF, DIC, PL	Tilting Trinocular 0-35°, EW10x/25mm, Manual Sextuple Nosepiece, NIS45 Infinity Plan Semi-APO BD 5x10x, LWD 20x50x100x, 4" Stage 210x170/102x102mm, LWD Condenser N.A.0.65, Coaxial, Reflect 12V100W	1 pcs / 2 ctn 79*58*52cm, 34/30Kgs			
245	A13.1093-TR	Metallurgical Microscope, Semi-Auto, Reflect, Transmit/Reflect, Semi-APO, BF/DF, DIC, PL	Tilting Trinocular 0-35°, EW10x/25mm, Auto Sextuple Nosepiece, NIS45 Infinity Plan Semi-APO BD 5x10x, LWD 20x50x100x, 4" Stage 210x170/102x102mm, LWD Condenser N.A.0.65, Coaxial, Reflect 12V100W	1 pcs / 2 ctn 79*58*52cm, 34/30Kgs			
246	A13.1095-TR	Metallurgical Microscope, Full Auto, Reflect, Transmit/Reflect, Semi-APO, BF+DF+DIC+PL	Tilting Trinocular 0-35°, EW10x/25mm, Auto Sextuple Nosepiece, NIS45 Infinity Plan Semi-APO BD 5x10x, LWD 20x50x100x, 4" Stage 210x170/125x75mm, LWD Condenser N.A.0.65, Coaxial, Reflect 12V100W Halogen	1 pcs / 2 ctn 79*58*52cm, 34/30Kgs			
247	A13.1096	Inverted Metallurgical Microscope, Semi-APO, BF/DF/PL/FL/DIC	Trinocular, Built-in Bertrand Lens, SW10x/22mm, Sextuple Nosepiece With DIC Slot, Infinity Plan Sem-APO BD 5x10x20x APO BD50x100x, Triple Layer Working Stage 340x230mm, Moving Range 130x85mm, With 3	1 pcs / 1 ctn 120*80*70cm, 60/50Kgs			
248	A13.1098	Full Motorized Inverted Metallurgical Microscope, Semi-APO, BF/DF/PL/FL/DIC	Trinocular, Built-in Bertrand Lens, SW10x/22mm, Motorized Sextuple Nosepiece With DIC Slot, Infinity Plan Sem-APO BD 5x10x20x APO BD50x100x, XYZ Motorized Working Stage With Joy Stick, Coaxial Coarse	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			



No.	Cata. No.	Product	Specification	Package		Picture
Inverted Metallurgical Microscope						
249	A15.1302-A	Polarizing Microscope	Monocular, WF10x Retical, Quadruple, Ach. 4x,10x,40x,100x, Round Stage Dia.120mm, Coaxial Focusing, Abbe N.A.1.25, 1W LED, Polarizer, Analyzer, Optical Compensator λ Slip (First Class Red), 1/4 λ Slip	1 pcs / 1 ctn 37*23*47cm, 5/4Kgs		
250	A15.1302-B	Polarizing Microscope	Binocular, WF10x/18 +WF10x Retical, Quadruple, Ach. 4x,10x,40x,100x, Round Stage Dia.120mm, Coaxial Focusing, Abbe N.A.1.25, 1W LED, Polarizer, Analyzer, Optical Compensator λ Slip (First Class Red), 1/4 λ Slip	1 pcs / 1 ctn 37*23*47cm, 6/4.5Kgs		
251	A15.1302-C	Polarizing Microscope	Trinocular, WF10x/18 +WF10x Retical, Quadruple, Ach. 4x,10x,40x,100x, Round Stage Dia.120mm, Coaxial Focusing, Abbe N.A.1.25, 1W LED, Polarizer, Analyzer, Optical Compensator λ Slip (First Class Red), 1/4 λ Slip	1 pcs / 1 ctn 37*23*47cm, 6/4.5Kgs		
252	A15.2604	Polarizing Microscope	Monocular, WF10x/18mm With Cross Reticule 0.1mm, Quadruple, Strain-free Achromatic 4x10x40x, Analyzer Unit 0-90°, Bertrand Lens, Compensator λ & 1/4 λ Slip, Polarizing Unit, Round Stage Dia. 142mm, Strain-free Abbe	1 pcs / 1 ctn 34*25*42cm, 6.5/5Kgs		
253	A15.1017-T	Polarizing Microscope	Trinocular, WF10x Retical + WF10x, Quadruple Center Adjustable, Strain Free Achromatic 4x10x40x, Round Stage, Abbe 1.25, Coaxial Coarse&Fine, Halogen 6V20W, Swing-Out Polarizer On Collector, Analyzer	1 pcs / 1 ctn 33*28*43cm, 10/0Kgs		
254	A15.0204-A	Polarizing Microscope, Transmit Light	Trinocular, WF10x/22mm, Retical WF10x/22mm, Quadruple, Strain-free Infinity Plan 4x10x40x60x, Round Stage Dia.150mm, Transmit Halogen 6V30W With Polarizer, Intermedia Analyzer With Bertrand Lens, λ , $\lambda/4$ And Quarts	1 pcs / 1 ctn 36*31*68cm, 13/18Kgs		
255	A15.0204-B	Polarizing Microscope, Transmit Light	Trinocular, WF10x/22mm, Retical WF10x/22mm, Quadruple, Strain-free Infinity Plan Coverless 5x10x40x60x, Round Stage Dia.150mm, Transmit Halogen 6V30W, Reflect Halogen 6V30W, With Polarizer, With Bertrand	1 pcs / 1 ctn 36*31*68cm, 13/18Kgs		
256	A15.0701-T-T	Polarizing Microscope, Transmit	Seidentopf Trinocular Viewing Head, Inclined at 30°, WF10x/20 With Cross Hair, WF10x/20 With Graticule, WF10x/20, Infinity Plan Achromatic Objective 4x,10x for Reflected/Transmitted System, Infinity	1 pcs / 1 ctn 39*28*64cm, 15/0Kgs		
257	A15.0701-TR-T	Polarizing Microscope, Transmit & Reflect	Seidentopf Trinocular Viewing Head, Inclined at 30°, WF10x/20 With Cross Hair, WF10x/20 With Graticule, WF10x/20, Infinity Plan Achromatic Objective 4x,10x for Reflected/Transmitted System, Infinity	1 pcs / 1 ctn 39*28*64cm, 15/0Kgs		
258	A15.2603-B	Polarizing Microscope, Transmit	Trinocular, WF10x/20/1pc, WF10x/20 Retical/1pc, Quadruple Center Adjustable, Strain-free Infinity Plan 4x10x20x40x, Round Stage Dia.142mm, λ + $\lambda/4$, Quartz Wedge, Transmint Light 5W LED, Low Voltage	1 pcs / 1 ctn 57*28*45cm, 8.6/7.5Kgs		

No.	Cata. No.	Product	Specification	Package			Picture
259	A15.2601-R-T	Polarizing Microscope, Reflect	Trinocular, WF10x/20mm 1pc, WF10x/20mm Rectical 1pc, Qunituple Nosepiece Center Adjustable, LWD Infinity Strain-free Plan 5x10x20x50x, Round Working Stage Dia.148mm, λ , $\lambda/4$, Quartz Wedge, Reflect Light 5W LED	1 pcs / 1 ctn 52*47*43cm, 16/12.5Kgs			
260	A15.2601-RT-T	Polarizing Microscope, Transmit & Reflect	Trinocular, WF10x/20mm 1pc, WF10x/20 Rectical 1pc, Qunituple Nosepiece Center Adjustable, Infinity Strain-free Plan 4x10x20x40x, LWD Infinity Strain-free Plan 50x, Round Working Stage Dia.174mm, λ , $\lambda/4$, Quartz Wedge, Reflect Light 5W LED	1 pcs / 1 ctn 75*36*45cm, 16/12Kgs			
261	A15.2601-T-T	Polarizing Microscope, Transmit	Trinocular, WF10x/22mm 1pc, WF10x/22 Rectical 1pc, Qunituple Center Adjustable, Infinity Strain-free Plan 4x10x20x40x, Round Working Stage Dia.174mm, Swing-Out Condenser, Transmit 5W LED, λ , $\lambda/4$, Quartz Wedge, Reflect Light 5W LED	1 pcs / 1 ctn 57*31*41cm, 12.5/10Kgs			
262	A15.2701 -B	Polarizing Microscope	Trinocular, 40X ~ 1000X, WF10X/20, Plan Polarizing reticule eyepiece 10x/20mm. Strain-free infinity achromatic 4x, 10x, 40x(S) and 100x (S/Oil) objectives, quadruple nosepiece, Swing-out condenser	1 pcs / 1 ctn 48*48*32cm, 19/10Kgs			
263	A15.2702-B	Polarizing Microscope	Trinocular, 40x-1000x, WF10X/20mm, plan polarizing reticule eyepiece 10x/20mm, strain-free infinity achromatic 4x, 10x, 40x(S) and 100x (S/Oil) objectives. coaxial coarse and fine focus. quadruple, swing-out condenser	1 pcs / 1 ctn 48*48.5*68cm, 20/0Kgs			
264	A15.0908-T	Polarizing Microscope, BF+PL, Transmit & Reflect	Binocular, PL10x/22mm, PL10x/22 Cross, Quaduple, Infinity Plan Strain-Free Objectives 4x10x40x60x, Intermedia Analyzer Set, Circular Rotatable Working Stage Dia.160mm, Swing-out Condenser NA0.9, With Polarizer, Reflect Light 5W LED	1 pcs / 1 ctn 68*44*49cm, 15/13Kgs			
265	A15.0908-R	Polarizing Microscope, BF+PL, Transmit & Reflect	Binocular, PL10x/22mm, PL10x/22 Cross, Quaduple, LWD Infinity Plan Destressing 5x10x20x50x, Intermedia Analyzer Set, Circular Rotatable Working Stage Dia.160mm, Swing-out Condenser NA0.9 With Polarizer, Reflect Light 5W LED	1 pcs / 1 ctn 68*44*49cm, 15/13Kgs			
266	A15.0908-RT	Polarizing Microscope, BF+PL, Transmit & Reflect	Binocular, PL10x/22mm, PL10x/22 Cross, Quaduple, LWD Infinity Plan Destressing 5x10x20x50x, Intermedia Analyzer Set, Circular Rotatable Working Stage Dia.160mm, Swing-out Condenser NA0.9 With Polarizer, Reflect Light 5W LED	1 pcs / 1 ctn 68*44*49cm, 15/13Kgs			
267	A15.1091-T	Polarizing Microscope, Semi-APO, Transmit & Reflect Light	Tilting Trinocular 0-35°, EW10x/22mm, EW10x/22 Scale Cross, Quintuple Nosepiece Center Adjustable, NIS60 Infinity Plan Strainless Semi-APO PL5x10x20x50x100x, Round Stage Dia.190mm, With Attachable Moving Stage	1 pcs / 1 ctn 76*55*50cm, 24/20Kgs			
268	A15.1091-TR	Polarizing Microscope, Semi-APO, Transmit & Reflect Light	Tilting Trinocular 0-35°, EW10x/22mm, EW10x/22 Scale Cross, Quintuple Nosepiece Center Adjustable, NIS60 Infinity Plan Strainless Semi-APO PL5x10x20x50x100x, Round Stage Dia.190mm, With Attachable Moving Stage	1 pcs / 2 ctn 79*58*52cm, 32/30Kgs			


No.	Cata. No.	Product	Specification	Package			Picture
269	A15.1091-R	Polarizing Microscope, Semi-APO, Transmit & Reflect Light	Tilting Trinocular 0-35°, EW10x/22mm, EW10x/22 Scale Cross, Quintuple Nosepiece Center Adjustable, NIS60 Infinity Plan Strainless Semi-APO PL5x10x20x50x100x, Round Stage Dia.190mm, With Attachable Moving	1 pcs / 2 ctn 79*58*52cm, 32/30Kgs			
Upright Fluorescent Microscope							
270	A16.0204	Flourescent Microscope	Trinocular, WF10x, Quarduple, Plan 4x10x100x, Fluorescence 25x, Mechanical Stage 135x125mm, Coaxial, Halogen 6V20W, Mercury Lamp 100W + Fluorescence B,G	1 pcs / 1 ctn 43.5*30*49.5cm , 14.5/0Kgs			
271	A16.0207-2	LED Flourescent Microscope	Trinocular, WF10x/22, Quarduple, Infinity Plan 4x10x100x, PLFL40x, Mechanical Stage 210x140mm/75x50mm, Transmit 3W LED, Incident 3W LED Fluorescence B,G,	1 pcs / 1 ctn 37*33*69cm, 15/0Kgs			
272	A16.2601-NL-1	LED Fluorescent Microscope	Trinocular, WF10x/20, Quadruple, Infinity E-Plan 4x10x40x100x, Double Layer Mechanical Stage 145x140/76x52mm, 3W LED, 5W LED Fluorescent, China B or G Filter, Safety Low Voltage Charger	1 pcs / 1 ctn 42*28*45cm, 18/15Kgs			
273	A16.2601-NL-2	LED Fluorescent Microscope	Trinocular, WF10x/20, Quadruple, Infinity E-Plan 4x10x40x100x, Double Layer Mechanical Stage 145x140/76x52mm, 3W LED, 5W LED Fluorescent, China B+G Filter, Safety Low Voltage Charger	1 pcs / 1 ctn 42*28*45cm, 18/15Kgs			
274	A16.2603-NL-1	LED Fluorescent Microscope, Semi-APO	Trinocular, WF10x/22, Quinduple, Infinity Plan 4x10x20x40x100x, Rackless Mechanical Stage 182x140mm, Swing Out Condenser, Halogen 6V30W, 5W LED Fluorescent, China B or G Filter	1 pcs / 1 ctn 51*57*29cm, 12.5/11Kgs			
275	A16.2603-NL-2	LED Fluorescent Microscope, Semi-APO	Trinocular, WF10x/22, Quinduple, Infinity Plan 4x10x20x40x100x, Rackless Mechanical Stage 182x140mm, Swing Out Condenser, Halogen 6V30W, 5W LED Fluorescent, China B+G Filter	1 pcs / 1 ctn 51*57*29cm, 12.5/11Kgs			
276	A16.0910	Research Fluorescence Microscope, B,G,UV, Semi-APO, Upgrade BF+DF+PL+PH+FI+DIC	Trinocular, PL10x/25, Sextuple BF Nosepiece With Slot For DIC/PL, Infinity Plan Semi-APO Fluorescent 4x,10x,20x,40x,100x, Double Mechanical Stage 187*166mm, Moving Range 85*55mm, Swing-Out	1 pcs / 1 ctn 89*51.5*81.5cm , 39/26Kgs			
277	A16.1062-T	LED Fluorescent Microscope, B,G	Trinocular, WF10x/22, Coding Quintuple, Infinity Plan NIS60 4x10x40x100x, Rackless Stage 230x150/78x54mm, N.A.1.25 Condenser With Slot, Coaxial, 3W LED, LCD Screen Base, Fluorescent Filter	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
278	A16.1091-H	Upright Fluorescent Microscope, Semi-APO, Manual	Trinocular, EW10x/25mm, Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler Illumination, 12V100W Halogen , ECO	1 pcs / 2 ctn 79*58*52cm, 34/30Kgs			











No.	Cata. No.	Product	Specification	Package			Picture
279	A16.1093-H	Upright Fluorescent Microscope, Semi-APO, Semi-Motorized	Tilting Trinocular 0-35°, EW10x/25mm, Auto Sextuple Nosepiece, Infinity Plan Semi-APO Fluorescent 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Auto Swing-Out Condenser NA0.9/0.25, Coaxial,	1 pcs / 2 ctn 79*58*52cm, 34/30Kgs			
280	A16.1095	Upright Fluorescent Microscope, Semi-APO, Full Motorized	Trinocular, EW10x/25mm, Motorized Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Motorized X/Y/Z Control Stage 190x152/ 125x75mm, Motorized Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler Illumination	1 pcs / 2 ctn 79*58*52cm, 34/30Kgs			
281	A16.2601-T2	Fluorescent Microscope	Trinocular, WF10x/20, Quaduple, Infinity EPlan 4x10x40x100x, 3W LED, Safety Low Voltage Charger, Disc 100W Mercury Fluorescent, China B/G Filter	1 pcs / 1 ctn 97*37*44cm, 22/18Kgs			
282	A16.2601-L-T2	LED Fluorescent Microscope	Trinocular, WF10x/20, Quaduple, Infinity E-Plan 4x10x40x100x, 3W LED, Safety Low Voltage Charger, 5W Disc LED Fluorescent, China B,G Filter	1 pcs / 2 ctn 42*28*45cm, 18/15Kgs			
283	A16.2603-T2	Fluorescent Microscope, Semi-APO	Trinocular, WF10x/22, Quinduple, Infinity Plan 4x10x20x40x100x, Rackless Mechanical Stage 182x140mm, Swing Out Condenser, Halogen 6V30W, Disc 100W Mercury Fluorescent, China B/G Filter	1 pcs / 1 ctn 99*34*44cm, 23/21Kgs			
284	A16.2603-L-T2	LED Fluorescent Microscope, Semi-APO	Trinocular, WF10x/22, Quinduple, Infinity Plan 4x10x20x40x100x, Rackless Mechanical Stage 182x140mm, Swing Out Condenser, Halogen 6V30W, 5W Disc LED Fluorescent China B/G	1 pcs / 1 ctn 52*47*49cm, 23/21Kgs			
285	A16.0208-2	LED Fluorescent Microscope, For Phthisis Checking	Trinocular, WF10x/22, Quintuple, Infinity Plan 4x10x100x, Infinity Plan Fluorescent 40x, Double Layer Mechanical 180x150mm/75x50mm, Transmit 3W LED, Epi Fluorescent 5W LED, Fluorescent Filter B+G	1 pcs / 1 ctn 37*33*69cm, 15/0Kgs			
286	A16.0907-L	LED Fluorescence Microscope, Phthisis Checking	Trinocular, PL10x/20mm, Quaduple, Infinity Plan 10x20x40x100x, Working Stage 150x140/76x50mm, Condenser N.A.1.25 With Slot For Phase Contrast & Dark Field, With Dark Field Slide, Kohler 3W LED, Plug Type Lamp	1 pcs / 1 ctn 56.4*50.7*56.6cm, 13/12Kgs			
287	A16.0911-P	LED Fluorescence Microscope, Phthisis Checking	Trinocular, Eyepiece PL10X/20mm, Infinity Plan Achromatic 10x20x40x100x, 5 Holes Nosepiece, Rackless Working Stage 150x162/76x50mm, Abbe Condenser N.A.1.25 With Slot For Dark Field,	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
288	A16.0908-L	LED Fluorescence Microscope, B,G,UV, Semi-APO, B4 Tuberculosis Phthisis Checking Optional	Trinocular, PL10x/22mm, Quintuple, Infinity Plan Semi-APO Fluorescent 4x10x20x40x100x, Swing Out Condenser NA1.2/0.22, Double Layers 175x145/76x50mm, 6V30W Halogen, Disc LED Fluorescent Unit (4 Holes),	1 pcs / 1 ctn 58*46*67cm, 23/20Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
289	A16.2603-P-L	LED Fluorescent Microscope, Phthisis Checking	Trinocular, WF10x/22, Quintuple, Infinity Plan 10x40x60x100x, Rackless Mechanical Stage 182x140mm, Swing Out Condenser, 5W LED, LED 5W Fluorescent 455nm (Chroma Filter) For Phthisis Checking	1 pcs / 1 ctn 75*34*44cm, 17/15Kgs			
Inverted Flourescent Microscope							
290	A16.0206	Inverted Flourescent Microscope, Phase Contrast	Trinocular, WF10x/22, Quintuple, LWD Infinity Plan 10x20x40x, LWD Infinity Plan Phase Contrast 10x20x40x, Fixed Stage Overall Size 227x208mm, Glass Round Stage Size Dia.118mm, Moving Range 77x114mm, Incident	1 pcs / 1 ctn 83*39*53cm, 26/24Kgs			
291	A16.0912-L	Inverted LED Fluorescence Microscope, Semi-APO, BF+PH+FL, ECO	Trinocular, WF10x/22mm, Quintuple, LWD Infinity Plan Fluorescent 4x10x, LWD Infinity Plan Semi-APO PH 20x40x, Phase Contrast Slide 4x10+20x40x, Centering Telescope, Fixed Working Stage 215x250mm,	1 pcs / 1 ctn 74*67*57cm, 24/20Kgs			
292	A16.1063	Inverted LED Fluorescent Microscope, B,G,U	Binocular (Photo Port On Side), WF10x/22, Quintuple, LWD Infinity Plan Semi-APO Fluorescent Objective 4x10x20x40x, Plain Stage 170x250mm, Attachable Moving Ruler 128x80mm, Glass Plate, Auxilliary Plate, Universal	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
293	A16.1064	Inverted LED Fluorescent Microscope, B,G,U	Binocular (Photo Port On Side), WF10x/22, Coded Quintuple Nosepiece, LWD Infinity Plan Semi-APO Fluorescent Objective 4x10x20x40x, Plain Stage 170x250mm, Attachable Moving Ruler 128x80mm, Glass Plate,	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
294	A16.1065	Inverted LED Fluorescent Microscope, B,G,U	Tilttable Binocular (Photo Port On Side), WF10x/22, Coded Quintuple Nosepiece, LWD Infinity Plan Semi-APO Fluorescent Objective 4x10x20x40x, Plain Stage 170x250mm, Attachable Moving Ruler 128x80mm, Glass Plate,	1 pcs / 2 ctn 69*39*65cm, 20/18Kgs			
295	A16.1096	Inverted Fluorescent Microscope, Manual, Semi-APO, BF+PH+FL, DF/PL/DIC	Trinocular, Built-in Bertrand Lens, SW10x/22mm, Sextuple Nosepiece With DIC Slot, Infintiy Plan Sem-APO PH10x20x40x, Triple Layer Working Stage Moving Range 130x85mm, With 3 Holder, LWD Condensor N.A.0.55	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
296	A16.1097	Inverted Fluorescent Microscope, Semi-Auto, Semi-APO, BF+PH+FL, DF/PL/DIC	Trinocular, Built-in Bertrand Lens, SW10x/22mm, LCD Touch Screen + Coded Sextuple Nosepiece With DIC Slot, Infintiy Plan Sem-APO PH10x20x40x, Triple Layer Working Stage Moving Range 130x85mm, With	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
297	A16.1098	Full Motorized Inverted Fluorescent Microscope, Full Auto, Semi-APO, BF+PH+FL, DE/PI/DIC	Trinocular, Built-in Bertrand Lens, SW10x/22mm, Motorized Sextuple Nosepiece With DIC Slot, Infintiy Plan Sem-APO PH10x20x40x, Motorized Triple Layer Working Stage 325x144mm Moving Range 130x85mm, With	1 pcs / 3 ctn 0*0*0cm, 52/40Kgs			
298	A16.2614-2	Inverted Fluorescent Microscope	Trinocular, WF10x/22, Quintuple, Infinity Plan LWD 10x20x40x, Infinity Phase Contrast LWD 20x Objective, 10x+20x/40x Phase Contrast Annular Spot, Transmit Halogen 6V30W, Reflect Fluorescent 100W Mercury	1 pcs / 2 ctn 66*59*33cm, 24.8/19.6Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
299	A16.2614-L-2	Inverted LED Fluorescent Microscope	Trinocular, WF10x/22, Quintuple, Infinity Plan LWD 10x20x40x, Infinity Phase Contrast LWD 20x Objective, 10x+20x/40x Phase Contrast Annular Spot, Transmit Halogen 6V30W, Reflect 5W LED Fluorescent	1 pcs / 2 ctn 66*59*33cm, 24.8/19.6Kgs			
300	A16.2615-L-2	Inverted LED Fluorescent Microscope	Trinocular + Dual Photo Output, WF10x/22, Quintuple, Infinity Plan LWD 4x, Infinity Phase Contrast LWD 10x20x40x, Turret Phase Contrast Condenser, Incl. Attachable Stage For 96 Holes Plate, ECO, Transmit Halogen	1 pcs / 1 ctn 66*60*33cm, 17/13.6Kgs			
301	A16.2702-2i	Inverted Fluorescent Microscope	Binocular + Side Photo Port, WF10x/20mm, Quintuple, LWD Infinity Plan 10x20x40x, LWD Infinity Plan Phase Contrast 20x, Double Layer Mechanical Stage 240x260mm/135x85mm, With Water Drop Holder +	1 pcs / 2 ctn 69*68*42cm, 43/0Kgs			

Multi-Viewing Microscope

302	A17.0940	Multi-Viewing Microscope, 10 People	1pc Trinocular Head, 9pcs Binocular Heads, Plan WF10X/22mm, Quintuple Nosepiece (5 Holes), Infinity Plan 4x10x20x40x100x, Double Layer Stage 175x145mm, 10W LED, Swing-out Condenser	1 pcs / 1 ctn 43.9*41.6*55cm 15.5/11Kgs			
303	A17.0950-5	Multi View Microscope	5 People, 1xTrinocular, 4xBinocular, PL10x/22mm, Sextuple BF Nosepiece With Slot For DIC/PL, Infinity Plan 4x, 10x, 20x, 40x, 100x, Double Mechanical Stage 187*166mm, Moving Range 80*55mm, Swing-Out	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
304	A17.0950-3	Multi View Microscope	3 People, 1xTrinocular, 2xBinocular, PL10x/22mm, Sextuple BF Nosepiece With Slot For DIC/PL, Infinity Plan 4x, 10x, 20x, 40x, 100x, Double Mechanical Stage 187*166mm, Moving Range 80*55mm, Swing-Out	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
305	A17.0950-10	Multi View Microscope	10 People, 1xTrinocular, 9xBinocular, PL10x/22mm, Sextuple BF Nosepiece With Slot For DIC/PL, Infinity Plan 4x, 10x, 20x, 40x, 100x, Double Mechanical Stage 187*166mm, Moving Range 80*55mm, Swing-Out	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
306	A17.1091-2	Multi Viewing Research Laboratory Microscope, 10 People, Manual	2 Heads Opsite, Trinocular, EW10x/25mm, Auto Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Auto Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			
307	A17.1091-2R	Multi Viewing Research Laboratory Microscope, 10 People, Manual	2 Heads Right Side, Trinocular, EW10x/25mm, Auto Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Auto Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			
308	A17.1091-2L	Multi Viewing Research Laboratory Microscope, 10 People, Manual	2 Heads Left Side, Trinocular, EW10x/25mm, Auto Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Auto Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			









No.	Cata. No.	Product	Specification	Package			Picture
309	A17.1091-3	Multi Viewing Research Laboratory Microscope, 10 People, Manual	3 Heads, Trinocular, EW10x/25mm, Auto Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Auto Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler Illumination, 24V100W	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			
310	A17.1091-5	Multi Viewing Research Laboratory Microscope, 10 People, Manual	5 Heads, Trinocular, EW10x/25mm, Auto Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Auto Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler Illumination, 24V100W	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			
311	A17.1091-10	Multi Viewing Research Laboratory Microscope, 10 People, Manual	Deca Heads, Trinocular, EW10x/25mm, Sextuple Nosepiece, Infinity Plan NIS60 4x10x20x40x100x, Rackless Stage 190x152/ 78x54mm, Swing-Out Condenser NA0.9/0.25, Coaxial, Kohler Illumination, 24V100W Halogen/3W LED	1 pcs / 1 ctn 79*58*52cm, 24/20Kgs			
312	A17.1112	Multi-Viewing Microscope, 10 People	1pc Trinocular Head, 9pcs Binocular Heads, WF10X/22mm, Quadruple, Infinity Plan 4x,10x,40x,100x, Mechanical Stage 180x160mm/ 80x50mm, Kohler Illumination, 12V100W Halogen, LED Cursor	1 pcs / 3 ctn 0*0*0cm, 0/0Kgs			
313	A17.1113	Multi-Viewing Microscope, 10 People	1pc Trinocular Head, 9pcs Binocular Heads, WF10X/22mm, Quadruple, Infinity Plan 4x,10x,40x,100x, Mechanical Stage 180x160mm/ 80x50mm, Kohler Illumination, 5W LED Cursor Brightness Adjustable	1 pcs / 3 ctn 0*0*0cm, 0/0Kgs			
Forensic Comparison Microscope							
314	A18.1816-A	Questioned Document Examination System	1~30x Motorized Zoom, 1/3" CCD 800 Lines Digital Camera, Diffuse Light 1x24V150W, UVC 254nm 1x10W, UVA 365nm 1x10W, Side Light 1x50W, Transmit UV 2x10W, Transmit White 2x20W, Exciting Filter 9 Bands,	1 pcs / 1 ctn 100*69*78cm, 76/70Kgs			
315	A18.1818	Questioned Document Examination System	1~30x Motorized Zoom, 1/3" CCD 800 Lines Digital Camera, Incl. Software+Laptop, Wooden Packing, Light Source : --Strong Light 1x24V/150W Halogen --Ultra Violet 365nm, 1x10W Check Filter	1 pcs / 1 ctn 80*60*70cm, 50/45Kgs			
316	A18.1825	Comparison Microscope, Step Zoom, 3.2~192x	3.8~230x, Binocular 45° Inclined, 48-75mm, WF10x/22, WF20x/12, Step Zoom Lens 0.8x,1.25x,2x,3.2x,4.8x, Auxiliary Lens 0.4x.2.0x, Bridge Magnification 1.2x, 5W LED+12V/50W Halogen, Coaxial Light, Polarizing Light, UV Filter, Filter Red/Blue/Yellow/Green, 25.0mm	1 pcs / 1 ctn 55*49*76cm, 25/30Kgs			
317	A18.1825-LCD	LCD Comparison Microscope, Step Zoom, 3.2~192x	3.8~230x, Binocular 45° Inclined, 48-75mm, WF10x/22, WF20x/12, Step Zoom Lens 0.8x,1.25x,2x,3.2x,4.8x, Auxiliary Lens 0.4x.2.0x, Bridge Magnification 1.2x, 5W LED+12V/50W Halogen, Coaxial Light, Polarizing Light, UV Filter, Filter Red/Blue/Yellow/Green, 25.0mm	1 pcs / 1 ctn 68*50*85cm, 40/35Kgs			
318	A18.1826	Comparison Microscope, 3.36~216x	3.36~216x, WF10x,20x, Zoom Lens 0.7-4.5x, View Field 2mm-60mm, Auxiliary Lens 0.4x.2.0x, Bridge Magnification 1.2x, 5W LED+12V/50W Halogen, Coaxial Light, Polarizing Light, UV Filter, Filter Red/Blue/Yellow/Green, 25.0mm	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			











No.	Cata. No.	Product	Specification	Package			Picture
319	A18.1828	Comparison Microscope, 2~240x	2~240x, Binocular, WF10x/22mm, WF20x/12mm, 6 Holes Nosepiec, LWD Infinity Plan Objective 0.5x,1x,1.5x,2x,4x,6x,0.5x, Auxiliary Lens 0.4x,2.0x, 5W LED+12V/50W Halogen, Coaxial Light, Polarizing Light,	1 pcs / 1 ctn 68*50*85cm, 45/40Kgs			
320	A18.1828-LCD	LCD Comparison Microscope, 2~240x	2x~240x, Binocular, WF10x/22mm, WF20x/12mm, 6 Holes Nosepiec, LWD Infinity Plan Objective 0.5x,1x,1.5x,2x,4x,6x,0.5x, Auxiliary Lens 0.4x,2.0x, 5W LED+12V/50W Halogen, Coaxial Light, Polarizing Light,	1 pcs / 1 ctn 68*50*85cm, 50/45Kgs			
321	A18.1829	Motorized Digital Comparison Microscope, 2x~240x	2x-240x, 0-45° Binocular, Eyepiece WF10x/22, WF20x/12, Sextuple Nosepiece, LWD Infinity Plan Objectives 0.5x,1x,1.5x,2x,4x,6x, Auxiliary Lens 4x,2.0x, Upper Light 5W LED,Side Light12V/50W Cooling Halogen, Bottom	1 pcs / 3 ctn 58*78*102cm, 85/50Kgs			
322	A18.1830	Motorized Digital Comparison Microscope, 2x~240x, With Motorized Table	2x-240x, 0-45° Binocular, Eyepiece WF10x/22, WF20x/12, Sextuple Nosepiece, LWD Infinity Plan Objectives 0.5x,1x,1.5x,2x,4x,6x, Auxiliary Lens 4x,2.0x, Upper Light 5W LED,Side Light12V/50W Cooling Halogen, Bottom	1 pcs / 3 ctn 0*0*0cm, 130/110Kgs			
323	A18.1831	Comparison Microscope, Step Zoom, 3.2~288x	3.2x~288x, Binocular, WF10x/22, WF20x/12, Step Zoom 0.8x,1.25x,2x,3.2x,4.8x, Auxiliary Lens 0.4x,2.0x, 5W LED+12V/50W Halogen+108 LEDs Ring Light, UV Filter, Filter Red/Blue/Yellow/Green,	1 pcs / 1 ctn 68*50*85cm, 45/40Kgs			
324	A18.4902	Dual Light Path Full Digital Comparison Microscope	Optical Magnfication 0.83x-10x, Total Magnfication 57x-715x (21"LCD), Digital Camera x2 8.3M, CMOS 1/2.5", 3840x2160, W.D.72mm, View of Field 0.65x0.36~8.2x4.6mm, Double Layer Mechanical Stage x2, Size	1 pcs / 1 ctn 87*61*41cm, 33.9/18.5Kgs			
Phase Contrast Microscope							
325	A19.0209-40	Phase Contrast Microscope, Achromatic Phase Contrast 40x	Binocular, WF10x, Quarduple, Achromatic 10x20x100x, Phase Contrast Objective 40x, Double Layer Mechanical Stage, 140*140mm, Mvoing Range 50*75mm, Phase Contrast Slide 40x, Telescope 10x, Phase Contrast	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
326	A19.0214	Phase Contrast Microscope, Achromatic Turret Phase Contrast	Trinocular, WF10x/18mm, Quarduple, Achromatic Phase Contrast Objective 10x20x40x100x, Double Layer Mechanical Stage 160x140mm/75x50mm, Turret Disc Phase Contrast Condenser, Centering	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
327	A19.0911	Phase Contrast Microscope, BF+PH, Coded Nosepiece	Trinocular, WF10x/20mm, Infinity Plan 4x, 5 Holes Coded Nosepiece, Double Layer Mechanical 150x162mm, Moving 76x50mm, 3W LED Kohler, With LCD Screen, Light Temperature Adjustable (3000-5600K), Brightness Memory	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
328	A19.1008-40	Phase Contrast Microscope	Binocular, WF10x16x, Quarduple, Achromatic 4x10x100x, 4 Holes, Double Layer 140x140mm, 3W LED, Phase Contrast Objective 40x, Abbe Condenser N.A.1.25, Phase Contrast	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			



No.	Cata. No.	Product	Specification	Package			Picture
329	A19.2701-B	Phase Contrast Microscope	Binocular, simple PH attachment, infinity phase contrast objective	1 pcs / 1 ctn 48*48*32cm, 9/0Kgs			
330	A19.2702-B	Phase Contrast Microscope, Infinity Plan, Turret Disc Condenser	Binocular, WF10x/22, Quaduple, Infinity Plan Phase Contrast 10x20x40x100x, Coaxial, Mechanical 156x138mm, Moving 76x52mm, Turret Disc PH Condenser, Centering Telescope, Kohler 3W LED	0 pcs / 0 ctn 50*50*34cm, 13.5/0Kgs			
Dark Field Microscope							
331	A10.0202	Dark Field Microscope, 12V50W Halogen Lamp	Trinocular, WF10x/20, Quintple, Plan 4x10x20x40x100x, Double Layers Mechanical 160x140mm/ 70x50mm, Abbe NA1.25, Coaxial Coarse & Fine Focusing, Halogen 12V50W Kohler, Dark Field Condenser Dry + Immersion	1 pcs / 1 ctn 45*30*35cm, 10.5/9Kgs			
332	A10.1007	Dark Field Microscope	Trinocular, WF10x, Quadruple, Achromatic 4x10x40x100x, Double Layer Mechanical 140x140mm, Coaxial Coarse&Fine, Abbe 1.25, LED	1 pcs / 1 ctn 39*35*52cm, 12/10Kgs			
333	A10.1018	Professional Dark Field Microscope	A55.1002-45 C-Mount 0.45x Trinocular, EW10x/22mm, Quintuple, Infinity Plan Objective 4x10x40x100x, Mechanical Stage 185x142/ 55x75mm, Swingout NA0.9/0.25, Coaxial Coarse & Fine, External Halogen 6V30W, Kohler Illumination Asperical Collector, Eyepiece	1 pcs / 1 ctn 43*39*61cm, 15/14Kgs			
Mono Zoom Microscope							
334	A21.1609	3D Motorized Rotate Zoom Mono Video Microscope	45°/60° Optional, Lens 0.28x-1.875x, W.D. 26mm/45°, 42mm/60°, LED Light, Mini Motor 5-27 Turns/Min, Stepless Adjustment, Stand SD11-D, Moving Stage Range 100x86mm, 0.5x CCD Adapter, Control Box	1 pcs / 1 ctn 43*41*36cm, 8/0Kgs			
335	A21.1610	2D/3D Manual Rotate Zoom Video Microscope	45° Lens, 360° Manual Rotatable, Zoom Body 0.6x-5.0x, C-Mount 0.5x, 2.0M 1080p HDMI Camera, USB For Mouse & U Disk, 12.5" LCD Monitor, Total Magnification 11-95x, View Field 24x13.5-2.88x1.62mm, W.D. 86mm, LED 4 Zones Control, Track Stand	1 pcs / 1 ctn 43*41*36cm, 8/0Kgs			
336	A21.1611	2D/3D Auto Rotate Zoom Video Microscope, Dual Light	45° Lens+LED Light, 360° Auto Rotate Control By HDMI Camera Mouse, Zoom Body 0.6x-5.0x, With 0.5x C-Mount, 8.0M 60fps HDMI Camera, USB For Mouse & U Disk, Total Magnification 14x-120x (15.6" LCD), View Field 24x13.5-2.88x1.62mm, W.D. 86mm, LED 4 Zones Control, Track Stand	1 pcs / 1 ctn 43*41*36cm, 8/0Kgs			
337	A21.1620	Calibration Free Motorized Zoom Video Microscope, 0.6x-5.0x	0.6x-5.0x, 2.0M 1080p HDMI Camera, USB Mouse+ U Disk, Total Magnification 14-120x (15.6" LCD), View Field 24x13.5-2.88x1.62mm, W.D. 86mm, LED 4 Zones Control, Track Stand Size 330x300mm With Coaxial	1 pcs / 1 ctn 43*41*36cm, 8/0Kgs			
338	A21.3601-B3	Monocular Zoom Video Microscope, 0.7-4.5x	0.7~4.5x, With CCD Mount, Lens 1x, W.D.100mm, XLB45+B3 Stand, No Light Source	1 pcs / 1 ctn 44.5*32.5*23cm, 4.5/3.7Kgs			

Step Zoom Stereo Microscope








339	A22.1001-13	Stereo Microscope	Binocular,1x/3x, EWF10x/23, Up & Bottom LED Light	1 pcs / 1 ctn 39*32*56cm, 7/6Kgs			
340	A22.1201-C	Stereo Microscope	Binocular 45°, WF10x, LWD Objective 1.0x, Boom Pole Stand, Incident Light Gooseneck LED	1 pcs / 1 ctn 41*20*43cm, 7/5Kgs			
341	A22.1233-TA	Stereo Microscope, Fixed Head	1x2x/1x3x/2x4x, Binocular,WF10x/20, Track Stand+Touch Switch, Up+Down LED, With 3 x AAA Rechargeable Batteries 1300mAh	2 pcs / 1 ctn 48*34*37cm, 6/4Kgs			
342	A22.1235-TA	Stereo Microscope, 360° Rotate Small Head	1x2x/1x3x/2x4x, Binocular,WF10x/20, Track Stand+Touch Switch, Up+Down LED, With 3 x AAA Rechargeable Batteries 1300mAh	2 pcs / 1 ctn 48*34*37cm, 6/4Kgs			
343	A22.1237-TA	Stereo Microscope, 360° Rotate Bigger Head	1x2x/1x3x/2x4x, Binocular,WF10x/20, Track Stand+Touch Switch, Up+Down LED, With 3 x AAA Rechargeable Batteries 1300mAh	2 pcs / 1 ctn 48*34*37cm, 6/4Kgs			
344	A22.1239-TA	Stereo Microscope, 360° Rotate Big Head With 3 Objectives 1x2x3x/1x2x4x	1x2x3x/1x2x4x, Binocular,WF10x/20, Track Stand+Touch Switch, Up+Down LED, With 3 x AAA Rechargeable Batteries 1300mAh	2 pcs / 1 ctn 48*34*37cm, 6/4Kgs			
345	A22.1307-13	Stereo Microscope	1x3x,Binocular,WF10x, Pole Stand, Up/Down LED Lamp	4 pcs / 1 ctn 67*40*41cm, 12/10Kgs			
346	A22.1308-13	Stereo Microscope	1x3x,Binocular,WF10x,Track Stand, Up/Down LED Lamp	4 pcs / 1 ctn 67*40*41cm, 12/10Kgs			
347	A22.1309-24	Stereo Microscope, 2x4x	2x4x, Binocular, WF10x/18mm, Track Stand, Forest Glass Plate Dia.95mm, Up Coaxial/Down LED Light, Working Distance 100mm	1 pcs / 1 ctn 34*41*25cm, 5/4Kgs			
348	A22.1504-C	Stereo Microscope	WF10x. Binocular, 2x/4x, Pole Stand, Upper & Bottom Light Source 12V10W Halogen	2 pcs / 1 ctn 31*29*36cm, 6.5/0Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
349	A22.1504-A	Stereo Microscope	WF10x. Binocular, 2x/4x, Pole Stand, No Light Source	4 pcs / 1 ctn 61*29*36cm, 9/0Kgs			
350	A22.1504-B	Stereo Microscope	WF10x. Binocular, 2x/4x, Pole Stand, Upper Light Source 12V10W Halogen	4 pcs / 1 ctn 61*29*36cm, 12.5/0Kgs			
351	A22.3660-B1	Step Zoom Stereo Microscope, 2x/4x	2x/4x, Binocular, WF10x/20, B1 Pole Stand + A1 Holder, Stand Size 200*255*22mm, Pole Height 240mm, Dia.32mm, No Light Source	1 pcs / 1 ctn 38.5*25*37cm, 4.4/3.5Kgs			
352	A22.3660-B2	Step Zoom Stereo Microscope, 2x/4x	2x/4x, Binocular, WF10x/20, B2 Track Stand + A1 Holder, Stand Size 200*255*60mm, Pole Height 240mm, Dia.32mm, Up/Bottom 12V15W Halogen Lamp	1 pcs / 1 ctn 50*26*37cm, 6.5/5Kgs			
353	A22.3660-B8L	Step Zoom Stereo Microscope, 2x/4x	2x/4x, Binocular, WF10x/20, B8L Pole LED Stand + A1 Holder, Stand Size 205*275*40mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED	1 pcs / 1 ctn 51*36*30cm, 5.4/4.2Kgs			
354	A22.3660-B9L	Step Zoom Stereo Microscope, 2x/4x	2x/4x, Binocular, WF10x/20, B9L Track LED Stand, Size 205*275*40mm, Track Height 300mm, Focusing Range 106mm, Dia.76mm, Up/Bottom 3W LED	1 pcs / 1 ctn 51*36*30cm, 5.5/4.3Kgs			
355	A22.3660N-B1	Step Zoom Stereo Microscope, 2x/4x	2x/4x, Binocular, WF10x/20, B1 Pole LED Stand + A1 Holder, Stand Size 200*255*22mm, Pole Height 240mm, Dia.32mm, No Light Source	1 pcs / 1 ctn 38.5*24*37cm, 5.4/4.4Kgs			
356	A22.3660N-B2	Step Zoom Stereo Microscope, 2x/4x	2x/4x, Binocular, WF10x/20, B2 Track Stand + A1 Holder, Stand Size 200*255*60mm, Pole Height 240mm, Dia.32mm, Up/Bottom 12V15W Halogen Lamp	1 pcs / 1 ctn 50*26*37cm, 5/4.4Kgs			
357	A22.3660N-B8L	Step Zoom Stereo Microscope, 2x/4x	2x/4x, Binocular, WF10x/20, B8L Pole LED Stand + A1 Holder, Stand Size 205*275*40mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED	1 pcs / 1 ctn 39*29*41cm, 5.5/4Kgs			
358	A22.3660N-B9L	Step Zoom Stereo Microscope, 2x/4x	2x/4x, Binocular, WF10x/20, B9L Track LED Stand, Size 205*275*40mm, Track Height 300mm, Focusing Range 106mm, Dia.76mm, Up/Bottom 3W LED	1 pcs / 1 ctn 39*29*41cm, 5.5/4Kgs			




No.	Cata. No.	Product	Specification	Package			Picture
359	A22.3670-B1	Stereo Microscope	2x/4x, Binocular, WF10x/23, B1 Pole Stand + A1 Holder, Stand Size 200*255*22mm, Pole Height 240mm, Dia.32mm, No Light Source	1 pcs / 1 ctn 38.5*24*37cm, 5.4/4.4Kgs			
360	A22.3670-B2L	Stereo Microscope	2x/4x, Binocular, WF10x/23, B2 Track Stand + A1 Holder, Stand Size 200*255*60mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED	1 pcs / 1 ctn 50*26*37cm, 5.5/4Kgs			
361	A22.3670-B8L	Stereo Microscope	2x/4x, Binocular, WF10x/23, B8L Pole LED Stand + A1 Holder, Stand Size 205*275*40mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED	1 pcs / 1 ctn 51*36*30cm, 6.4/5.2Kgs			
362	A22.3670-B9L	Stereo Microscope	2x/4x, Binocular, WF10x/23, B9L Track LED Stand, Size 205*275*40mm, Track Height 300mm, Focusing Range 106mm, Dia.76mm, Up/Bottom 3W LED	1 pcs / 1 ctn 51*36*30cm, 6.5/5.3Kgs			
Continous Zoom Stereo Microscope							
363	A23.3645-B1T	Zoom Stereo Microscope, 0.7~4.5x	0.7-4.5x, Trinocular, WF10x/20, B1 Pole Stand + A1 Holder, Stand Size 200*255*22mm, Pole Height 240mm, Dia.32mm, No Light Source, Incl. Eyepiece Adapter	2 pcs / 1 ctn 50*40*38cm, 11/9Kgs			
364	A23.3645-B2T	Zoom Stereo Microscope, 0.7~4.5x	0.7~4.5x, Trinocular, WF10x/20, B2 Pole Stand + A1 Holder, Stand Size 200*255*60mm, Pole Height 240mm, Dia.32mm, Up/Bottom 12V15W Halogen Lamp, Incl. Eyepiece Adapter	1 pcs / 1 ctn 50*26*37cm, 7.2/5.7Kgs			
365	A23.3645-B8LT	Zoom Stereo Microscope, 0.7~4.5x	0.7-4.5x, Trinocular, WF10x/20, B8L Pole LED Stand + A1 Holder, Stand Size 205*275*40mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED, Incl. Eyepiece Adapter	1 pcs / 1 ctn 51*36*30cm, 6.1/4.9Kgs			
366	A23.3645-B9LT	Zoom Stereo Microscope, 0.7~4.5x	0.7-4.5x, Trinocular, WF10x/20, B9L Track LED Stand, Size 205*275*40mm, Track Height 300mm, Focusing Range 106mm, Dia.76mm, Up/Bottom 3W LED, Incl. Eyepiece Adapter	1 pcs / 1 ctn 51*36*30cm, 6.2/5Kgs			
367	A23.3645N-B1T	Zoom Stereo Microscope, 0.7~4.5x	0.7-4.5x, Trinocular, WF10x/20, B1 Pole Stand + A1 Holder, Stand Size 200*255*22mm, Pole Height 240mm, Dia.32mm, No Light Source, Incl. Eyepiece Adapter	1 pcs / 1 ctn 38.5*24*37cm, 4.2/3.3Kgs			
368	A23.3645N-B2T	Zoom Stereo Microscope, 0.7~4.5x	0.7-4.5x, Trinocular, WF10x/20, B2 Track Stand + A1 Holder, Stand Size 200*255*60mm, Pole Height 240mm, Dia.32mm, UP/Bottom 12V15W Halogen Lamp, Incl. Eyepiece Adapter	1 pcs / 1 ctn 50*26*37cm, 6.3/4.8Kgs			











No.	Cata. No.	Product	Specification	Package			Picture
369	A23.3645N-B8LT	Zoom Stereo Microscope, 0.7~4.5x	0.7-4.5x, Trinocular, WF10x/20, B8L Pole LED Stand + A1 Holder, Stand Size 205*275*40mm, Pole Height 240mm, Dia.32mm, Up/Bottom 3W LED, Incl. Eyepiece Adapter	1 pcs / 1 ctn 51*36*30cm, 5.2/4.1Kgs			
370	A23.3645N-B9LT	Zoom Stereo Microscope, 0.7~4.5x	0.7-4.5x, Trinocular, WF10x/20, B9L Track LED Stand, Size 205*275*40mm, Track Height 300mm, Focusing Range 106mm, Dia.76mm, Up/Bottom 3W LED, Incl. Eyepiece Adapter	1 pcs / 1 ctn 51*36*30cm, 5.3/4.1Kgs			
371	A23.3667-B1T	Zoom Stereo Microscope, 0.67~4.5x	0.67-4.5x, Trinocular, WF10x/22, B1 Pole Stand + A1 Holder, Stand Size 200*255*22mm, Pole Height 240mm, Dia.32mm, No Light Source, Incl. Eyepiece Adapter	1 pcs / 1 ctn 38.5*24*37cm, 4.4/3.4Kgs			
372	A23.3667-B2T	Zoom Stereo Microscope, 0.67~4.5x	0.67-4.5x, Trinocular, WF10x/22, B2 Track Stand + A1 Holder, Stand Size 200*255*60mm, Pole Height 240mm, Dia.32mm, UP/Bottom 12V15W Halogen Lamp, Incl. Eyepiece Adapter	1 pcs / 1 ctn 50*26*37cm, 6.5/4.9Kgs			
373	A23.3667-B8LT	Zoom Stereo Microscope, 0.67~4.5x	0.67-4.5x, Trinocular, WF10x/22, B8L Pole LED Stand + A1 Holder, Stand Size 205*275*40mm, Pole Height 240mm, , Incl. Eyepiece Adapter Dia.32mm, Up/Bottom 3W LED	1 pcs / 1 ctn 51*36*30cm, 5.4/4.2Kgs			
374	A23.3667-B9LT	Zoom Stereo Microscope, 0.67~4.5x	0.67-4.5x, Trinocular, WF10x/22, B9L Track LED Stand, Size 205*275*40mm, Track Height 300mm, Focusing Range 106mm, Dia.76mm, Up/Bottom 3W LED, Incl. Eyepiece Adapter	1 pcs / 1 ctn 51*36*30cm, 5.5/4.3Kgs			
375	A23.3667L-B1T	Zoom Stereo Microscope, 0.67~4.5x	0.67-4.5x, Trinocular, WF10x/22, B1 Pole Stand + A1 Holder, Stand Size 200*255*22mm, Pole Height 240mm, Dia.32mm, No Light Source, Incl. Eyepiece Adapter	1 pcs / 1 ctn 39*29*41cm, 5.5/4Kgs			
376	A23.3667L-B2T	Zoom Stereo Microscope, 0.67~4.5x	0.67-4.5x, Trinocular, WF10x/22, B2 Track Stand + A1 Holder, Stand Size 200*255*60mm, Pole Height 240mm, Dia.32mm, UP/Bottom 12V15W Halogen Lamp, Incl. Eyepiece Adapter	1 pcs / 1 ctn 39*29*41cm, 5.5/4Kgs			
377	A23.3667L-B8LT	Zoom Stereo Microscope, 0.67~4.5x	0.67-4.5x, Trinocular, WF10x/22, B8L Pole LED Stand + A1 Holder, Stand Size 205*275*40mm, Pole Height 240mm, , Incl. Eyepiece Adapter Dia.32mm, Up/Bottom 3W LED	1 pcs / 1 ctn 39*29*41cm, 5.5/4Kgs			
378	A23.3667L-B9LT	Zoom Stereo Microscope, 0.67~4.5x	0.67-4.5x, Trinocular, WF10x/22, B9L Track LED Stand, Size 205*275*40mm, Track Height 300mm, Focusing Range 106mm, Dia.76mm, Up/Bottom 3W LED, Incl. Eyepiece Adapter	1 pcs / 1 ctn 39*29*41cm, 5.5/4Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
379	A23.0203-T	Zoom Stereo Microscope, 0.7-6.3x, 1:9	Trinocular, WF10x/22, 0.7-6.3x, 1:9, W.D.110mm, Track Stand, Reflected & Transmitted Illumination 3W LED, Stage 310x280mm	1 pcs / 1 ctn 59*43*42cm, 9/7Kgs			
380	A23.0301-31	Zoom Stereo Microscope, 0.7x-4.5x	0.7-5x, trinocular, big track stand, up & bottom light source.	0 pcs / 0 ctn 42.5*35.5*59.5cm, 9/0Kgs			
381	A23.0910-B1	Parallel Zoom Stereo Microscope, 0.63x-8x, 1:12.5, APO	0.63x-8x, Zoom Ratio 1:12.5, Parallel, 5-45° Titable Trinocular, WF10x/22mm, 1.0x APO Lens W.D.78mm, B1 Track Stand, Without Bottom Light, Upper LED Ring Light	1 pcs / 1 ctn 58.5*55*81cm, 14.8/14Kgs			
382	A23.0910-B2	Parallel Zoom Stereo Microscope, 0.63x-8x, 1:12.5, APO	0.63x-8x, Zoom Ratio 1:12.5, Parallel, 5-45° Titable Trinocular, WF10x/22mm, 1.0x APO Lens W.D.78mm, B2 Track Stand, Bottom Light Mirror 360° Rotatable, Incline Light 5W LED External Power Box + Single Optical Filter	1 pcs / 1 ctn 58.5*55*81cm, 16/14Kgs			
383	A23.0910-B10	Parallel Zoom Stereo Microscope, 0.63x-8x, 1:12.5, APO	0.63x-8x, Zoom Ratio 1:12.5, Parallel, 5-45° Titable Trinocular, WF10x/22mm, 1.0x APO Lens W.D.78mm, B10 Track Stand, Bottom Light 5W LED Color Temperature Adjustable, Upper LED Ring Light	1 pcs / 1 ctn 58.5*55*81cm, 14.8/14Kgs			
384	A23.1001-10	Parallel Zoom Stereo Microscope, Max 0.8-8x 1:10	0.8-8x, 1:10, Parallel, Binocular, EW10x/22, Track Stand, Transmit & Reflect LED Light Source	1 pcs / 1 ctn 37*32*67cm, 9/8Kgs			
385	A23.1001-F-10	Flourescent Zoom Stereo Microscope	0.8-8x, Binocular, EW10x/22, Transmit & Reflec LED, 100W Ultra hi-voltage spherical mercury lamp,	1 pcs / 1 ctn 37*31*66cm, 10/0Kgs			
386	A23.1001-T-10	Zoom Stereo Microscope, Tilttable Head, Max 0.8-8x, 1:10	0.8-8x, Binocular, EW10x/22, Transmit & Reflec LED,	1 pcs / 1 ctn 37*31*66cm, 10/0Kgs			
387	A23.1008	Parallel Zoom Flourescent Stereo Microscope, 0.75x-13.5x, 1:18, APO	0.75x-13.5x, Zoom Ratio 18:1, Trinocular 20° Incline, Light Split E100:P0/E0:P100, EW10X/23mm, Parallell Light Source, Plan APO 1.0x N.A.0.15, W.D.60mm, Zoom Range 60+90mm, LED Light Stand, 1x C-	1 pcs / 1 ctn 37*31*66cm, 10/0Kgs			
388	A23.1008-F	Parallel Zoom Flourescent Stereo Microscope, 0.75x-13.5x, 1:18, APO	0.75x-13.5x, Zoom Ratio 18:1, Trinocular 20° Incline, Light Split E100:P0/E0:P100, EW10X/23mm, Parallell Light Source, Plan APO 1.0x N.A.0.15, W.D.60mm, Zoom Range 60+90mm, LED Light Stand,	1 pcs / 1 ctn 37*31*66cm, 10/0Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
389	A23.1501-T1	Zoom Stereo Microscope, 0.8~5x	0.8~5x, Trinocular, WF10x/22mm, Track Big Fan Stand, No Light, 0.5x C Mount	1 pcs / 1 ctn 36*35*41cm, 8/7Kgs			
390	A23.1501-T2	Zoom Stereo Microscope, 0.8~5x	0.8~5x, Trinocular, WF10x/22mm, Track Big Fan Stand, Upper/Bottom LED, 0.5x C Mount	1 pcs / 1 ctn 36*35*41cm, 8/7Kgs			
391	A23.1501-T3	Zoom Stereo Microscope, 0.8~5x	0.8~5x, Trinocular, WF10x/22mm, Track Big Square Stand, No Light, 0.5x C Mount	1 pcs / 1 ctn 36*35*41cm, 8/7Kgs			
392	A23.1501-T4	Zoom Stereo Microscope, 0.8~5x	0.8~5x, Trinocular, WF10x/22mm, Track Big Square Stand, Upper/Bottom LED, 0.5x C Mount	1 pcs / 1 ctn 36*35*41cm, 8/7Kgs			
393	A23.2603-BL	Zoom Stereo Microscope, 0.7-4.5x, 1:6.5	0.7-4.5x, 1:6.5, Binocular, WF10x/22, Track Stand 33x30cm, Refelect & Transmit 3W LED Light, Safety Low Voltage Charger	1 pcs / 1 ctn 38*38*48cm, 6.8/5Kgs			
394	A23.2604-TL	Zoom Stereo Microscope, 0.68~4.7x, 1:6.8	0.68~4.7x, 1:6.8, Trinocular, WF10x/23mm, Track Stand 33x30cm, Up/Down 3W LED, Safety Low Voltage Charger	1 pcs / 1 ctn 38*38*48cm, 7.8/5.7Kgs			
395	A23.2605-TL	Zoom Stereo Microscope, 0.65~5.3x, 1:8.1	0.65~5.3x, 1:8.1, Trinocular, WF10x/24, Track Stand, Up/Down 3W LED 33x30cm, Safety Low Voltage Charger	1 pcs / 1 ctn 38*38*48cm, 7.8/5.7Kgs			








Naked Eye 3D Stereo Microscope

396	A22.0302-A	Naked Eye 3D Stereo Microscope	Magnification 2x4x6x, G Clamp Type, Built-in 3W LED	1 pcs / 1 ctn 67.5*43.5*69.5cm, 12.5/11Kgs			
397	A22.0302-B	Naked Eye 3D Stereo Microscope	Magnification 2x4x6x, Square Stand, Built-in 3W LED --Base Size 400x300x20mm (15.74x11.81x0.78 Inch) --Vertical Arm Height 500mm (19.68 Inch)	1 pcs / 1 ctn 67.5*43.5*69.5cm, 13/11Kgs			
398	A22.0302-C	Naked Eye 3D Stereo Microscope	Magnification 2x4x6x, Double Square Stand, Built-in 3W LED --Double Square Stand, Base Size 400x300x20mm (15.74x11.81x0.78 Inch) --Vertical Arm Height 500mm (19.68 Inch)	1 pcs / 1 ctn 67.5*43.5*69.5cm, 13/11Kgs			




No.	Cata. No.	Product	Specification	Package			Picture
399	A22.0302-D	Naked Eye 3D Stereo Microscope	Magnification 2x4x6x, Universal Base, Built-in 3W LED --Universal Base Total Height: 1800mm --Observation Head Rotatory Freely, Extension Range 850mm --Vertical Up/Down Range 580mm,	1 pcs / 1 ctn 67.5*43.5*69.5cm, 13/11Kgs			
400	A23.4501	3D Zoom Stereo Microscope, 0.7-5.6x	3D Touch-Screen 7", 0.7-5.6x, Ratio 1:8, Reflected & transmitted 3W LED, Brightness Adjustable, W.D. 110mm	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
Gem Microscope							
401	A24.0901-B	Gem Microscope, 0.67x-4.5x	0.67~4.5x, Binocular, WF10x/22, Upper Light 7W Fluorescent or 1W LED, Bottom Light 6V30W Halogen	1 pcs / 1 ctn 54*33*62cm, 12/8Kgs			
402	A24.1001-B	Gem Microscope, 0.8-5.0x	0.8-5.0x, Binocular, WF15x/16, Track Stand 0-38° Tilttable, Base 0-325° Rotatable, Working Distance 115mm, Reflect 9W Fluorescent, LED Point Light, Transmit 6V20W Halogen, With Dark Field Stage, Bright/Dark Field	1 pcs / 1 ctn 40*36*48cm, 6/0Kgs			
403	A24.1201-A24.1201	Jewellery Gem Microscope, 2x/4x	2x/4x, Binocular, WF10x/20mm, Track Stand, Reflect & Transmit Halogen 12V10W, Jewel Tweezer, Darkfield Working Stage	4 pcs / 1 ctn 77*31*40cm, 17/0Kgs			
404	A24.1202-A	Jewery Microscope, 1x-4x, 0.75x-3.6x	1~4x, Binocular, SWF10x/22mm, Pole Stand, Reflect 12V10W Halogen + 10W Sunshine Fluorescent Light, Transmit 12V10W Halogen, Jewellery Tweezer, Dark Field Working Stage	2 pcs / 1 ctn 50*38*40cm, 13/11Kgs			
405	A24.1203-T	Jewery Microscope, 0.75~5x	0.75~5x, Trinocular, SWF10x/23mm, Track Stand Tilttable, Reflect 12V10W Halogen + 10W Sunshine Fluorescent Light, Transmit 12V10W Halogen, Jewellery Tweezer, Dark Field Working Stage	1 pcs / 1 ctn 50*38*48cm, 13/10Kgs			
406	A24.1204	Jewellery Gem Microscope, 2x	Monocular, 2x, WF10x, Transmit LED, With Dark Field Stage, With Jewellery Twizzer	10 pcs / 1 ctn 70*39*29cm, 6/0Kgs			
Digital Biological Microscope							
407	A31.5121-M	Digital Biological Microscope, 5.0M	5.0M, Teaching Head, WF10x, WF16x, Quarduple, Ach. 4x10x40x100x, Mechanical Stage 115x125mm, Coaxial Focusing, LED, With 10 Prepared Slides, 0.3x C-Mount, With Experiment Kit, Pack In Hand Carry Case+Carton	2 pcs / 1 ctn 0*0*0cm, 7.5/0Kgs			
408	A31.5121-T	Digital Biological Microscope, 5.0M	5.0M, Trinocular Head, WF10x, WF16x, Quarduple, Ach. 4x10x40x100x, Mechanical Stage 115x125mm, Coaxial Focusing, LED, With 10 Prepared Slides, 0.3x C-Mount, With Experiment Kit, Pack In Hand Carry Case+Carton	2 pcs / 1 ctn 0*0*0cm, 7.5/0Kgs			











No.	Cata. No.	Product	Specification	Package			Picture
409	A31.6603-5.0M	Digital Biological Microscope, 5.0M	5.0M, 1/2" CMOS, Resolution 2952x1944, 15-30fps, USB2.0, Built-in 0.5x C-Mount Adapter, With Software, Binocular, WF10x/20, Quadruple, Infinity E-Plan 4x10x40x100x, Mechanical Stage 145X140/76x52mm, Coaxial, 3W LED	1 pcs / 1 ctn 37*26*61cm, 7.5/6Kgs			
410	A31.6603-5.0MH	Digital Biological Microscope, 5.0M	5.0M, 1/2" CMOS, Resolution 2952x1944, 15-30fps, USB2.0, Built-in 0.5x C-Mount Adapter, With Software, With Pad Holder, Binocular, WF10x/20, Quadruple, Infinity E-Plan 4x10x40x100x, Mechanical Stage 145X140/76x52mm, Coaxial, 3W LED	1 pcs / 1 ctn 37*26*61cm, 7.5/6Kgs			
411	A31.1006-A	Digital Microscope, 1.3M	1.3M, USB2.0, Monocular, WF10x/18mm, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 115x125mm/ 70x20mm, Abbe NA1.25, Coaxial, 1W LED	1 pcs / 1 ctn 28*22*39cm, 3.5/3Kgs			
412	A31.1006-B	Digital Microscope, 1.3M	1.3M, USB2.0, Binocular, WF10x/18mm, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 115x125mm/ 70x20mm, Abbe NA1.25, Coaxial, 1W LED	1 pcs / 1 ctn 28*22*39cm, 3.7/3Kgs			
413	A31.1007-3.0M	Digital Microscope	3.0M, USB2.0, Binocular, WF10x/18mm, Quadruple, Achromatic 4x10x40x100x, Mechanical Stage 140x140mm/ 75x50mm, Abbe N.A.1.20, Coaxial, Halogen 6V20W	1 pcs / 1 ctn 37*29*47cm, 13/0Kgs			
414	A31.1030	Digital Laboratory Microscope, Coded Nosepiece & LCD, 8.0M	8.0M Binocular, WF10x/20, Coded Quadruple, Infinity Plan NIS45 4x10x40x100x, Rackless Stage 180x130/ 74x30mm, N.A.1.25 Condenser With Objective Position Scales, Coaxial Focusing, 3W LED	1 pcs / 1 ctn 47*32*59cm, 12.1/11Kgs			
415	A31.1062-5.0M	Digital WIFI Laboratory Microscope, 5.0M	Digital Binocular 5.0M, USB2.0, WF10x/22, Coding Quintuple, Infinity Plan NIS60 4x10x40x100x, Rackless Stage 230x150/ 78x54mm, N.A.1.25 Condenser With Slot, Coaxial, 3W LED, LCD Screen Base, Brightness Auto	1 pcs / 1 ctn 47*32*59cm, 12.1/11Kgs			
416	A31.0911	Digital Biological Microscope, 4.0M, HDMI+USB	Digital Binocular Head, Built-in 4.0M 1/3" CMOS Camera, Output HDMI++MiniUSB+SD Card, Interpupillary Distance 50-75mm, ±5 Diopter Adjustable, Light Split 8:2, Eyepiece PL10X/20mm, Infinity Plan	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
417	A31.0921-2M	WIFI Digital Biological Microscope, 2.0M, Android	Monocular, Built-in 2M WIFI Camera, With 6-13" Pad Holder, PL10X/18mm, Triple, Achromatic, 4x10x40, Plain Stage, Condenser 0.65+ Disk Diaphragm, LED. WIFI Support Max 2 Devices, Support Android	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
418	A31.0923-M2M	WIFI Digital Biological Microscope, 2.0M, Android	Monocular, Built-in 2M WIFI Camera, With 6-13" Pad Holder, PL10X/18mm, Quadruple, Achromatic, 4x10x40x100x, Mechanical Stage 115x125mm, Abbe Condenser 1.25+ Iris Diaphragm, LED. WIFI Support Max 2 Devices, Support	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
419	A31.0923-B2M	WiFi Digital Biological Microscope, 2.0M, Android	Binocular, Built-in 2M WiFi Camera, With 6-13" Pad Holder, PL10X/18mm, Quarduple, Achromatic, 4x10x40x100x, Mechanical Stage 115x125mm, Abbe Condenser 1.25+ Iris Diaphragm, LED. WiFi Support Max 2 Devices, Support	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
420	A31.0925-5M	WiFi Digital Biological Microscope, 5.0M, Win, iMac, iOS, Android	Binocular, Built-in 5M WiFi Camera, With 6-13" Pad Holder, PL10X/18mm, Quarduple, Plan 4x10x40x100x, Mechanical Stage 132x140mm, Abbe Condenser 1.25+ Iris Diaphragm, LED. WiFi Support Max 6 Devices	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
421	A31.0927-5M	WiFi Digital Biological Microscope, 5.0M, 12M, Win, iMac, iOS, Android	Monocular, Built-in 5M WiFi Camera, With 6-13" Pad Holder, PL10X/20mm, Quarduple, Plan 4x10x40x100x, Mechanical Stage 150x140mm, Abbe Condenser 1.25+ Iris Diaphragm + Slot for Dark Field and PH, LED. WiFi	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
422	A31.0927-12M	WiFi Digital Biological Microscope, 5.0M, 12M, Win, iMac, iOS, Android	Binocular, Built-in 2M WiFi Camera, With 6-13" Pad Holder, PL10X/18mm, Quarduple, Achromatic, 4x10x40x100x, Mechanical Stage 115x125mm, Abbe Condenser 1.25+ Iris Diaphragm, LED. WiFi Support Max 2 Devices	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
423	A31.3701-14M	Digital Biological & Stereo Dual Light Cordless Microscope	14.0M, 30fps/1080p, USB3.0, Binocular, WF10x/18, Quintple, Infinity Achromatic 4x10x40x100x + Stereo 1x, Double Layer Mechanical Stage 140x130/76x36mm, Coaxial Focusing, LED Up/Side Light Source	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
Digital Stereo Microscope							
424	A32.1201	Digital Stereo Microscope, 2x/4x, 1.3M	2x4x, 1.3M, Binocular, WF10x/20mm, Track Stand, Coarse Focusing, Reflect & Transmit Halogen 12V10W, Brightness Adjustable	4 pcs / 1 ctn 68*31*37cm, 16/0Kgs			
425	A32.1239-3.0M	Digital Stereo Microscope, 1.3M	3.0M, 1x2x4x, Binocular, WF10x/20mm, Track Stand Dial Swich, Glass Plate + B/W Board Dia.6cm, Up/Down LED, With 3 x AAA Rechargeable Batteries 1300mAh	2 pcs / 1 ctn 48*34*37cm, 6/4Kgs			
426	A32.1501	Digital Stereo Microscope	2x/4x, 1.3M CMOS, USB2.0	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
427	A32.3610	Digital Stereo Microscope	0.7~4.5x, 2M, LED Ring Light, Track Stand, Output HDMI/ U Disk/ TF, USB Mouse Control , Measure Function, Total Magnification 17.2~110.5x, Optional Polarzing, Auxiliary Lens 0.5x, 2.0x	1 pcs / 1 ctn 38.5*24*37cm, 4.9/4Kgs			
428	A32.3645-B2	Digital Stereo Microscope, 0.7x~4.5x, 5.0M	5.0M, Binocular, WF10x/20mm, 7x-45xZoom Body, B2 Pole Stand + A1 Holder, Stand Size 200*255*60mm, Pole Height 240mm, Dia.32mm, Up/Bottom 12V15W Halogen Lamp	1 pcs / 1 ctn 38.5*24*37cm, 4.9/4Kgs			





No.	Cata. No.	Product	Specification	Package			Picture
429	A32.3645-B8LD	Digital Stereo Microscope, 0.7x-4.5x, 5.0M	5.0M, Binocular, WF10x/20mm, 7x-45xZoom Body, B8LD Pole Stand + A1 Holder, Stand Size 205*275*40mm, Pole Height 240mm, Dia.32mm, Up 30 LED & Bottom 92 LED	1 pcs / 1 ctn 38.5*24*37cm, 4.9/4Kgs			
430	A32.3645-B9LD	Digital Stereo Microscope, 0.7x-4.5x, 5.0M	5.0M, Binocular, WF10x/20mm, 7x-45xZoom Body, B9LD Track LED Stand, Size 205*275*40mm, Track Height 300mm, Focusing Range 106mm, Dia.76mm, Up 30 LED & Bottom 92 LED	1 pcs / 1 ctn 38.5*24*37cm, 4.9/4Kgs			
431	A32.0902	Digital Stereo Microscope, 0.8x-3.5x, 5.0M	Binocular, 60°Inclined, WF10x/22, Zoom Objective 0.8x-3.5x,W.D.95mm, Track Stand,Up & Bottom 5W LED, Separate Light Control, 5.0M CMOS, Output HDMI, Video, USB2.0,With SD Solt, With Photo Capture Button	1 pcs / 1 ctn 46*26*42cm, 7/0Kgs			
432	A32.0907	Auto Focus Digital Stereo Microscope, 0.7x-4.5x, 2.0M, HDMI Measure	Auto Focus, 2.0M 1/2.8" CMOS Resolution 1920*1080@60fps, HDMI Output, W.D.100mm, USB Disk Storage, Optical Lens 0.7-4.5x, Digital Zoom 1-6x, Mouse Control Measure Under HDMI Interface, LED Ring Light	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
433	A32.0914	Auto Focus Digital Stereo Microscope, 1x-14x, 2.0M, HDMI Measure	Auto Focus, 2.0M 1/2.8" CMOS Resolution 1920*1080@60fps, HDMI Output, W.D.100mm, TF Card Storage, Optical Lens 1-14x, Digital Zoom 1-6x, Mouse Control Measure Under HDMI Interface, LED Ring Light 5500K??,	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
434	A32.0920	Digital Stereo Microscope	0.7-5.0x, 4.0M, 1/3" CMOS, HDR Image, OSD Manual, 16G Memory, Output HDMI+USB2.0, 1.0x Lens W.D.105mm, LED 56pcs Brightness Adjustable, Track Stand Lifting 225mm, DC12V	1 pcs / 1 ctn 49*35*43cm, 6/5Kgs			
435	A32.6401-12M	Digital Stereo Microscope, 0.7x-4.5x, 12M, HDMI+SD+WIFI(Optional)		1 pcs / 1 ctn 44*31*18cm, 4.5/4Kgs			

LCD Digital Biological Microscope






436	A33.5102	7" LCD Biological Digital Microscope, 2.0M	Mono Straight Head, WF10x, WF16x, Triple, Ach. 4x10x40x, Mechanical Stage, Separate Coarse & Fine Focusing, 7" LCD Resolution 800x480, 0.6M 1/3" CMOS, 976x592 Pixels, AV Output, Power by 5V DC USB Adapter,	1 pcs / 1 ctn 30*19.5*39.5cm , 3.9/3Kgs			 NEW
437	A33.5111	7" LCD Biological Digital Microscope, 2.0M, USB+AV+TF Card	7" IPS LCD, Resolution 1024x600, 2.0M 1/2.75" CMOS, 1920*1080 Pixels, 30 Fps, Max Support 32G TF Card, AV Output, Digital Zoom 2x/4x, USB2.0 Output to PC, With Software StrangeView, Chinese/Eng/Rus/Urk,	0 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
438	A33.5112	9.7" LCD Biological Digital Microscope, 5.0M, HDMI+TF Card	5.0M, 9.7" LCD, Resolution 2048*1536, Android 5.1, CPU 4 Core 1.8G, RAM 2G ROM 8G, HDMI+TF Card Output, Touch Screen USB Mouse Control, Installed Measure Software, Chinese/English Manual,	1 pcs / 1 ctn 30*40*50cm, 6/4Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
439	A33.5121-TH	7" LCD Biological + USB Portable Dual Lens Digital Microscope, 2.0M	Trinocular, WF10x, WF16x, Quaduple, Ach. 4x10x40x100x, Mechanical Stage 115x125mm, Coaxial Focusing, LED, With Prepared Slides, 0.3x C-Mount, With Experiment Kit, Pack In Hand Carry Case+Carton	1 pcs / 1 ctn 30*20*40cm, 5/4.5Kgs			 HOT
440	A33.5121-MH	7" LCD Biological + USB Portable Dual Lens Digital Microscope, 2.0M	Teaching Head, WF10x, WF16x, Quaduple, Ach. 4x10x40x100x, Mechanical Stage 115x125mm, Coaxial Focusing, LED, With Prepared Slides, 0.3x C-Mount, With Experiment Kit, Pack In Hand Carry Case+Carton	1 pcs / 1 ctn 30*20*40cm, 5/4.5Kgs			 HOT
441	A33.1009	9" LCD Digital Biological Microscope, 5.0M HDMI+USB	Sliding Trinocular Black, WF10x16x, Quaduple, Achromatic 4x10x40x100x, Double Layer Mechanical, Coaxial Coarse&Fine, Abbe 1.25, LED, Foam +Carton Packing, 1.0x C-Mount, 5.0M, 9" LCD, 1/2.8" CMOS, Screen	1 pcs / 1 ctn 33*28*55cm, 9/7Kgs			 HOT
442	A33.1501	3.5" LCD Digital Microscope, 2.0M	3.5" LCD, 2.0M, Monocular, WF10x, 16x, Triple, Achromatic 4x10x40x, Vernial Scale Stage, Up/Bottom LED, With Accessories: Plastic Microtom, Hatching Box, Tweezer, Sea salt 1 Bottle, Yeast 1 Bottle,	2 pcs / 1 ctn 42*36*28cm, 7/6Kgs			
443	A33.1502	10.5" LCD Touch Screen Digital Microscope, 8.0M, Android Pad	10.5" LCD, 8.0M, Touch Screen Android Pad, Installed Measuring Software, Trinocular, WF10x/20, Quaduple, Achromatic HC 4x,10x,40x100x, Double Layer Mechanical Graphite 143*132mm,	1 pcs / 1 ctn 53*32.5*53.5cm 11.5/9Kgs			 HOT
444	A33.3712-16M	10.1" LCD Digital Biological Microscope	10.1" LCD, 16.0M, 30fps/1080p, Output TF Card + HDMI, Binocular, WF10x/18, Quaduple, Achromatic 4x10x40x100x + Stereo 1x, Double Layer Mechanical Stage 140x130/76x36mm, Coaxial Focusing, LED Up/Side Light Source	1 pcs / 1 ctn 41*32*51.5cm, 10/8Kgs			
445	A33.3713-16M	10.1" LCD Digital Biological Microscope, Internet Version	10.1" LCD, 16.0M, 30fps/1080p, Output Internet/WIFI/USB2.0 + Built-in 8G + U Disk, Binocular, WF10x/18, Quinduple, Achromatic 4x10x40x100x + Stereo 1x, Double Layer Mechanical Stage 140x130/76x36mm, Coaxial Focusing,	1 pcs / 1 ctn 41*32*51.5cm, 10/8Kgs			
446	A33.3714-14M	10.1" LCD Digital Biological Microscope, Win10 Pad +Internet Version	10.1" LCD, 16.0M, 30fps/1080p, Output HDMI+VGA, Built-in PC 4G Memory, 64G Hard Disk, Win10, Binocular, WF10x/18, Quinduple, Achromatic 4x10x40x100x + Stereo 1x, Double Layer Mechanical Stage	1 pcs / 1 ctn 41*32*51.5cm, 10/8Kgs			
447	A33.3733-16M	10.1" LCD Digital Fluorescent Microscope, Ethernet Version	10.1" LCD, 16.0M, 30fps/1080p, Output Internet/WIFI/USB2.0/Built-in 8G/U Disk, Quinduple, Infinity Plan 4x10x40x100x + Stereo 1x, Double Layer Mechanical Stage 140x130/76x36mm, Coaxial Focusing,	1 pcs / 1 ctn 33.5*32*58cm, 10/8Kgs			
LCD Digital Stereo Microscope							
448	A36.5001	3.5"LCD Digital Microscope, 500x,5.0M	3.5" LCD, 5.0M, 20~500x, 8 LED	6 pcs / 1 ctn 44*42*35cm, 9/8Kgs			


No.	Cata. No.	Product	Specification	Package			Picture
449	A36.5006	4.3" LCD Digital Microscope, 8.0M, 600x	4.3"LCD, 600x, Capture Resolution 8M, 5M, 3M,1.3M, VGA, Manual Focus Range 5~80mm, Storage Micro-SD卡 (Included), Measuring Function When USB output, LED Adjustable, For Windows 7/8/10, Mac, Portable Capture	1 pcs / 1 ctn 21*16.5*12.5cm 0.7/0Kgs			
450	A36.5049	4" LCD Digital Microscope With Polarizer, 10x-300x/1200x, 12M	4" LCD Resolution 800*400 IPS, Photo 12M, Video 1080p, 4" LCD Magnification 3x-44x, Output HDMI, Optical Magnification 10~300x (On 12M Photo), Digital Zoom Wheel 4x, MicroSD Slot (Max Support 32G),	20 pcs / 1 ctn 42.5*40.5*35.5cm, 14/12Kgs			 NEW
451	A36.5058	5" LCD Digital Microscope With Polarizer, 10x-300x/1200x, 12M	5" LCD Resolution 960*540 IPS, Photo12M, Video 1080p, 5" LCD Magnification 3x~48x, Optical Magnification 10~300x (On 12M Photo), Digital Zoom Wheel 4x, MicroSD Slot (Max Support 32G), Output HDMI, AV,	4 pcs / 1 ctn 42.5*40*33cm, 6.3/5Kgs			 NEW
452	A36.5501	7" LCD Stereo + Biological Digital Microscope, 20M, 500x/2000x	7.0" LCD, Main Body 2.0M Camera, 5x-53x Zoom, Portable USB2.0 1.3M, 1.5x-25x Zoom, One Key Switch Between Main/Portable, Support Max 32G TF Card, USB Output To Win/Mac, With StrangeView Measuring Software, LED	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			
453	A36.0920	LCD Digital Stereo Microscope, 10.6"	10.1" IPS Touch Screen, 5~15" Adjustable, 0.7~5.0x, 4.0M, 1/3" CMOS, HDR Image, OSD Manual, 16G Memory, Output HDMI+USB2.0, 1.0x Lens W.D.105mm, LED 56pcs Brightness Adjustable, Track Stand	1 pcs / 1 ctn 49*35*43cm, 6/5Kgs			
454	A36.1210	9" LCD Digital Measure Stereo Microscope, 46x	1x-4x, 3.0M, 9" LCD, Output USB2.0+SD Card, Up LED Ring Light, Bottom LED Light, Touch Control Total Mag. 11.5 ~ 46x, Mouse Control Measure Function	1 pcs / 1 ctn 41*39*42cm, 6/0Kgs			
455	A36.6401-12MH	10.6" LCD Digital Stereo Microscope, 14x-90x, 12M	10.6" LCD, 14x-90x, 12M, Video 1080p/60fps, 0.7x-4.5x, Up & Side LED Brightness Adjustable, Black Pole Stand 38*25cm, DC12V1A, SD Slot, HDMI	1 pcs / 1 ctn 38*40*44cm, 5.8/5Kgs			
456	A36.3601	11.6" LCD Digital Measure Stereo Microscope, 110.5x 5.0M	0.7~4.5x, 5.0M, 11.6" LCD 1080p, Output HDMI/ USB Disk/TF Card, Up LED Ring Light, Track Stand Base 260x320x15mm, Total Magnification 17.2~110.5x, USB Mouse Control, Measure Function,	1 pcs / 1 ctn 43*40*24cm, 6/5Kgs			
457	A36.2201-A36.2201	11.6" LCD Auto Focus Microscope, 1.296x~25.938x	2.0M, 1/2.8", 1080P CMOS, HDMI Output, USB Mouse Control Measure, SD Card Storage, Auto Focus, LED Ring Light, Built-in Zoom Lens 0.028-0.56x, W.D.195mm, Incl. 76mm Bracket, 11.6" LCD Screen, Total Mag.	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
458	A36.4970	15.6" LCD HDMI+USB Digital Microscope, 4.5x-49.2x, 2.0M	4.5x-49.2x, 15.6" LCD Full HD 1080p, 2.0M, 1/2.8" CMOS, Back Focus, TF Card, HDMI Output, USB Output With Measuring Software, W.D.44-230mm, Depth of View 0.8-6.5mm, View Field 7x4mm-45.7x25.8mm, Universal Hinge	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
459	A36.4971	12.5" LCD HDMI+USB Digital Microscope, 4.5x-49.2x, 2.0M	4.5x-49.2x, 15.6" LCD Full HD 1080p, 2.0M, 1/2.8" CMOS, Back Focus, TF Card, HDMI Output, USB Output With Measuring Software, W.D.44-230mm, Depth of View 0.8-6.5mm, View Field 7x4mm-45.7x25.8mm, Stainless Pole	0 pcs / 0 ctn 36*21*29cm, 3.2/0Kgs			
460	A36.5101	7" LCD Mono Zoom+ USB Portable Dual Lens Digital Microscope, 2.0M+1.3M	7.0" LCD, Main Body 2.0M Camera, 5x-53x Zoom, Portable USB2.0 1.3M, 1.5x-25x Zoom, One Key Switch Between Main/Portable, Support Max 32G TF Card, USB Output To Win/Mac, With StrangeView Measuring Software, LED	1 pcs / 1 ctn 33*25.5*25cm, 1.8/0Kgs			 NEW
461	A36.5102	7" LCD Stereo + USB Portable Dual Lens Digital Microscope, 2.0M+1.3M	7.0" LCD, 2.0M, 8.17x-53x zoom, handheld USB2.0 1.3M, 1.5x-25x zoom, one-key switch host/handheld, support up to 32G TF card photo and video, USB output to computer, including measurement Software StrangeView	1 pcs / 1 ctn 28*38*48cm, 6.5/5Kgs			 HOT
462	A36.6401-B8L-16MX	LCD Digital Stereo Microscope	0.7~4.5x, 16M, 10.6" LCD, 1080p, USB3.0, HDMI, LED Brightness Adjustable, Pole Stand B8L, DC12V, SD Slot, HDMI Output	1 pcs / 1 ctn 34*44*44cm, 4/0Kgs			

Portable USB/WIFI Digital Microscope

463	A34.4191	USB Digital Microscope, 50x~1000x	VGA, 640*480, 50x~1000x, 8LEDs, USB2.0	54 pcs / 1 ctn 51.5*27*37cm, 15/14Kgs			
464	A34.5001	USB Digital Microscope, 300x, 5.0M	5.0M, 20~300x, 8 LED, USB 2.0, With Stand	12 pcs / 1 ctn 45*36*43cm, 12/11Kgs			
465	A34.4970	HDMI+USB Digital Microscope, 5.3x-39.4x, 2.0M	5.3x-39.4x, 2.0M, 1/2.8" CMOS, Back Focus, TF Card, HDMI Output, USB Output With Measuring Software, W.D.44-230mm, Depth of View 0.8-6.5mm, View Field 7x4mm-45.7x25.8mm, Universal Hinge Stand	1 pcs / 1 ctn 37*28*18cm, 2.1/1.5Kgs			
466	A34.4199	WIFI+USB Miroscope, 1000x, 2.0M, iPad/Android/Win/Mac	2.0M, 50x~1000x, 8 LEDs, WIFI 10m Range, WIFI Support Android, IOS, USB Support Win, Mac, Output MJPEG/AVI, Batteries For 3 Hours Working, Black/Grey/Blue Optional	54 pcs / 1 ctn 51.5*27*37cm, 15/14Kgs			
467	A34.5017	WIFI Digital Miroscope, 200x, 2.0M, iPad/Android	2.0M, 1080p, 10x~200x, 6 LED, Micro SD Card, With Stand, Focus 0-50mm, iPhone/iPad/Android, App: Tiny Capture, Rechargeable Battery	6 pcs / 1 ctn 44*42*35cm, 9/8Kgs			




Digital Camera For Microscope

468	A59.4902	HDMI+USB Digital Camera 3.5M, Manual Focus	3.5M, 1080p 30 fps, Output USB/HDMI+SD Card, Manual Focus, C-Mount	40 pcs / 1 ctn 55*31*39cm, 18.5/17Kgs			
-----	----------	--	--	---	--	--	---

No.	Cata. No.	Product	Specification	Package			Picture
469	A59.4910	USB2.0 CMOS Digital Camera, 5.0M, Fix USB Cable	5.0M, 1/2.5" CMOS, USB2.0, C-Mount, Fixed USB Cable	40 pcs / 1 ctn 54*38*31cm, 11.5/10Kgs			
470	A59.4910-M05	USB2.0 CMOS Digital Camera, 5.0M, Fix USB Cable	5.0M, 1/2.5" CMOS, USB2.0, C-Mount, Fixed USB Cable, With 0.5x Eyepiece Tube Adapter Ring Set (23.2mm, 30.0mm, 30.5mm)	40 pcs / 1 ctn 54*38*31cm, 16/15Kgs			
471	A59.4920-W	5G WIFI+USB Digital Camera, 5.0M, 30Fps, Win,Android,iOS	5.0M, 30fps, 1/2.5" CMOS, 5G WIFI+USB2.0, Biological/Stereo Modes, LED/Halogen Modes, Win/Android/iOS	12 pcs / 1 ctn 52*44*39cm, 12/10Kgs			
472	A59.4972	5G WIFI+4K HDMI+USB+WAN Digital Camera, 12.0M, C-Mount, USB/WIFI Measuring	12.0M, 1/2.3" CMOS, 5G WIFI+4K HDMI+USB2.0+Ethernet Full Output Port, 4000x3000, 60fps@1928*1080p, USB+WIFI Measure Software Support Win/Android/iOS, View Mode Biological/Stereo, LED/Halogen	0 pcs / 0 ctn 30.5*16.5*10.5cm, 0.9/0Kgs			
473	A59.4974	5G WIFI+4K HDMI+USB Disk Digital Camera, 12.0M, C-Mount, HDMI Measuring	12.0M, 1/2.3" CMOS, 5G WIFI+4K HDMI+USB3.0 Disk, HDMI Built-in Operating System+Office, USB2.0 Mouse Keyboard Control, Built-in 32G Memory, WIFI Output Support Win/iOS/Android, 30FPS@4000*30000, 30FPS@1928*1080, View Mode	0 pcs / 0 ctn 30.5*16.5*10.5cm, 0.9/0Kgs			
474	A59.2211-10MPA	USB3.0 CMOS Digital Camera, U3CMOS	10M, 1/2.3" CMOS, USB3.0, C-Mount	20 pcs / 1 ctn 52*32*33cm, 15/14Kgs			
475	A59.2213-12MPA	USB3.0 CMOS Digital Camera, High Fps, Fluorescent Image, 1.5M~45M	12.0M, 1/1.7" CMOS, USB3.0, C-Mount, 50FPS@2048x1080	20 pcs / 1 ctn 50*47*35cm, 18/0Kgs			
476	A59.2225-20MPA	USB3.0 CMOS Cooled Digital Camera, Fluorescent Image, GS Shutter, 1.7M~45M	20M, 1" CMOS, USB3.0, 30@1824x1216, Cooling 42 Degree Temperature Drop, Anti-fog Structure, 1 Hour Exposure, C-Mount	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
477	A59.3515	USB3.0 CMOS Digital Camera, 16M	USB3.0, CMOS, 16MP, 4632x3488@8, DC5V	1 pcs / 1 ctn 12*23*28cm, 2/1.5Kgs			
478	A59.3830-16M	USB3.0 CMOS Digital Camera, 16M	16MP, CMOS 1/2.33", USB3.0, Windows/Mac/Linux	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
479	A59.2226-6MPA	USB3.0 CCD Cooled Digital Camera, Fluorescent Image	6M, 1" CCD, USB3.0, 14@2748x1092, Cooling 50 Degree Temperature Drop, Anti-fog Structure, 1 Hour Exposure, C-Mount	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
480	A59.3516	USB3.0 CCD Digital Camera, 2.8M	USB3.0, Sony CCD, 280万像素, 1940x1460@25, DC5V	1 pcs / 1 ctn 12*23*28cm, 2/1.5Kgs			
481	A59.2205-5.1MPA	USB2.0 CMOS Digital Eyepiece Camera, 0.92M~12M	5.0M, 1/2.5", USB2.0, 30@2592x1944, Dia. 23.2mm Eyepiece Tube, With Adapter Ring 30mm, 30.5mm	50 pcs / 1 ctn 30*30*50cm, 12/0Kgs			
482	A59.2205-E-2.0M	USB2.0 CMOS Digital Eyepiece Camera	2.0M, 1/3.2", USB2.0, Dia. 23.2mm Eyepiece Tube, With Adapter Ring 30mm, 30.5mm	50 pcs / 1 ctn 30*30*50cm, 12/0Kgs			
483	A59.2208-5.1MPA	USB2.0 CMOS Digital Camera, 1.3M~14M	5.1M, 1/2.5", USB2.0, C-Mount	20 pcs / 1 ctn 50*47*35cm, 18/0Kgs			
484	A59.2208-14MPA	USB2.0 CMOS Digital Camera, 1.3M~14M	14M, 1/2.3", USB2.0, C-Mount	20 pcs / 1 ctn 50*47*35cm, 18/0Kgs			
485	A59.3808-12M	USB2.0 CMOS Digital Camera	12MP, Panasonic 1/2.33", 4000x3000, 12fps, USB2.0, PC Win7/8/10	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
486	A59.3808-5M	USB2.0 CMOS Digital Camera	5.0MP, Sony 1/2.8", 2560x1920, 22fps, USB2.0, PC Win7/8/10	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
487	A59.4920	USB2.0 CMOS Digital Camera, 5.0M 30Fps	5.0M, 30fps, 1/2.5" CMOS, USB2.0, Biological/Stereo Modes, LED/Halogen Modes, Win/Android/IOS	40 pcs / 1 ctn 54*38*31cm, 11.5/10Kgs			
488	A59.5101-2.0MA	USB2.0 CMOS Digital Eyepiece Camera, 1.3M/2.0M	2.0M, USB2.0, Dia. 23.2mm Eyepiece Tube, USBCable, Software For Win&Mac	60 pcs / 1 ctn 61*26*49.5cm, 12/10Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
489	A59.5102-5.0M	USB2.0 CMOS Eyepiece Camera, 5.0M, Eyepiece Type	5.0M, 1/2.5" CMOS, Eyepiece Tube Adapter, Photo Format MJPG, Video Format YUY2, 16fps/2592*1944, 30fps/640*480, USB Cable, Software For Windows & Mac.	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
490	A59.5103-5.0M	USB2.0 CMOS Eyepiece Camera, 5.0M, C-Mount	5.0M, 1/2.5" CMOS, C-Mount, Photo Format MJPG, Video Format YUY2, 16fps/2592*1944, 30fps/640*480, USB Cable, Software For Win&Mac	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
491	A59.2221-2M	HDMI Digital Camera, 720p, Mouse Measure	720P, 2M, 1/2.8", 30@1280x720, HDMI Control Panel By Mouse, Small Size (50*50*54mm), C-mount, SD Card Slot	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
492	A59.2222-2MPA	HDMI+USB+SD Digital Camera, 2.0M~8.0M, Mouse Measure	1080P, 2.0M, 1/2.8", HDMI+USB+SD Card, FPS 60@1920x1080 (HDMI), 50@1920x1080 (USB), HDMI + USB2.0 Output At Same Time, HDMI Control Panel For Mouse Measure, C-Mount	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
493	A59.2241-8MPA	4K HDMI+USB+SD Digital Camera, 8.0M, Mouse Measure	4K, 8.0M, 1/1.8", HDMI+USB+SD Card, FPS 30@3840x2160 (HDMI), 20@3840x2160 (USB), HDMI UI + USB Mouse With Measure Function, C-Mount	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
494	A59.3507	HDMI+Type-C Digital Camera, 2.0M, Mouse Measure	2.0M, 1/2.8" CMOS, 1920*1080@60FPS, HDMI+Type C UVC Output To PC, With Measuring Software S-Eye, USB2.0 Connect To U Disk + Mouse Control, HDMI Interface Measure Menu, Power By Type C 5-12V DC	1 pcs / 1 ctn 33*24.5*9cm, 13.5/0Kgs			
495	A59.3513	4K HDMI+USB+U Disk, Mouse Control Measuring Digital Camera, 12M	12M, 3840x2160@60FPS, HDMI+UVC Output To PC With S-Eye Software, USB3.0 To U Disk, USB2.0 To Mouse + HDMI Measuring Function Menu, DC12V/2A Power Adapter 100~240V	1 pcs / 1 ctn 33*24.5*9cm, 13.5/0Kgs			
496	A59.3815	HDMI+USB+SD Digital Camera, 5.0M, Mouse Measure	5.0M, 1/2.5", 1080P, HDMI+USB2.0+SD, HDMI Panel Control By Mouse, C-Mount	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
497	A59.4231	HDMI+USB+TF Digital Camera, 48M, USB Measure	48M, Sony 1/2.9", Output HDMI+USB+TF Card, HDMI Resolution Photo 48M 8000x6000, Video 4K@25fps, 2K@25fps 1080p@60fps, Digital Zoom 5.0x, USB2.0 Resolution 2.0M, Photo 1920*1080, HDMI+USB Output At Same Time, USB Measuring ImageView, Input DC12V, C/CS Mount, HDMI Control Panel By Mouse	1 pcs / 1 ctn 16*10*8cm, 0.5/0.2Kgs			
498	A59.4235	HDMI+USB+TF Digital Camera, 60M, USB Measure	60M, 1/2.3", Output HDMI+USB+TF Card, HDMI Resolution Photo 60M, Video 2K@30FPS, 1080p@60FPS, USB Resolution 2.0M, HDMI+USB Output At Same Time, USB Measuring ImageView, Input DC12V, C/CS Mount, HDMI Control Panel By Mouse	1 pcs / 1 ctn 16*10*8cm, 0.5/0.2Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
499	A59.2223-5M	HDMI+WIFI+SD Digital Camera, 1080p/5.0M, HDMI Measure	5M, 1/1.8", 1080P, HDMI+WIFI+SD, Incl. USB Mouse/USB WiFi Adapter, HDMI & WIFI Output Same Time, HDMI Panel Control By Mouse, C-Mount	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
500	A59.2243-8MA	4K HDMI+USB+WIFI +LAN+SD Digital Camera, 8.0M, Mouse Measure	4K, 8.0M, 1/1.8", WIFI+HDMI+USB+WAN+SD卡, FPS 30@3840x2160 (HDMI/LAN/USB), HDMI UI + USB Mouse Measure, USB/WIFI Measure, C-Mount	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
501	A59.3508	WIFI+HDMI+USB +WAN, Mouse Control Measuring Digital Camera, 8.0M	8.0M, 1/1.8" CMOS, 1920*1080@60FPS, WIFI+HDMI+Type C UVC Output To PC+WAN, With Measuring Software S-Eye, USB2.0 Connect To U Disk + Mouse Control, HDMI Interface Measure Menu, Power	1 pcs / 1 ctn 33*24.5*9cm, 13.5/0Kgs			

LCD & Digital Camera For Microscope






502	A59.2259	11.6" HDMI LCD Monitor, For A59.2240.41.42	11.6" Panasonic IPS LCD, 1920x1080P, Contrast 1000:1, For A59.2217-A59.2223	10 pcs / 1 ctn 53*33*37cm, 18/16Kgs			
503	A59.2260	13.3" HDMI LCD Monitor, For A59.2240.41.42	13.3" Panasonic IPS LCD, 1920x1080P, Contrast 1000:1, For A59.2217-A59.2223	10 pcs / 1 ctn 53*33*37cm, 18/16Kgs			
504	A59.3509	9" LCD Digital Camera, USB2.0 +Mouse Measure, 5.0M	5.0M, 9" LCD, 1/2.8" CMOS, Screen Resolution 1024*600, USB2.0 Real Time Image Output, TF Card, USB HDMI Mouse Control Menu For Measure, DC12V/2A, 100~240V.	1 pcs / 1 ctn 27*20*7.5cm, 0.6/0.5Kgs			
505	A59.3550	11.6" LCD Digital Camera, HDMI+Mouse Measure, 5.0M	5.0M, 11.6" LCD, 1/2.8" CMOS, Screen Resolution 1920*1080 @60FPS, HDMI+USB Mouse Control Menu For Measure, DC12V/2A, 100~240V	1 pcs / 1 ctn 33*24.5*9cm, 13.5/0Kgs			
506	A59.5105	5" LCD Digital Camera, 0.3M, AV Output	5" LCD Resolution 800x480, 0.3M 1/3" CMOS, 648x488 Pixels, AV Output, Power by 5V DC USB Adapter, Mount to 23.2mm Eyepiece Tube	40 pcs / 1 ctn 59*38.5*30.5cm, 12/10Kgs			
507	A59.5107	7" LCD Digital Camera, 0.6M, AV Output	7" LCD Resolution 800x480, 0.6M 1/3" CMOS, 976x592 Pixels, AV Output, Power by 5V DC USB Adapter, Mount to 23.2mm Eyepiece Tube	20 pcs / 1 ctn 51*36*48.5cm, 12/10Kgs			
508	A59.5109	9" LCD Digital Camera, 0.6M, AV Output	9" LCD Resolution 800x480, 0.6M 1/3" CMOS, 976x592 Pixels, AV Output, Power by 5V DC USB Adapter, Mount to 23.2mm Eyepiece Tube	10 pcs / 1 ctn 57*31*36.5cm, 10/9Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
509	A59.5117	7" LCD Digital Camera, 2.0M, USB+AV+TF Card	7" IPS LCD Resolution 1024x760, 2.0M 1/2.75" CMOS, 1920x1080 Pixels, 30 Fps, TF Card Max 32G, AV Output, USB2.0 Output, With Measuring Software, Support Win/iMac, Eyepiece Tube 23.2mm Adapter	20 pcs / 1 ctn 51*36*48.5cm, 12/10Kgs			
510	A59.5119	9" LCD Digital Camera, 2.0M, USB+AV+TF Card	9" LCD Resolution 1024x760, 2.0M 1/2.75" CMOS, 1920x1080 Pixels, 30 Fps, TF Card Max 32G, AV Output, USB2.0 Output, With Measuring Software, Support Win/iMac, Eyepiece Tube 23.2mm Adapter	20 pcs / 1 ctn 51*36*48.5cm, 12/10Kgs			
LCD Pad Digital Camera For Microscope							
511	A59.3520	9.7" Touch Screen Android Pad Digital Camera, 5.0M	5.0M, 9.7" LCD, Resolution 2048*1536, Android 5.1, CPU 4 Core 1.8G, RAM 2G ROM 8G, HDMI+TF Card Output, Touch Screen USB Mouse Control, Installed Measure Software, Chinese/English Manual,	1 pcs / 1 ctn 34*27*11cm, 1.5/1Kgs			
512	A59.3521	10.5" Touch Screen Android Pad Digital Camera, 8.0M	8.0M, 10.5" Touch Screen Android Pad, 1/1.8" CMOS, 4K 3840*2158@30FPS, Built-in Memory RAM 2G ROM 16G, Built-in Android 11 System With S-Eye Measure Software, Control by Touch Screen & USB Mouse+Keyboard, Built-in 1080P AV Output	1 pcs / 1 ctn 34*27*11cm, 1.5/1Kgs			
513	A59.3551	11.6" Touch Screen Win10 Tablet Digital Camera, 16.0M	11.6" IPS LCD, 16M 1/2.3" CMOS, LCD Resolution 1920*1080, Win10 64 Tablet, 2.4G/5G WIFI, Bluetooth, USB3.0, HDMI, Ethernet, Type C Power Supply DC12V/2A, Incl. S-EYE Measure Software	1 pcs / 1 ctn 33*24.5*9cm, 13.5/0Kgs			
514	A59.0921	10.1" Touch Screen Android Tablet Digital Camera, 4.0M, HDMI+USB/HDMI+WIFI	10.1" Touch Screen Android Pad, 4.0M COMS, 2688x1520@15fps, 1920x1080@30fps, Built-in Android System With Measure Software, 16G SD Card, HDMI+USB/HDMI+WIFI Synchronous Output, Support Win/iMac	1 pcs / 1 ctn 41.6*30.2*23.6cm, 5/2.5Kgs			
Auto Focus Digital Camera For Microscope							
515	A59.0912	Auto Focus 4K HDMI+WIFI+USB Camera, Mouse Measure	4K 12.0M 1/1.7" CMOS, HDMI+USB2.0 输出, 60@4000x3000, TF Card Slot, C-Mount, Auto Focus, With Measurement Function	0 pcs / 0 ctn 26*20*10cm, 1.5/1Kgs			
516	A59.2239	Auto Focus 4K HDMI+WIFI+USB Camera, Mouse Measure	4K 8.0M, 1/1.8", HDMI+WIFI+USB, 30@3840x2160(HDMI), 30@3840x2160(WLAN), 30@3840x2160(USB), Incl. USB Mouse/USB WiFi Adapter, HDMI & WIFI Output Same Time, HDMI Mouse Control	0 pcs / 0 ctn 26*20*10cm, 1.5/1Kgs			
517	A59.2240-5M	Auto Focus HDMI+WIFI+SD Camera, Mouse Measure	5M, 1/1.8", 1080P, HDMI+WIFI+SD, 30@1920x1080 (HDMI), 25@1920x1080 (WiFi), Incl. USB Mouse/USB WiFi Adapter, HDMI & WIFI Output Same Time, HDMI Panel Control By Mouse, C-Mount, AutoFocus	0 pcs / 0 ctn 26*20*10cm, 1.5/1Kgs			
518	A59.3811	Auto Focus HDMI+USB+SD Digital Camera, 6.0M, Mouse Control	6.0M, 1/2.8", 1080P, HDMI+USB2.0+SD, HDMI Panel Control By Mouse, C-Mount, AutoFocus	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			






No.	Cata. No.	Product	Specification	Package			Picture
519	A59.3820	3D Super EDF Auto-Focus Measuring HDMI Digital Camera	2.0M, 1080P 60fps, 1/2", HDMI+Ethernet Output, USB3.0 Output To U Disk, Built-in 2D+3D Measure Software, Control By USB Mouse & Keyboard, C-Mount, DC12V	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
Light Source For Microscope							
520	A56.5110	Attachable LED Light Source, Rechargeable	LED 0.75W, Brightness Adjustable, Recharge Battery 3.7V, 400mAh, Power By DC3.5 Charger	100 pcs / 1 ctn 42*47*35cm, 0/0Kgs			
521	A56.2103-A	LED Ring Light, 56 LEDs, Dia.63mm	56 LEDs, 4.5W, Inner Dia.63mm, Light Path Dia.58mm, AC 90~240V	0 pcs / 0 ctn 19*11*5cm, 0.3/0Kgs			
522	A56.1211-B-56	LED Ring Light, 56/64/144 LEDs, Dia.61mm	56 bulbs, Continous Intensity Adjustable	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
523	A56.1211-B-144	LED Ring Light, 56/64/144 LEDs, Dia.61mm	144 bulbs, Continous Intensity Adjustable	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
524	A56.3902-C	LED Ring Light, 40 LEDs, Install Dia.46.6mm, For Mono Zoom Microscope	40 LEDs, Inner Dia. 28mm, Install Dia. 28-46.5mm, Separate Brightness Adjustment	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
525	A56.2418	LED Polarizing Ring Light, 78 LEDs, 7W	78 LEDs, Including Polarizing & Analyzer Lens, Brightness@100mm >24000lx @6400K, Connect Dia.70mm, 7W, Packing List: Adapter 1 pc, LED Ring Light 1 pc, M48X1.5Plastic Ring Adapter 1 pc	1 pcs / 1 ctn 22*26*7cm, 0.5/0.4Kgs			
526	A56.3910	Polarizing LED Ring Light	120 LEDs, Light Path Dia.35mm, Out Dia. 90mm, Install Inner Dia. 62mm	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
527	A56.2617-10	LED Cool Light Illuminator	10W-LED,6000K 60000Lx	1 pcs / 1 ctn 28.5*23*25.5cm , 2.9/2.4Kgs			
528	A56.2617-20	LED Cool Light Illuminator	20W-LED,6000K,90000Lx,Cooling Fan	1 pcs / 1 ctn 28.5*23*25.5cm , 2.9/2.4Kgs			










No.	Cata. No.	Product	Specification	Package			Picture
529	A56.2617-30	LED Cool Light Illuminator	30W-LED,6000K,150000Lx,Cooling Fan	1 pcs / 1 ctn 28.5*23*25.5cm , 2.9/2.4Kgs			
530	A56.2618	Cold Light Source, Halogen 150W	21V150W/Philips13629,3000K, 100000Lx, Cooling Fan	1 pcs / 1 ctn 28.5*23*25.5cm , 2.7/2.1Kgs			
531	A56.4801	Dual Pipe Desktop LED Cold Light, 3Wx2	1.Voltage:AC90~240V 50Hz-60Hz 2.Power of bulb:LED 3W*2 3.Color temperature:3500~4500K 4.Life of bulb:50000h 5.Illumination(200mm):6000Lux*2	1 pcs / 0 ctn 35*17*26cm, 1.5/0.7Kgs			
532	A56.4821	LED Cold Light Source, 20W	LED 20Wx1, Color Temperature 5000±500K, Brightness 1600000 Lux, LED Lamp Life 20000H, Power AC90~240V 50Hz-60Hz, Size 95x175x155mm	1 pcs / 0 ctn 27*20*19cm, 2.5/1.6Kgs			
533	A56.4822-80W	LED Cold Light Source, 80W	1. Voltage: AC90~240V 50HZ-60Hz 2. Bulb: LED 80W 3. Bulb life: 20000h 4. Plane Illumination: ≥4,000,000Lux 5.C.T.:5000±500K 6.Size: 285*270*120mm	0 pcs / 0 ctn 40*40*26cm, 5/3.9Kgs			
534	A56.4832	Halogen Cold Light Source, 150W, Single Bulb, Single Hole	1. Voltage: AC110V or 220V 50Hz~60Hz 2. Bulb: halogen 24V150W 3. Bulb life: 100h 4. Plane Illumination:≥1,600,000Lux 5.C.T.:3600K 6.Size: 285*270*120mm	0 pcs / 0 ctn 35*17*26cm, 3.8/3Kgs			
535	A56.4833-150W	Halogen Cold Light Source, 150W/250W, Double Bulb, Double Hole	1. Voltage: AC110V or 220V 50Hz~60Hz 2. Bulb: halogen 24V150W 3. Bulb life: 100h 4. Plane Illumination:≥1,600,000Lux 5.C.T.:3600K 6.Size: 285*270*120mm	0 pcs / 0 ctn 38*38*26cm, 5/4.2Kgs			
536	A56.4844-250	XenonLight Source 250W	1. Voltage: AC110V or 220V 50Hz~60Hz 2. Bulb: xenon 250W 3. Bulb life: 500h 4. Plane Illumination:5,800,000Lux 5.C.T.:5600K 6.Size: 285*270*120mm	0 pcs / 0 ctn 41*39*25cm, 7.5/6.5Kgs			
Micrometer, Retical, Cross							
537	A58.2001-901	Stage Micrometer	100 DIV*0.01mm=1mm	200 pcs / 1 ctn 32*34*38cm, 8/6Kgs			
538	A58.2002-935	Eyepiece Micrometer	X 10*10 DIV * 0.1mm=10mm, Y Line	200 pcs / 1 ctn 32*34*38cm, 8/6Kgs			












No.	Cata. No.	Product	Specification	Package			Picture
Heating Stage							
539	A54.6001-2	Temperature Controlled Heating Stage, 10°C~320°C	Room Temp. +10 oC—320 oC Heating Plate Φ65mm Outside Dia.Φ 95mm 220V input, 36V DC output	1 pcs / 1 ctn 28*35*40cm, 3.5/2.5Kgs			
540	A54.9602	Temperature Controlled Heating Stage, 10°C~70°C	Temperature Range : 15 to 70 °C . Digital Display, Preset Temperature Control, Measuring Range ≤± 0 .2° C , Working Condition: 15° C - 28° C	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
541	A54.9605	Temperature Controlled Heating Stage, 10°C~50°C	Temperature Range : 10 to 50° C . LED Digital Display, Preset Temperature Control, Heating Stage: 5.5*9.5cm, Main Body 10*14*16.8cm	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
Eyepiece & Mobile Phone Adapter							
542	A55.2001	Eyepiece Adapter Set, 1.0x	C-Mount To Eyepiece Tube, Including 1.0x Adapter Tube 23.2mm, Ring 30mm,30.5mm	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
543	A55.2002	Eyepiece Adapter Set, 0.5x	C-Mount To Eyepiece Tube, Including 0.5x Adapter Tube 23.2mm, Ring 30mm,30.5mm	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
544	A55.2003	Eyepiece Adapter Set, 0.5x	C-Mount To Eyepiece Tube, Including 0.5x Adapter Tube 23.2mm, Ring 30mm,30.5mm	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
545	A55.2202-A50	Eyepiece Tube Adapter, for Microscope	Fit to 1/2" ~ 2/3" size sensor 0.5X Magnification Adjustable Focal length C-Mount to Dia.23.2mm eye tube	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
546	A55.2202-A37	Eyepiece Tube Adapter, for Microscope	Fit to 1/4" ~ 1/3" size sensor 0.37X Magnification Adjustable Focal length C-Mount to Dia.23.2mm eye tube	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
547	A55.2202-A75	Eyepiece Tube Adapter, for Microscope	Fit to 1/1.8" ~ 1" size sensor 0.75X Magnification Adjustable Focal length C-Mount to Dia.23.2mm eye tube	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
548	A55.2204-A37	Eyepiece Tube Adapter, For Telescope	Fit to 1/4" ~ 1/3" size sensor 0.37X Magnification Adjustable Focal length C-Mount to Dia.31.75mm Telescope eye tube	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
549	A55.2204-A50	Eyepiece Tube Adapter, For Telescope	Fit to 1/2" ~ 2/3" size sensor 0.5X Magnification Adjustable Focal length C-Mount to Dia.31.75mm Telescope eye tube	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
550	A55.2204-A75	Eyepiece Tube Adapter, For Telescope	Fit to 1/1.8" ~ 1" size sensor 0.75X Magnification Adjustable Focal length C-Mount to Dia.31.75mm Telescope eye tube	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			
551	A55.9011	Univeral Holder for Mobile Phone to Microscope	Plastic Board 100*190*5mm to Hold Mobile Phone, Magnet Slide, Eyepiece 10x, Eyepiece Tube Dia.23.2mm For Biological Microscope, Adapt Ring 30, 30.5 For Stereo Microscope	20 pcs / 1 ctn 58*31*45cm, 15/12Kgs			
552	A55.9009	Universal Holder For Mobile Phone	Mobile Phone Width 52~100mm, Eyepiece Dia25~48mm of All Kinks Microscope, Binocular, Telescope Connection	100 pcs / 1 ctn 28*48*43cm, 9/0Kgs			
553	A55.9010	Universal Holder For Mobile Phone	Mobile Phone Width 55~98mm, Eyepiece Dia26.4~46.4mm of All Kinks Microscope, Binocular, Telescope Connection	100 pcs / 1 ctn 28*48*43cm, 9/0Kgs			

Operating Surgical Microscope

554	A41.1901-A	Operating Microscope, One Head 90°, Manual Step Zoom, 3x~12x, For ENT, Surgery	Straight Binocular 6x, Zoom Body 0.6x, 1.0x, 1.6x, Objective F=200mm, 300mm, Total Magnification 3x~12x, Built-in Blue, Orange Filter, Fine Focusing Range 10mm, 10W LED, Column + Five Star Base With Caster	1 pcs / 2 ctn 67*48*33cm, 49/45Kgs			
555	A41.1901-B	Operating Microscope, One Head 45°, Manual Step Zoom, 3x~12x, For Dental, Ophthalmology	45° Binocular 6x, Objective F=200mm, 300mm, Total Magnification 3x~12x, Built-in Blue, Orange Filter, Fine Focusing Range 10mm, 10W LED, Column + Five Star Base With Caster	1 pcs / 2 ctn 67*48*33cm, 49/45Kgs			
556	A41.1901-C	Operating Microscope, One Head 90°, Manual Step Zoom, 3x~12x, Table Mount, For ENT, Surgery, Table	Straight Binocular 6x, Objective F=200mm, 300mm, Total Magnification 3x~12x, Built-in Blue, Orange Filter, Fine Focusing Range 10mm, 10W LED, Table Mount Without Stand	1 pcs / 1 ctn 67*48*33cm, 16/12Kgs			
557	A41.1901-CZ	Operating Microscope, One Head 90°, Manual Step Zoom, 3x~12x, Pole Stand, For ENT, Surgery, Column	Straight Binocular 6x, Objective F=200mm, 300mm, Total Magnification 3x~12x, Built-in Blue, Orange Filter, Fine Focusing Range 10mm, 10W LED, Customize Column Mount, Without Stand	1 pcs / 1 ctn 67*48*33cm, 15/12Kgs			
558	A41.1901-D	Operating Microscope, One Head 45°, Manual Step Zoom, 3x~12x, Table Mount, For ENT, Surgery, Table	45° Binocular 6x, Objective F=200mm, 300mm, Total Magnification 3x~12x, Built-in Blue, Orange Filter, Fine Focusing Range 10mm, 10W LED, Table Mount Without Stand	1 pcs / 1 ctn 67*48*33cm, 15/12Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
559	A41.1902-C	Operating Microscope, One Head 90°, Manual Step Zoom, 3x-12x, For ENT, Surgery	Straight Binocular 6x, Zoom Body 0.6x, 1.0x, 1.6x, Objective F=200mm, F=300mm, With Column Stand, 10W LED Lamp, With 50:50 Beamsplitter, With C-Mount	1 pcs / 2 ctn 87*66*24cm, 63/55Kgs			
560	A41.1902-D	Operating Microscope, One Head 45°, Manual Step Zoom 3x-12x, For Dental, Ophthalmology	45° Inclined Binocular 6x, Objective F=200mm, F=300mm, With Column Stand, 10W LED Lamp, With 50:50 Beamsplitter, With C-Mount	1 pcs / 2 ctn 87*66*24cm, 59/55Kgs			
561	A41.1903	Operating Microscope, One Head 0-180°, Manual Step Zoom, 4.8x-16x, For Dental, ENT, Ophthalmology	8x Tilting 0-180° Binocular, Eyepiece 12.5x, Zoom Body 0.6x, 1.0x, 1.6x, Objective F=200mm, 250mm, Total Magnification 4.8x-16x, Built-in Blue, Orange Filter, Fine Focusing Range 10mm, 10W LED, Column + Five Star Base	1 pcs / 2 ctn 86*67*24cm, 59/50Kgs			 HOT
562	A41.1921	Operating Microscope, One Head 0-200°, Manual 6 Step, 2.7x, 4.5x, 7.2x, 10.8x, 18x, 27x, For Dental, ENT	One Binocular Head, Tilttable 0-200°, Interpupillary 55-75mm, Diopter +/-6D, Eyepiece 12.5x, 6 Steps Manual Zoom, Manification 2.7x, 4.5x, 7.2x, 10.8x, 18x, 27x, 12.6x, View Field Dia. 74.1-7.4mm, Objective F=250mm, LED Lamp, 10W, 50:50 Beamsplitter	1 pcs / 3 ctn 595*460*230cm , 111/100Kgs			
563	A41.1925	Operating Microscope, One Head 0-200°, Manual Zoom, 2.5x-21x, For Dental, ENT	One Binocular Head, Tilttable 0-200°, Interpupillary 55-75mm, Diopter +/-6D, Eyepiece 12.5x, Manual Continous Zoom 1:6, Manification 2.5x-21x, View Field Dia. 78.6-9mm, Objective F=250mm, Head X/Y Swing +/-45°, LED Lamp, 10W, 50:50 Beamsplitter	1 pcs / 3 ctn 595*460*230cm , 111/100Kgs			
564	A41.1931	Operating Microscope, Dual Head, Manual 3 Step, 6x10x16x, 2 Function Foot Switch, For Ophthalmology	Dual Head, Both Inclined 45°, Interpupillary 55-75mm, Diopter +/-6D, Eyepiece 12.5x, 3 Steps Manual Zoom, Manification 6x10x16x, View Field Dia. 34, 20, 13mm, Objective F=200mm, Assistant Head 6x10x16x, 2nd Objective F=125mm, 3.2x, 5.0x, 7.5x, 10.0x, 15.0x, 20.0x, View Field Dia. 24.5-12.2mm, Objective F=200mm, LED Lamp, 10W, 50:50 Beamsplitter	1 pcs / 3 ctn 595*460*230cm , 104/90Kgs			
565	A41.1932	Operating Microscope, Dual Head, Manual 3 Step, 6x10x16x, 6 Functions Foot Switch, For Ophthalmology	Dual Head, Both Inclined 45°, Interpupillary 55-75mm, Diopter +/-6D, Eyepiece 12.5x, 3 Steps Manual Zoom, Manification 6x10x16x, View Field Dia. 34, 20, 13mm, Objective F=200mm, Assistant Head 6x10x16x, 2nd Objective F=125mm, 3.2x, 5.0x, 7.5x, 10.0x, 15.0x, 20.0x, View Field Dia. 24.5-12.2mm, Objective F=200mm, LED Lamp, 10W, 50:50 Beamsplitter	1 pcs / 3 ctn 595*460*230cm , 104/90Kgs			
566	A41.1935	Operating Microscope, Dual Head, Motor Zoom, 4.5x-27.3x, 8 Function Foot Switch, For Ophthalmology	Dual Head, Main Head Tilttable 30-90°, Interpupillary 55-75mm, Diopter +/-5D, Eyepiece 12.5x, Motorized Zoom, Manification 1:6 4.5x-27.3x, View Field Dia. 44-7.7mm, Objective F=200mm, Assistant Head 30°, Manual 3 Step Zoom, 3x-12x, View Field Dia. 34, 20, 13mm, Objective F=200mm, LED Lamp, 10W, 50:50 Beamsplitter	1 pcs / 4 ctn 595*460*230cm , 165/160Kgs			
567	A41.1941	Operating Microscope, Dual Head, Step Zoom, 6x10x16x, 2 Function Foot Switch, For Orthopedics	Dual Head 180° Opposite, Both Inclined 45°, Interpupillary 55-75mm, Diopter +/-6D, Eyepiece 12.5x, 3 Steps Manual Zoom, Manification 6x10x16x, View Field Dia. 34, 20, 13mm, Objective F=200mm, Assistant Head 6x10x16x, 2nd Objective F=125mm, 3.2x, 5.0x, 7.5x, 10.0x, 15.0x, 20.0x, View Field Dia. 24.5-12.2mm, Objective F=200mm, LED Lamp, 10W, 50:50 Beamsplitter	1 pcs / 0 ctn 595*460*230cm , 104/100Kgs			
568	A41.1942	Operating Microscope, Dual Head, Manual Step Zoom, 6x10x16x, 6 Function Foot Switch, For Orthopedics	Dual Head, Both Inclined 45°, Interpupillary 55-75mm, Diopter +/-6D, Eyepiece 12.5x, 3 Steps Manual Zoom, Manification 6x10x16x, View Field Dia. 34, 20, 13mm, Objective F=200mm, Assistant Head 6x10x16x, 2nd Objective F=125mm, 3.2x, 5.0x, 7.5x, 10.0x, 15.0x, 20.0x, View Field Dia. 24.5-12.2mm, Objective F=200mm, LED Lamp, 10W, 50:50 Beamsplitter	1 pcs / 0 ctn 595*460*230cm , 104/100Kgs			
569	A41.1945	Operating Microscope, Dual Head, Motor Zoom, 2x-30x, 8 Function Foot Switch, For Orthopedics	Dual Head 180° Opposite, Main Head Tilttable 30-90°, Interpupillary 55-75mm, Diopter +/-5D, Eyepiece 12.5x, Motorized Zoom, Manification 1:6 4x-24x, View Field Dia. 50-8.2mm, Objective F=200mm, Assistant Head 45°	1 pcs / 0 ctn 595*460*230cm , 165/160Kgs			

No.	Cata. No.	Product	Specification	Package			Picture
570	A41.1950	Operating Microscope, Dual Head, Motor Zoom, 3x~16x, 8 Function Foot/Handle, For Neurosurgery	Dual Head Binocular, Eyepiece 12.5x, Main Binocular Variable Angle 30-90°, Diopter Adjustable +/-6D, IPD 55~75mm, Motorized Continuous 1:6 Zoom 3x~16x, Built-in Beam Splitter, Assistant's Binocular 0° Can Swing 90°	1 pcs / 4 ctn 594*460*230cm 169/140Kgs			 HOT
571	A41.1951	Operating Microscope, Dual Head, Motor Zoom, 1.8x~21x, 12 Function Handle, For Neurosurgery	Dual Head Binocular, Eyepiece 12.5x, Main Binocular Variable Angle 0-200°, Diopter Adjustable +/-6D, IPD 55~75mm, Motorized Continuous Zoom 1.8x~21x, Built-in Beam Splitter, Assistant's Binocular 30° Can Swing 90°	1 pcs / 5 ctn 595*460*230cm 169/160Kgs			
572	A41.1952	Operating Microscope, Dual Head, Motor Zoom, 1.8x~19x, 10 Function Handle, For Neurosurgery	Dual Head Binocular, Eyepiece 12.5x, Main Binocular Variable Angle 0-200°, Diopter Adjustable +/-6D, IPD 55~75mm, Motorized Continuous Zoom 1.8x~19x, Built-in Beam Splitter, Assistant's Binocular Can Swing 90° or 0°	1 pcs / 2 ctn 67*48*33cm, 49/45Kgs			
573	A41.3403	Operating Microscope, One Head 45°, Manual Step Zoom, 5.3x,8x,12x, 2 Directions Foot Switch, For Neurosurgery	One Binocular, Eyepiece 12.5x/18mm, IPD 50-70mm, Diopter ±6, Objective F=200mm, Main Microscope 5.3x/37mm, 8x/25mm,12x/16.7mm, W.D.190mm, Halogen 12V100W, 6°+0° Coaxial Illumination ≥35000 Lux, Arm	1 pcs / 1 ctn 0*0*0cm, 41/0Kgs			
574	A41.3405	Operating Microscope, Dual Binocular, 5 Steps 4x,6x,10x,16x,25x, 6 Function Foot Switch, For Neurosurgery	Dual Binocular, Eyepiece 12.5x/18mm, IPD 50-70mm, Diopter ±6, Objective F=200mm, Main Microscope 5 Step 4x,6x,10x,16x,25x, Assistant Microscope 6x,10x,16x, W.D.170mm, Halogen 12V100W, 6°+0° Coaxial	1 pcs / 6 ctn 0*0*0cm, 207/190Kgs			
575	A41.3421	Operating Microscope, Dual Head, Motor Zoom, 4.6x~27x, 8 Directions Foot Switch, For Phththalmology	Dual Binocular, Eyepiece 12.5x/18mm, IPD 50-70mm, Diopter ±6, Objective F=200mm, Main Microscope Motorized/Manual Zoom 4.5x-25x, Assistant Microscope 6x,10x,16x, W.D.170mm, Halogen 12V100W, 6°+0°	1 pcs / 6 ctn 0*0*0cm, 205/190Kgs			
Motorized Measuring Microscope							
576	M11.5805	XY Motorized Measuring Microscope, Mono Zoom Body	0.6x-5x, 1:8.3, W.D.145mm, Focus Range 50mm, 5W LED Ring Light, XY Motorized Stage, Size 255x205mm, Moving 100x100mm, 9.0M USB3.0 GS Shutter Digital Camera, PC Dell i5 64G 256G+1T, 2G Graphic, 27" 4K,	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			 NEW
577	M11.5810	XY Motorized Measuring Microscope, Binocular Head Zoom Body	Binocular, WF10x/22, 0.8x-5x, 1:6.3, W.D.115mm, Focus Range 50mm, 5W LED Ring Light, XY Motorized Stage, Size 255x205mm, Moving 100x100mm, 9.0M USB3.0 GS Shutter Digital Camera, PC Dell i5 64G 256G+1T, 2G	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			 NEW
578	M12.5810	Motorized Biological Microscope, BF, XYZ Motorized	Biological + XYZ Motor Stage Trinocular, WF10X/22 Dia.30mm, Quintuple, Infinity Plan Semi-APO FL4x10x20x40x100x, XYZ Motorized Stage 255x210mm. Moving 100x100x42mm, Swing-out Condenser,	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			 NEW
579	M12.5820	Motorized Biological Microscope, BF, XYZ Motorized	Biological + XYZ Motor Stage + Motor Nosepiece Trinocular, EW10X/25 Dia.30mm, 6 Holes Motor Noespiece, Infinity Plan NIS60 Semi-APO 4x10x20x40x100x, XYZ Motorized Stage 255x210mm,	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			 NEW
580	M12.5850	Motorized Biological Microscope, BF, XYZ Motorized	Biological + XYZ Motor Stage + Motor Nosepiece Trinocular, PL10X/25, Motorized 6 Holes Motor Noespiece With DIC/PL Slot, Infinity Plan Semi-APO FL4x10x20x40x100x, XYZ Motorized	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			 NEW

No.	Cata. No.	Product	Specification	Package			Picture
581	M13.5820	Motorized Metallurgical Microscope, BF+DF+PL, Semi-APO, DIC, XYZ & Nosepiece	Metallurgical + XYZ Motor Stage + Motor Nosepiece Tilting Trinocular 0-35°, EW10x/25mm, Manual Sextuple Nosepiece, NIS45 Infinity Plan Semi-APO BD 5x10x, LWD 20x50x100x, XYZ Motorized Stage	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			 NEW
582	M13.5850	Motorized Metallurgical Microscope, BF=DF+PL, Semi-APO, DIC, XYZ & Nosepiece	Upright TR Metallurgical B+PL+XYZ Motor + Motor Nosepiece Trinocular, PL10x/25, 6 Holes Motor Noespiece With DIC/PL Slot, Infinity Plan Semi-APO BD DIC5x10x20x50x100x, XYZ Motorized	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			 NEW
583	M16.5810	Motorized Fluorescent Microscope, BF/FL, XYZ Motorized	Fluorescent + XYZ Motor Stage Trinocular, WF10x/22 Dia.30mm, Quintuple, Infinity Plan Semi-APO FL4x10x20x40x100x, XYZ Motorized Stage 255x210mm, Moving 85x70x42mm, Swing-out Condenser,	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			 NEW
584	M16.5820	Motorized Fluorescent Microscope, BF/FL, XYZ+Nosepiece	Fluorescent + XYZ Motor Stage + Motor Nosepiece Trinocular, EW10X/25 Dia.30mm, 6 Holes Motor Noespiece, Infinity Plan NIS60 Semi-APO 4x10x20x40x100x, XYZ Motorized Stage 255x210mm,	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			 NEW
585	M16.5850	Motorized Fluorescent Microscope, BF/FL, XYZ Motorized	BF/FL + XYZ Motor Stage + Motor Nosepiece Trinocular, PL10x/25, Motor 6 Holes Noespiece With DIC/PL Slot, Infinity Plan Semi-APO FL4x10x20x40x100x, XYZ Motorized Stage 255x210mm,	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			 NEW
586	M20.3820	3D Auto-Focus Super EDF Measuring Digital Microscope	16x-160x, 2.0M, 1080P 60fps, 1/2", HDMI+Ethernet Output, USB3.0 Output To U Disk, Built-in 2D+3D Measure Software, Control By USB Mouse & Keyboard, C-Mount, DC12V	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			 NEW
587	M20.5850	All-in-One Industry Microscope, Motorized XYZ +Motorized Nosepiece +Marco Camera	6 Holes Motor Nosepiece, Infinity Plan LWD Metallurgical Semi-APO 5x10x20x50x100x, XYZ Motorized Stage Moving Range 50x50x20mm, Built-in LED Illumination, 9.0M USB3.0 GS Shutter Digital Camera, 1" CMOS,	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			 NEW
588	M20.5880	Gantry type Industry Microscope, Motorized XYZ +Motorized Nosepiece +Marco Camera	6 Holes Motor Nosepiece, Infinity Plan LWD Metallurgical 2x5x10x20x50x100x, XYZ Motorized Stage 550x570x65mm, Moving Range 300x300x20mm, Built-in LED Transmit/Reflect Illumination, 9.0M	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			 NEW
589	M20.5860	All-in-One Industry Microscope, Motorized XYZ +Motorized Nosepiece + Inclined Stand	12x-7000x, Built-in Motorized Zoom Body 0.5x-6x, Zoom Range 1:12, Intelligent Magnification Display, Manual + Auto Focus, Motor Z Move Range 120mm, Max Specimen Height 150mm, Vertical & +/-90° Inclined Stand, Stage	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			 NEW
Motorized Slide Scanner							
590	M30.1001	Digital Microscope Slide Scanner	6 Slides Scan Holder, 20x Plan Apochromatic Objective, N.A.=0.75. Automatic Focus, Automatic Search For Scanning Samples. LED Light With Fly-eye Lens, High Precision Full Closed Loop Drive Control System: X/Y/Z Axis	1 pcs / 1 ctn 0*0*0cm, 0/0Kgs			 NEW
591	M30.5810-A04	Full Auto Microscope Slide Scanner	1 Objective, 1 Camera, 4 Slides Scan --Infinity Plan 20x Objective --1 Micro Camera --2D Standard Software --2C Computer	0 pcs / 0 ctn 0*0*0cm, 0/0Kgs			 NEW